OSYnet-SafetyController, Best.-Nr. 77 2100 00

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OSYnet-SafetyController, order no. 77 2100 00

Safety controller for direct connection to the group bus (OSYnet). Usable in parallel with OSY-ControlCenter, per physical group. Intended use as redundancy for maintaining the system operation in case of OSY-ControlCenter fail.

- Controlling the base light call functions of one physical group (DIN VDE 0834)
- Control of room lamps and displays
- Autonomous and automatic switching between normal and safety operation
- 1 potential free digital output (change-over contact) for signalling the operation mode
- 1 potential free digital output (change-over contact) for fault signalling
- Digital output: internal: 24 VDC/60 mA, external: 24 VAC/1 A
- Max. current consumption: 230 mA

Note! The complete installation of the system is described in the technical manual.

Attention! The printed circuit board includes electrostatic sensitive components. Avoid touching.

A Mounting

Mounting on supporting rail or wall mounting in the electrical room of the ward.

A1: Mounting on supporting rail (35 mm)

Click the OSYnet-SafetyController 1 onto the supporting rail 3.

A2: Wall mounting

1. Push out the three fixing clips 2 until the holes for the screws 4 appear.
2. Screw the OSYnet-SafetyController 1 to the wall.

B Setting of address (126, 125)

For setting the address on the printed circuit board please use the code switches 1 to 7 with the values 1, 2, 4, 8, 16, 32, 64. Set the appropriate switch to ON.

At the OSYnet-SafetyController you have to set **address 126**. If you install a second OSYnet-SafetyController in the same physical group, you have to set **address 125** at the second OSYnet-SafetyController.

Do not change code switches 8 and 9! Factory set: OFF.

Reset

To reset the OSYnet-SafetyController you have to set code switch 10 to ON for one second, then back to OFF. A reset has to be made, if a code switch setting has been changed. Otherwise the new setting will not be activated.

B Control LED (red)

Light OFF:	OSYnet-SafetyController without power supply or defect
Flashing slowly: (1 s ON - 1 s OFF)	OSY-ControlCenter controls the group. OSYnet-SafetyController ready-to-operate (standby)
Flashing quickly: (200 ms ON - 200 ms OFF)	OSYnet-SafetyController controls the group.
Light ON for more than 10 s:	OSYnet-SafetyController not ready to operate.


C ESD protection cover (order no. 00 0276 53)


After you have connected the cables (see reverse side), you have to put the ESD protection cover on the OSYnet-SafetyController. (ESD = Electro Static Discharge)

D Dismantling from supporting rail

Put a screw driver into the visible hole of the black fixing clip. Then pull the clip upwards, until the OSYnet-SafetyController 1 comes loose from the supporting rail 3.

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Output for signalling the operation mode

OUT 1 closed & OUT 2 open = Normal operation

OUT 1 open & OUT 2 closed = Safety operation (SafetyController controls the group)

Jumper OUT 1/2 set = Internal power supply (24V DC, maximum 60 mA for triggering a relay)

Jumper OUT 1/2 not set = External power supply (24V AC or DC, maximum 1 A)

Output for fault signalling

Faults of group bus users

OUT 3 closed & OUT 4 open = no fault

OUT 3 open & OUT 4 closed = Fault!

Jumper OUT 3/4 set = Internal power supply (24V DC, maximum 60 mA for triggering a relay)

Jumper OUT 3/4 not set = External power supply (24V AC or DC, maximum 1 A)

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A Mounting

Optionally mounting on supporting rail or wall mounting in the electrical room of the ward.

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At the OSYnet-SafetyController you have to set address 126. If you install a second OSYnet-SafetyController in the same physical group, you have to set address 125 at the second OSYnet-SafetyController.

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Reset

To reset the OSYlink you have to set code switch 10 to ON for one second, then back to OFF. A reset has to be made, if a code switch setting has been changed. Otherwise the new setting will not be activated.

B Control LED (red)


Light OFF:	OSYnet-SafetyController without power supply or defect
Flashing slowly: (1s ON - 1s OFF)	OSY-ControlCenter controls the group. OSYnet-SafetyController ready-to-operate (standby)
Flashing quickly: (200 ms ON - 200 ms OFF)	OSYnet-SafetyController controls the group.
Light ON for more than 10 s:	OSYnet-SafetyController not ready to operate.


C ESD protection cover (order no. 00 0276 53)

After you have connected the cables (see reverse side), you have to put the ESD protection cover on the OSYlink. (ESD = Electro Static Discharge)

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Put a screw driver into the visible hole of the black fixing clip. Then pull the clip upwards, until the OSYlink 1 comes loose from the supporting rail 3.

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© Tunstall GmbH, Okoliten 66, D-48291 Telgte, www.tunstall.de