Steuermodul ZLB, Bestell-Nr. 19 0701 20

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Control module ZLB, order no. 19 0701 20

The control module ZLB is intended for connection to the ward bus of a CONCENTO^{PLUS} nurse call system. The room bus RAN is connected to the control module ZLB. Additionally a room lamp universal, order no. 77 0180 10 or 77 0181 10 can be connected to the control module ZLB.

System requirements of the CONCENTO^{PLUS} nurse call system for the installation of control module ZLB:

- Group electronic, order no.: 19 0700 10, firmware revision: 1.3.xxxx or higher.
- Terminals SD6 P, order no.: 19 0700 20, firmware revision: 1.3.xxxx or higher.
- ConLog^{PLUS} Management Software, order no. 19 0803 0x, software revision: 1.0.3.x or higher.

The control module ZLB controls all room functions in compliance with the German standard DIN VDE 0834. At the same time, the control module ZLB serves as a room lamp including light sections for all call types (red), staff presence (green), and a additional light section for calls in the WC area (white).

NOTE! The complete installation of the system is described in the technical manual.

CAUTION! The printed circuit board and the LED modules include electrostatic sensitive components. Avoid touching.

A Product description


For details refer to fig. C.

- | | |
|--------------------------------------|--|
| 1 - *Back box, 1-gang | 7 - Three LED modules |
| 2 - Connector, 8-pole, pluggable | 8 - Fastening screw for the PCB |
| 3 - Casing of the control module ZLB | 9 - Connector, 2-pole, pluggable |
| 4 - *Back box screws | 10 - Light dome with isolation insert |
| 5 - Printed circuit board PCB | * Not included with product delivery, please order separately. |
| 6 - DIP switch, 8-pole | |


B Dismantling

At the point of delivery the control module ZLB is assembled and must be dismantled as follows:

1. Remove the 8-pole connector [2] from the rear of the control module ZLB, see fig. C.
2. Lightly press the light dome [10] together from top and bottom.
3. Pull off the light dome [10].
4. Loosen the fastening screw [8] of the PCB and remove the PCB (first bottom, then top), see fig. C.

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
C Mounting

Wall mounting on a 1-gang back box.

1. Strip the connection cables to a suitable length. Remove the cable shield and drain wire up to the cable sheath. Strip the wires.
2. Connect the wires to the removed connectors [2] and [9] according to fig. E.
3. Set the 8-pole DIP switch according to section D and fig. D.
4. Guide the two wired connectors ([2], [9]) from behind through the opening in the rear of the casing [3].
5. Bolt the casing [3] to the back box [1] with the screws [4]. The arrow in the rear of the casing points toward the top.
6. Plug the wired 8-pole connector into the socket on the rear of the PCB [5].
7. Route the two wires of the wired 2-pole connector [9] through the slot in the PCB [5], so that the wires cannot be pinched, when the PCB [5] is inserted into the casing [3]. Now plug the 2-pole connector [9] into the socket on the PCB [5].
8. Insert the PCB [5] into the casing [3] (first top, then bottom) and screw in the fastening screw [8].
9. If the three LED modules [7] are not attached to the PCB [5], attach them now. The colour of an LED module can be identified by a colour marking on the surface of the LED module. For positioning of the colours refer to fig. D.
10. Perform a RAN test before commissioning the room installation, see section D DIP switch 3.
11. Press the light dome with isolating insert [10] onto the casing [3] until it latches in audibly.

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- 1 - *Back box, 1-gang
- 2 - Connector, 8-pole, pluggable
- 3 - Casing of the control module ZLB
- 4 - *Back box screws
- 5 - Printed circuit board PCB
- 6 - DIP switch, 8-pole

- 7 - Three LED modules
- 8 - Fastening screw for the PCB
- 9 - Connector, 2-pole, pluggable
- 10 - Light dome with isolation insert
- * Not included with product delivery.

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
D Setting the 8-pole DIP switch

Factory setting = All DIP switches are set to OFF.


DIP switch 1: Bus termination	
OFF	Control module ZLB is neither the first nor the last ward bus user.
ON	Control module ZLB is the first or the last ward bus user (= bus termination).
DIP switch 2: WC call cancelling	
OFF	WC calls can be cancelled only with a WC cancel button.
ON	WC calls can be cancelled with a WC cancel button or a presence button.
DIP switch 3: RAN test	
Test whether all devices on the room bus RAN have been connected correctly and are ready for operation. Perform this test before commissioning the room installation:	
1. Set DIP switch 3 to ON.	
✓	The LEDs of all devices on the room bus must flash. Devices not flashing have been incorrectly connected or are defective. NOTE: The LEDs of pear push switches and room lamp universal do not flash with this test.
2. Set DIP switch 3 back to OFF.	
3. Replace defective devices. Correctly connect devices that have not been correctly connected.	
DIP switches 4 – 8: RAN number (0 – 30)	
Here you have to set the number of devices on the room bus RAN. Devices on the room bus RAN = Switches, con- nection sockets, RAN interfaces. No pear push switch or other pluggable call device, no room lamp universal.	
DIP switch 4 ON = RAN number 1.	
DIP switch 5 ON = RAN number 2.	
DIP switch 6 ON = RAN number 4.	
DIP switch 7 ON = RAN number 8.	
DIP switch 8 ON = RAN number 16.	
Other RAN numbers are set by mathematical addition of these DIP switches. Example:	
RAN number 3 = DIP switch 4 (for RAN number 1) ON and DIP switch 5 (for RAN number 2) ON.	
RAN number 5 = DIP switch 4 (for RAN number 1) ON and DIP switch 6 (for RAN number 4) ON.	
Factory setting: RAN number 0 = all DIP switches OFF.	

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Technical data

Power supply	24 V DC
Standby current consumption	25 mA
Current consumption per light section	30 mA
8-pole connector: - Connection capacity - Insulation strip length	0.08 – 1.50 mm ² 6 mm
2-pole connector: - Connection capacity - Insulation strip length	0.33 – 2.00 mm ² 6 mm
Switching current per output RL, AW, WL	100 mA at 24 V DC
Dimensions (HxWxD)	110 x 150 x 40 mm
Degree of protection	IP 20
Ambient temperature	+5 °C – +40 °C
Relative humidity	0 % – 85 %

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© Tunstall GmbH, Orkotten 66, 48291 Teigte, Germany, www.tunstall.de