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**Tunstall**

## Intended use

IP-SystemManager 2 (76210100) and IP-SystemManager (76210000) are decentralized control units for the nurse call systems of the Flamenco system family and Flamenco SECURE. All control units of a nurse call system are connected to each other via an IP network

IP-SystemManagers (76210000) have been installed in nurse call systems since 2015. IP-SystemManager 2 (76210100) is the second generation of this control unit. IP-SystemManagers 2 (76210100) have been installed in nurse call systems since 2025. The two models cannot be operated together in one nurse call system. Please note the compatibility information in the next section.

Table 1. Products

Article no.	Product name
76210100	IP-SystemManager 2
76210000	IP-SystemManager

## Notes on compatibility

### IP-SystemManager 2 (76210100)

The following control units can be operated with IP-SystemManager 2 (76210100) in one nurse call system:

- IP-SystemManager 2 (76210100)
- IP-SystemManagers (76210000), if these have previously been updated to IP-SystemManager 2 firmware by Tunstall via a paid update.

Incompatibility: A connection to an OSY-ControlCenter or OSY-ControlCenter<sup>SE</sup> is not possible.

### IP-SystemManager (76210000)

The following control units can be operated with IP-SystemManager (76210000) in one nurse call system:

- IP-SystemManager (76210000)
- OSY-ControlCenter or OSY-ControlCenter<sup>SE</sup>

Upgrade option: An IP-SystemManager (762100 00) can be upgraded to the IP-SystemManager 2 (76210100) firmware for a fee so that it can be operated together in one nurse call system. However, it will then no longer be compatible with the OSY-ControlCenter or OSY-ControlCenter<sup>SE</sup>.

## Safety notes

### NOTICE

The complete installation of the system is described in the Technical Manual.

### CAUTION

The printed circuit board behind the service cover includes electrostatic sensitive components. Avoid touching.

## Product description


All IP-SystemManagers 2 (76210100) are equipped with identical hardware. All IP-SystemManagers (76210000) are equipped with identical hardware.

The individual IP-SystemManager and IP-SystemManager 2 devices differ in terms of their operating mode settings. These determine which functions are active and which connections are enabled.

Devices for the operating modes 'System', 'System + Local' and 'Local' are delivered unconfigured. The required operating mode for the devices must be set during configuration using the SystemOrganizer software. The operating mode, location and IP address must be entered on the device label during configuration.


The operating modes 'BMA', 'Voice Gateway' and 'Voice Gateway Master' can only be set by Tunstall. Devices for these operating modes are therefore preconfigured at the factory and must be installed exactly at the intended location.

Figure 1. Product views

- [1] Housing
  - [2] Service cover
  - [3] a. DIN rail
  - [4] a. 2 DIN rail clamps
- a Not included, order separately

Figure 2. Service cover

The service cover is removed by levering off. After use reattach the service cover.

- [1] Housing
- [2] Service cover
- [5] 6 grommets for fixing the service cover
- [6] Li battery, 3 V, CR2032
- [7] Reset button
- [8] DIP switch for controller ID

## Setting controller ID

Each IP-SystemManager or IP-SystemManager 2 is assigned a controller ID, which is used to identify this device in the SystemOrganizer software. This means that the controller ID specified for this device in SystemOrganizer must be set on this device.

The 6-pole DIP switch behind the service cover is used to set the controller ID, see Figure 2: Coding switches 1 to 6 with the values 1, 2, 4, 8, 16, 32. Set the switch to ON. Examples: Controller ID 1 is set by switching on coding switch 1, controller ID 4 by switching on coding switch 3. Controller ID 24 is formed by addition: set switches 4 and 5 to ON ( $8+16=24$ ).

If the controller ID is changed during operation, a reset must be performed afterwards by pressing the reset button below the DIP switch. Otherwise, the new setting will not be applied.

## Mounting


Installation in dry rooms only. Always easily accessible to authorised persons (access space at least 60 cm wide).

The device must be protected against unintentional contact, e.g. by mounting it in a 19" wall mounting cabinet, 7HU, GD (article no. 77351907).

In case of installation in switch cabinets or similar, dissipated heat shall be discharged by means of forced ventilation, if necessary.


Mounting on DIN rail according to fig. B. e.g. with 19" mounting set (article no. 76090001).

Figure 3. Mounting on a DIN rail

## Dismantling

Figure 4. Dismantling from the DIN rail

## Connecting

- Strip the connection cables to a suitable length. Remove the cable shield and drain wire up to the cable sheath. Strip off the wires 7 mm.

### CAUTION

If the IP-SystemManager/IP-SystemManager 2 is operated without insertion bridges at the power connections, high operating currents from the system flow through the IP-SystemManager/IP-SystemManager 2. This may damage the device. The two insertion bridges must be inserted according to Figure 6!

- Connect the wires to the connectors (pluggable screwtype connectors) according to Figure 5.
- Insert the patch cable's RJ45 plug into the Ethernet socket.

Figure 5. Connection diagram



Figure 6. Bridges

- [a] Insertion bridge  
[b] Wire bridge

## Connections

<p><b>Power:</b>            Connection to the 24 V power supply ring of the ward or directly to a 24 V power supply unit.            Insert insertion bridges!            Conductor cross-section: 2.5 mm<sup>2</sup>            Connector, 4-pole: spare part article no. 00021133.            Two insertion bridges: spare part article no. 00022052.</p>
<p><b>UPS Input:</b>            Input for the signalling contact of a connected uninterruptible power supply (UPS), which signals buffer operation. This signalling contact is available at:</p> <ul style="list-style-type: none"> <li>• UPS control unit 10A (77341100): <b>Buffering</b> contact.</li> <li>• DC -PS with capacitor storage 24V, 6kWS (77341300): <b>Buffering</b> contact.</li> </ul> <p>Conductor cross-section: 0.14 mm<sup>2</sup> – 1.5 mm<sup>2</sup>            Connector, 2-pole: spare part article no. 00021147.</p>
<p><b>Fault Relay:</b>            Only enabled for devices with "Local", "System" or "System + Local" operating mode.            The signalling of the output is in parallel with the <b>Fault LED</b> (below <b>Fault Relay</b>), for details refer to Section : "LED displays" (page 7).            Fault message relay output.            Change over contact, potential free.            Maximum contact load: 2 A, 24 V            Conductor cross-section: 0.14 mm<sup>2</sup> – 1.5 mm<sup>2</sup>            Connector, 3-pole: spare part article no. 00021145.</p>
<p><b>OSYnet 1 and NF 1:</b>            Only enabled for devices with "Local" or "System + Local" operating mode.            Connection of a group bus OSYnet. For the maximum number of bus users that can be connected, please refer to the Planning manual.  <b>CAUTION!</b> IP-SystemManager or IP-SystemManager 2 must not be installed in the middle of the group bus line. It must be connected at one end of the group bus line, as it includes an integrated bus terminating resistor.            In nurse call systems without speech communication, connection <b>NF 1</b> is not used.            For permitted cable types, refer to the Technical manual.            Connector, 3-pole, for OSYnet: spare part article no. 00021145            Connector, 4-pole, for NF: Spare part article no. 00021136.</p>
<p><b>OSYnet 2 and NF 2:</b>            Only enabled for devices with "Local" or "System + Local" operating mode.            As for <b>OSYnet 1</b> and <b>NF 1</b>, but for a second group bus line.            The users on both group buses must belong to the same ward; see the planning manual.</p>
<p><b>Ethernet:</b>            Connection to the Ethernet of the nurse call system. This connects the IP-SystemManager/IP-SystemManager 2 to the other IP-SystemManagers/SystemManagers 2 as well as the IP-based devices of the nurse call system.  <b>NOTE!</b> The Ethernet port is internally equipped with a separation point according to EN 60601-1 (2 x MOPP). An external separation device is not required.            Ethernet, 10/100 Mbit LAN connection, RJ45 socket.            Maximum cable length: 90 m. Patch cable, at least cat. 5e.</p>
<p><b>RS232 A:</b>            ESPA 4.4.4 protocol.            For example, connection of a radio paging system, DECT system or fire alarm system.  <b>NOTES!</b> If room terminals are connected to the OSYnet port and an RS232 port is used a separator device according to EN 60601-1 (2 x MOPP) must be connected between the connected device and the RS232 port, e.g. RS-232 isolator 2xMOPP, 4kV (article no. 76500010).            Maximum cable length: 10 m            Conductor cross-section: 0.14 mm<sup>2</sup> – 1.5 mm<sup>2</sup>            Connector, 5-pole, spare part order no.: 00021137.</p>
<p><b>RS232 B:</b>            As <b>RS232 A</b>, but for a second system.</p>

## LED displays

<p><b>Green LED 'Power':</b>            LED is lit continuously. = This device is powered.</p>
<p><b>Green LED 'Status':</b>            LED flashes once a second = This device is ok.            LED is permanently lit or is off. = Fault! This device is not operational.</p>
<p><b>Red LED 'Fault' below 'Fault Relay':</b>            LED is lit continuously. = There is at least one fault.            Signalling in parallel to <b>Fault Relay</b> output.            This is a collective display. The assigned fault locations depend on the operating mode of this device:</p> <p><b>Operating mode 'System' or 'System + Local':</b>            Collective indicator for faults in the entire nurse call system. Indicated types of fault:</p> <ul style="list-style-type: none"> <li>• An OSYnet bus user is defective or disconnected.</li> <li>• A RAN bus user is defective or disconnected.</li> <li>• The power supply on a ward is interrupted. This fault is recognized indirectly by a failure in OSYnet bus user polling.</li> <li>• The power supply at an IP-SystemManager/IP-SystemManager 2 is interrupted. This fault is recognized indirectly by a failure in OSYnet bus user polling.</li> <li>• OSYnet bus line is interrupted.</li> <li>• LAN connection to an IP-SystemManager/IP-SystemManager 2 is interrupted.</li> <li>• An IP-SystemManager/IP-SystemManager 2 is not operational.</li> <li>• The power supply unit is in buffer mode. Indication in parallel to the red LED <b>Fault</b> below <b>UPS Input</b>.</li> </ul>
<p><b>Operating mode 'Local':</b>            Collective indicator for faults in the ward connected to this device. Indicated fault types:</p> <ul style="list-style-type: none"> <li>• An OSYnet bus user is defective or disconnected.</li> <li>• A RAN bus user is defective or disconnected.</li> <li>• The power supply unit is in buffer mode. Indication in parallel to the red LED <b>Fault</b> below <b>UPS Input</b>.</li> </ul>
<p> <b>NOTICE</b>            The red LED <b>Fault</b> (below <b>Fault Relay</b>) and the <b>Fault Relay</b> connection of IP-SystemManagers/IP-SystemManagers 2 with operating mode "BMA", "Voice Gateway" and "Voice Gateway Master" are inactive.</p>
<p><b>Red LED 'Fault' below 'UPS Input':</b>            LED is lit continuously. = The uninterruptible power supply is operating in buffer mode.</p>
<p><b>'OSYnet 1': Green LED 'Rx', yellow LED 'Tx':</b>            LEDs are flashing. = Data traffic on the group bus, connected to <b>OSYnet 1</b>.            Yellow LED is lit continuously. = Group bus OSYnet 1 has not yet been configured with the SystemOrganizer.</p>
<p><b>'OSYnet 2': Green LED 'Rx', yellow LED 'Tx':</b>            As <b>OSYnet 1</b>, but for <b>OSYnet 2</b>.</p>
<p><b>Yellow LED 'NF1':</b>            LED is lit continuously. = There is a speech connection or an announcement in the group connected to <b>OSYnet 1</b>.</p>
<p><b>Yellow LED 'NF2':</b>            As <b>NF1</b>, but for <b>OSYnet 2</b>.</p>
<p><b>'Ethernet', green LED 'Link', yellow LED 'Speed':</b>            The LEDs are signalling "Link" and "Speed" of the Ethernet connection.</p>
<p><b>'RS232 A', green LED 'Rx', yellow LED 'Tx':</b>            LEDs are flashing. = Data traffic between this device and the device connected to <b>RS232 A</b>.</p>
<p><b>'RS232 B', green LED 'Rx', yellow LED 'Tx':</b>            As <b>RS232 A</b>, but for <b>RS232 B</b>.</p>

## Replace battery

Behind the service cover (see Figure 2) is a lithium button cell (3 V, CR2032) for buffering the internal clock in the event of a power failure.

In normal operation of the nurse call system, i.e. the device is supplied with power, the battery should be replaced every 2 years. If the device is not supplied with power, e.g. when the device is stored, an earlier replacement may be necessary due to the battery consumption.

Replacement battery: high quality Li battery, 3 V, CR2032

### Behaviour of the system when the battery is empty

An empty battery leads to the following situation: After a failure of the supply voltage of the IP-SystemManager/IP-SystemManager 2, the internal clock restarts at 00:00. An incorrect time is displayed on the display devices of the nurse call system (ComStation, corridor display, ComTerminals, RoomTerminals, CareTerminals). The time is correct again when a connected TimeServer sends the time. If no TimeServer is connected, the time must be reset manually via the SystemOrganizer.

### Battery disposal

The used battery must be returned to a collection point for used batteries or to Tunstall.

## Technical data

Dimensions (HxWxD)	132 x 216 x 48 mm
Housing material	PS UL94-HB
Service cover material	ABS
Weight	415 g
IP protection	IP20
Nominal voltage	24 V DC
Quiescent current consumption	120 mA
Permitted voltage range	20 – 26 V DC
Ambient temperature	+5 °C to +40 °C
Relative humidity	0 % to 85 % (non-condensing)

Tunstall GmbH, Orkotten 66, 48291 Telgte, Germany  
Tel: +49 (0) 2504/701-0, E-Mail: DE.info@tunstall.com  
[www.tunstall.de/en-de/](http://www.tunstall.de/en-de/)

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Tunstall GmbH, Orkotten 66, 48291 Telgte, Allemagne