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Flamenco Systemfamilie

Rufanlagen Flamenco^{IP+}, Flamenco^{IP}, Flamenco

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The diagram illustrates the network architecture of the HMI system. At the bottom, the **LAN** is connected to the **Ethernet** port of the HMI rack. The rack contains several modules: a **PoE+** module, a **PoE+** module, a **PoE+** module, a **PoE+** module, a **PoE+** module, a **PoE+** module, a **PoE+** module, and a **PoE+** module. These modules are connected to various components:

- Top Section:**
 - 76 0510 00** (HMI) connected to the first **PoE+** module.
 - 76 0510 00** (HMI) connected to the second **PoE+** module.
 - 76 0510 00** (HMI) connected to the third **PoE+** module.
 - 76 0510 00** (HMI) connected to the fourth **PoE+** module.
 - 76 0510 00** (HMI) connected to the fifth **PoE+** module.
 - 76 0510 00** (HMI) connected to the sixth **PoE+** module.
 - 76 0510 00** (HMI) connected to the seventh **PoE+** module.
 - 76 0510 00** (HMI) connected to the eighth **PoE+** module.
- Middle Section:**
 - 76 0510 00** (HMI) connected to the first **PoE+** module.
 - 76 0510 00** (HMI) connected to the second **PoE+** module.
 - 76 0510 00** (HMI) connected to the third **PoE+** module.
 - 76 0510 00** (HMI) connected to the fourth **PoE+** module.
 - 76 0510 00** (HMI) connected to the fifth **PoE+** module.
 - 76 0510 00** (HMI) connected to the sixth **PoE+** module.
 - 76 0510 00** (HMI) connected to the seventh **PoE+** module.
 - 76 0510 00** (HMI) connected to the eighth **PoE+** module.
- Bottom Section:**
 - 76 0510 00** (HMI) connected to the first **PoE+** module.
 - 76 0510 00** (HMI) connected to the second **PoE+** module.
 - 76 0510 00** (HMI) connected to the third **PoE+** module.
 - 76 0510 00** (HMI) connected to the fourth **PoE+** module.
 - 76 0510 00** (HMI) connected to the fifth **PoE+** module.
 - 76 0510 00** (HMI) connected to the sixth **PoE+** module.
 - 76 0510 00** (HMI) connected to the seventh **PoE+** module.
 - 76 0510 00** (HMI) connected to the eighth **PoE+** module.

Additional components and connections include:

- 76 0510 00** (HMI) connected to the first **PoE+** module.
- 76 0510 00** (HMI) connected to the second **PoE+** module.
- 76 0510 00** (HMI) connected to the third **PoE+** module.
- 76 0510 00** (HMI) connected to the fourth **PoE+** module.
- 76 0510 00** (HMI) connected to the fifth **PoE+** module.
- 76 0510 00** (HMI) connected to the sixth **PoE+** module.
- 76 0510 00** (HMI) connected to the seventh **PoE+** module.
- 76 0510 00** (HMI) connected to the eighth **PoE+** module.
- 76 0510 00** (HMI) connected to the first **PoE+** module.
- 76 0510 00** (HMI) connected to the second **PoE+** module.
- 76 0510 00** (HMI) connected to the third **PoE+** module.
- 76 0510 00** (HMI) connected to the fourth **PoE+** module.
- 76 0510 00** (HMI) connected to the fifth **PoE+** module.
- 76 0510 00** (HMI) connected to the sixth **PoE+** module.
- 76 0510 00** (HMI) connected to the seventh **PoE+** module.
- 76 0510 00** (HMI) connected to the eighth **PoE+** module.
- 76 0510 00** (HMI) connected to the first **PoE+** module.
- 76 0510 00** (HMI) connected to the second **PoE+** module.
- 76 0510 00** (HMI) connected to the third **PoE+** module.
- 76 0510 00** (HMI) connected to the fourth **PoE+** module.
- 76 0510 00** (HMI) connected to the fifth **PoE+** module.
- 76 0510 00** (HMI) connected to the sixth **PoE+** module.
- 76 0510 00** (HMI) connected to the seventh **PoE+** module.
- 76 0510 00** (HMI) connected to the eighth **PoE+** module.

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3. Wall mounting (→ fig. C)

Wall mounting using the keyhole openings on the back of the unit.

- Unobstructed view of the display must be ensured at a distance of up to 20 m.
- A distance of at least 50 cm to the wall must be maintained on the left and right of the corridor display^{IP}, because loudspeakers [2] are built in both side housing covers.
- The readability of the display must not be impaired by external light sources.
- Note the position of the ambient brightness sensor [4]! The sensor is used to automatically adjust the brightness of the display to the ambient brightness.
- Mounting height above floor level = 150 – 250 cm.

4. Dismantling (→ fig. C)

1. Disconnect the connection plugs.
2. Remove the corridor display^{IP} from the wall.

5. Connecting (→ fig. D)

1. Plug the RJ45 plug of the patch cable into the RJ45 socket of the corridor display^{IP}.
2. If power is not supplied by PoE+: Connect the corridor display^{IP} to the 24 V DC power supply ring of the ward or directly to a 24 V DC power supply unit via the 2-pole connector.

IP network for the nurse call system	
Connection socket	RJ45 socket
Cable type for PoE+	Min. CAT5e, shielded, Ø min. 0.64 mm (22 AWG)
Cable type without PoE+	Min. CAT5e, shielded
Max. cable length	90 m

24 V DC power supply	
Connector	Pluggable screw terminal, 2-pole; Spare part: 70 0807 08.
Screw	M3
Cable type	NYM 2x2.5 mm ²
Stripping length	8 mm
Connection capacity	0.08 – 4.00 mm ² , max. 1 wire

6. Start-up of the unit

When the corridor display^{IP} has been connected to the power supply, it starts up. Sequence of indications in the display: Hardware revision > Firmware revision > All LEDs light up. > Setting* of the function “call tone” (Call vol.) > Setting* of the function “announcements” (Call vol.) > Setting* of the function “time display” (Time) > MAC address.

*: “auto” = Enabled. “off” = Disabled.

At the end of the start-up, the display shows one of the following indications:

- **Time.** Meaning: Corridor display^{IP} is in operation.
- **LED at the bottom left of the display is flashing.** Meaning: Corridor display^{IP} is in operation. “Time display” is disabled.
- **MAC address of the unit.** Meaning: ERROR! Corridor display^{IP} is not operational. Corridor display^{IP} is not connected to the IP-SystemManager control unit. Correct the error!

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7. Functions that can be disabled (option)

The following main functions are active in the factory settings of the corridor display^{IP}, but can be disabled. Disabling the functions is done in the system menu, see chapter 8.

Call tone (Call vol.)

The “call tone” function, i.e. the acoustic signalling of calls shown in the display may only be disabled, if it is ensured that staff on duty are made aware of calls by other signalling.

Announcements (Ann. vol.)

The “announcements” function, i.e. the transmission of announcements may only be disabled, if it is ensured that the reception of the announcements is not required in the vicinity of the corridor display^{IP}.

Time display (Time)

If the display of the time is disabled, an LED flashes at the bottom left of the display in standby mode to indicate that the corridor display^{IP} is in operation.

8. Using the system menu (option)

The system menu serves to:

- display the revisions and IP address of the unit
- test the loudspeakers and the LEDs in the display.
- disable the functions that can be disabled

NOTE! Before using the system menu, be sure to read chapter “7. Functions that can be disabled”.

This is how you use the System menu:

1. Remove the protective cap from the programming button [5], see fig. D.
2. Press and hold the programming button for 10 seconds, e.g. with a small screwdriver or pen.
 - ✓ The display shows: “System-Menu”.
3. Press the programming button briefly.
 - ✓ Display shows hardware revision, e.g. “Rev. A2”.
4. Press the programming button briefly.
 - ✓ Display shows firmware revision, e.g. “FW R1.02”.
5. Press the programming button briefly.
 - ✓ The display shows the own IP address, e.g. “IP: 172.40.2.5”.
6. Press the programming button briefly.
 - ✓ Loudspeakers are tested by successively emitting a soft, a medium and a loud tone.
7. Press the programming button briefly.
 - ✓ All LEDs are tested by lighting up.
8. Press the programming button briefly.
 - ✓ Display shows the setting of the “call tone” (Call vol.) function: “auto” = Enabled. “off” = Disabled.
9. Set the desired option: Press the programming button briefly to switch between “auto” and “off”.

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Continued from previous page:

10. Press and hold the programming button for approx. 1.5 seconds.
✓ Display shows the setting of the “announcements” (Ann. vol.) function : “auto” = Enabled. “off” = Disabled.
11. Set the desired option: Press the programming button briefly to switch between “auto” and “off”.
12. Press and hold the programming button approx. 1.5 seconds.
✓ Display shows the setting of the “Time display” (Time) function: “auto” = Enabled. “off” = Disabled.
13. Set the desired option: Press the programming button briefly to switch between “auto” and “off”.
14. To close the system menu and save the settings, press and hold the programming button for 3 seconds.
✓ The system menu is closed. The settings are saved.
15. Replace the protective cap of the programming button [5].
NOTE! You can close the system menu at any menu item and save the settings by pressing and holding the programming button for 3 seconds.

Perform function tests

If you have used the System menu, you must check that all three functions that can be disabled are set as desired: Call tone, Transmission of announcements, time display.

9. Technical data

Standby current consumption	130 mA at 24 V DC
Current consumption on call	440 mA at 24 V DC
Standby power consumption	3.8 W at PoE+
Power consumption on call	12.0 W at PoE+
Dimensions (HxWxD)	125 x 785 x 55 mm
Character height	50 mm
Weight	approx. 1800 g
Housing material	Aluminium, lacquered, with end caps made of ABS, lacquered
Front window material	Acrylic glass
Degree of protection	IP20
Ambient temperature	+5°C – +40°C
Relative humidity	0 % – 85 % (non-condensing)

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Corridor display^{IP} Alpha 16, double-sided, order no. 76 0160 00

1. Intended use

The corridor display^{IP} is designed for connection to the IP network of a Flamenco^{IP+} nurse call system. Power is supplied either via PoE+ or via the 24 V DC power supply of the nurse call system. The corridor display^{IP} is used for displaying calls and for the transmission of announcements.

NOTE! The complete installation of the system is described in the Technical Manual.

2. Product description (→ fig. C)

The corridor display^{IP} shows the call type and room number of the call location alphanumerically. The display has 16 digits; longer texts are displayed as a ticker.

The calls are also signalled acoustically. The three call categories (calls, emergency calls, alarms) differ in their tone sequence.


The acoustic call signalling as well as announcements are transmitted via two integrated loudspeakers.

In standby mode, the time is displayed (factory setting).

The corridor display^{IP} is configured with the SystemOrganizer software. Three functions can be disabled on the corridor display^{IP} itself.

3. Mounting location requirements

- Unobstructed view of the display must be ensured at a distance of up to 20 m.
- A distance of at least 50 cm to the wall must be maintained on the left and right of the corridor display^{IP}, because loudspeakers [2] are built in both side housing covers.
- The readability of the display must not be impaired by external light sources.
- Note the position of the ambient brightness sensor [3]! The sensor is used to automatically adjust the brightness of the display to the ambient brightness. A Sensor [3] is only present on the side of the corridor display^{IP} shown in fig. A.
- Mounting height above floor level = 150 – 250 cm.

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- [1] Connection sockets
[2] Two loudspeakers
[3] Ambient brightness sensor

- [4] Programming button (2 mm under protective cap)
[12] Wire rope holder

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4. Ceiling suspension (→ fig. C – F)


1. Drill two holes in the ceiling (distance: 525 mm). Insert the dowels [5].
2. Each wire rope: Bolt together the wire rope [10] with fitting [8], ceiling fastener [7], screw cap [9] and hanger bolt [6] according to fig. E.
3. Each wire rope: Slide the wire rope [10] into the wire rope holder [12] up to the desired position. To lock the wire rope [10], pull it. A gripping mechanism is activated. Fig. F.

NOTE! To release an inadvertently activated locking, press the head element [11] of the wire rope holder [12]. The gripping mechanism will be released.

4. Screw the two hanger bolts [6] of the pre-mounted corridor display^{IP} into the dowels [5] in the ceiling.

Optional accessories

Extension set ceiling suspension, order no. 19 0780 05, to extend the ceiling suspension by 50 cm.

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- | | |
|-------------------------------|---|
| [1] Connection sockets | [8] Fitting |
| [2] Two loudspeakers | [9] Screw cap |
| [3] Ambient brightness sensor | [10] Wire rope (approx. 500 mm) |
| [5] Dowel (Ø 8 mm) | [11] Head element of the wire rope holder |
| [6] Hanger bolt | [12] Wire rope holder |
| [7] Ceiling fastener | |

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5. Dismantling (→ fig. F)

1. Disconnect the connection plugs.
2. A second person must hold the corridor display^{IP}.
3. Release the lockings of the two ropes one after the other by pressing down the respective head element [11].

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6. Connecting (→ fig. B)

1. Plug the RJ45 plug of the patch cable into the RJ45 socket of the corridor display^{IP}.
2. If power is not supplied by PoE+: Connect the corridor display^{IP} to the 24 V DC power supply ring of the ward or directly to a 24 V DC power supply unit via the 2-pole connector.

IP network for the nurse call system	
Connection socket	RJ45 socket
Cable type for PoE+	Min. CAT5e, shielded, Ø min. 0.64 mm (22 AWG)
Cable type without PoE+	Min. CAT5e, shielded
Max. cable length	90 m

24 VDC power supply	
Connector	Pluggable screw terminal, 2-pole; Spare part: 70 0807 08.
Screw	M3
Cable type	NYM 2x2.5 mm ²
Stripping length	8 mm
Connection capacity	0.08 – 4.00 mm ² , max. 1 wire

7. Start-up of the unit

When the corridor display^{IP} has been connected to the power supply, it starts up. Sequence of indications in the display: Hardware revision > Firmware revision > All LEDs light up. > Setting* of the function “call tone” (Call vol.) > Setting* of the function “announcements” (Call vol.) > Setting* of the function “time display” (Time) > MAC address.

*: “auto” = Enabled. “off” = Disabled.

At the end of the start-up, the display shows one of the following indications.

- **Time.** Meaning: Corridor display^{IP} is in operation.
- **LED at the bottom left of the display is flashing.** Meaning: Corridor display^{IP} is in operation. “Time display” is disabled.
- **MAC address of the unit.** Meaning: ERROR! Corridor display^{IP} is not operational. Corridor display^{IP} is not connected to the IP-SystemManager control unit. Correct the error!

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8. Functions that can be disabled (option)

The following main functions are active in the factory settings of the corridor display^{IP}, but can be disabled. Disabling the functions is done in the system menu, see chapter 9.

Call tone (Call vol.)

The “call tone” function, i.e. the acoustic signalling of calls shown in the display may only be disabled, if it is ensured that staff on duty are made aware of calls by other signalling.

Announcements (Ann. vol.)

The “announcements” function, i.e. the transmission of announcements may only be disabled, if it is ensured that the reception of the announcements is not required in the vicinity of the corridor display^{IP}.

Time display (Time)

If the display of the time is disabled, an LED flashes at the bottom left of the display in standby mode to indicate that the corridor display^{IP} is in operation.

9. Using the system menu (option)

The system menu serves to:

- display the revisions and IP address of the unit
- test the loudspeakers and the LEDs in the display.
- disable the functions that can be disabled

NOTE! Before using the system menu, be sure to read chapter “8. Functions that can be disabled”.

This is how you use the System menu:

1. Remove the protective cap from the programming button [4], see fig. D.
2. Press and hold the programming button for 10 seconds, e.g. with a small screwdriver or pen.
 - ✓ The display shows: “System-Menu”.
3. Press the programming button briefly.
 - ✓ Display shows hardware revision, e.g. “Rev. A2”.
4. Press the programming button briefly.
 - ✓ Display shows firmware revision, e.g. “FW R1.02”.
5. Press the programming button briefly.
 - ✓ The display shows the own IP address, e.g. “IP: 172.40.2.5”.
6. Press the programming button briefly.
 - ✓ Loudspeakers are tested by successively emitting a soft, a medium and a loud tone.
7. Press the programming button briefly.
 - ✓ All LEDs are tested by lighting up.
8. Press the programming button briefly.
 - ✓ Display shows the setting of the “call tone” (Call vol.) function: “auto” = Enabled. “off” = Disabled.
9. Set the desired option: Press the programming button briefly to switch between “auto” and “off”.

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Continued from previous page:

10. Press and hold the programming button for approx. 1.5 seconds.
✓ Display shows the setting of the “announcements” (Ann. vol.) function : “auto” = Enabled. “off” = Disabled.
11. Set the desired option: Press the programming button briefly to switch between “auto” and “off”.
12. Press and hold the programming button approx. 1.5 seconds.
✓ Display shows the setting of the “Time display” (Time) function: “auto” = Enabled. “off” = Disabled.
13. Set the desired option: Press the programming button briefly to switch between “auto” and “off”.
14. To close the system menu and save the settings, press and hold the programming button for 3 seconds.
✓ The system menu is closed. The settings are saved.
15. Replace the protective cap of the programming button [4].
NOTE! You can close the system menu at any menu item and save the settings by pressing and holding the programming button for 3 seconds.

Perform function tests

If you have used the System menu, you must check that all three functions that can be disabled are set as desired: Call tone, Transmission of announcements, time display.

10. Technical data

Standby current consumption	140 mA at 24 V DC
Current consumption on call	450 mA at 24 V DC
Standby power consumption	4.3 W at PoE+
Power consumption on call	13.0 W at PoE+
Dimensions (HxWxD)	125 x 785 x 55 mm
Character height	50 mm
Weight	approx. 2500 g
Housing material	Aluminium, lacquered, with end caps made of ABS, lacquered
Front windows material	Acrylic glass
Degree of protection	IP20
Ambient temperature	+5°C – +40°C
Relative humidity	0 % – 85 % (non-condensing)

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Group lamp, 3 groups, order no. 77 0113 02

Group lamp, 4 groups, order no. 77 0114 02

Signal lamps for several groups for group related display of calls: 2 groups, 3 groups or 4 groups. Two single sections (red and green) per group with connection possibilities for signalling of calls and call status incl. LED modules.

Dimensions (HxWxD) without mounting frame:

2 groups: 160 x 86 x 70 mm

3 groups: 240 x 86 x 70 mm


4 groups: 320 x 86 x 70 mm

Operating current per light section: 75 mA

Note! The complete installation of the system is described in the technical manual.

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
B

Anschlüsse / Connections:


Plus-Pol: Ruf (rot)
positive terminal: call (red)

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
Corridor Display Alpha 16, order no. 77 0150 00

Display for alphanumeric indication of calls and general system information.

- 16 digits, alphanumeric
- Longer texts are automatically displayed as ticker
- Buzzer for signalling forwarded calls
- Loudspeaker, e.g. for announcements or signal tones
- Supply current: 300 mA (average)

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A Setting of address (1 - 110)

For setting the address please use the code switches 1 to 7 with the values 1, 2, 4, 8, 16, 32, 64. Set the appropriate switch to ON.



Examples: Address 1 is set by using the code switch 1. Address 4 is set by using the code switch 3. Address 24 is selected through addition: coding switch 4 and 5 ($8 + 16 = 24$).

After the address is changed while system operation a **reset** has to be made. For the reset set code switch 10 for one second to ON.

B Mounting

The corridor display is designed for mounting in corridors or duty rooms.

The standard accessory is used for ceiling mounting.

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

C Connection




The corridor display is directly connected to the group bus OSYnet. Use the free wire ends of the connection cable for on site connection.

D Functional test

When the corridor display is connected, a functional test should be carried out.

- To start the functional test set code switch 9 "Test" to ON.
- The corridor display will run through the test steps described in the following table. Each step takes approx. 2 seconds.

Display shows:	Test ok, if...	What is being tested?
Start Test	<i>Functional test starts.</i>	
Tunstall GmbH Display Check 0123456789abcdef	The 3 texts are displayed one after the other.	Optical display.
Buzzer test	<i>Buzzer test starts.</i>	
Buzzer: low	Low signal tone sounds.	Buzzer setting "low".
Buzzer: medium	Medium signal tone sounds.	Buzzer setting "medium".
Buzzer: loud	Loud signal tone sounds.	Buzzer setting "loud".
Buzzer OFF	<i>Buzzer is switched off.</i>	
Speech relay 1 ON	Low relay click can be heard.	Speech channel 1 of system.
Speech relay 2 ON	Low relay click can be heard.	Speech channel 2 of system.
Speech relays OFF	Relay click can be heard.	Speech channels are closed.
 (displayed as ticker)	Display without "holes". Holes show defective LEDs.	LED function.
To stop the test: DIP switch 9 -> OFF	<i>To stop the test set code switch 9 to OFF.</i>	
End of Test	<i>Test ended.</i>	
00:00 (time)	<i>If OSYnet is connected: System time is displayed. Otherwise internal clock of corridor display starts at 00:00 h.</i>	

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
Corridor Display Alpha 16, double-sided, order no. 77 0160 00

Double-sided display for indication of calls and general system information.

- 16 digits, alphanumeric
- Longer texts are automatically displayed as ticker
- Tone generator for signalling forwarded calls
- Loudspeaker, e.g. for announcements or signal tones
- Supply current: 600 mA (average)

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

A Mounting

The corridor display is designed for mounting in corridors or duty rooms.

The standard accessory is used for ceiling mounting.

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C Setting of address (1 - 110)

Behind the protection cover at the side of one display there is a 10 part code switch. For setting the address please use the code switches 1 to 7 with the values 1, 2, 4, 8, 16, 32, 64. Set the appropriate switch to ON.

Examples: Address 1 is set by using the code switch 1. Address 4 is set by using the code switch 3. Address 24 is selected through addition: coding switch 4 and 5 ($8+16=24$).

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D Functional test


After the corridor display has been connected, a functional test should be carried out.

- To start the functional test set code switch 9 "Test" to ON.

The corridor display will run through the test steps described in the following table. Each step takes approx. 2 seconds. You have to carry out the test twice, i.e. check each display separately.

Display shows:	Test ok, if...	What is being tested?
Start Test	<i>Functional test starts.</i>	
Tunstall GmbH Display Check 0123456789abcdef	The 3 texts are displayed one after the other.	Optical display.
Buzzer test	<i>Buzzer test starts.</i>	
Buzzer: low	Low signal tone sounds.	Buzzer setting "low".
Buzzer: medium	Medium signal tone sounds.	Buzzer setting "medium".
Buzzer: loud	Loud signal tone sounds.	Buzzer setting "loud".
Buzzer OFF	<i>Buzzer is switched off.</i>	
Speech relay 1 ON	Low relay click can be heard. ¹⁾	Speech channel 1 of system.
Speech relay 2 ON	Low relay click can be heard. ¹⁾	Speech channel 2 of system.
Speech relays OFF	Relay click can be heard. ¹⁾	Speech channels are closed.
LEDs (displayed as ticker)	Display without "holes". Holes show defective LEDs.	LED function.
To stop the test: DIP switch 9 -> OFF	To stop the test set code switch 9 to OFF.	
End of Test	<i>Test ended.</i>	
00:00 (time)	If OSYnet is connected: System time is displayed. Otherwise internal clock of corridor display starts at 00:00 h.	

¹⁾ The relays for the speech channels are only in the display with the coding switch.

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
Red control LED flashes slowly (1 sec ON, 1 sec OFF), if connection to OSYnet is ok. Otherwise there is a fault.

77 0160 00
Flurdisplay Alpha 16, doppelseitig
Corridor Display Alpha 16, double-sided

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EN - Installation Instructions

Room lamp, 3 sections, order no. 77 0170 00

Signalling of all call types and staff presences.

Room lamp, 3 sections, glass decor order no. 77 0175 00

Same as 77 0170 00, but decorative glass frame.

Room lamp cardiac alarm, WC, order no. 77 0170 01

Complement to room lamp 77 0170 00 for explicit signalling of cardiac alarm and/or WC calls. Mounting either in the corridor next to 77 0170 00 or next to the WC door in the patient room.

Room lamp cardiac alarm, WC, glass decor, 77 0175 01

Same as 77 0170 01, but as complement to 77 0175 00, decorative glass frame.

Room lamp, 4 sections, order no. 77 0170 10

Signalling of all call types and staff presences as well as additional display for WC call.

Room lamp, 4 sections, glass decor, order no. 77 0175 10

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A Connection

1. Strip the connecting cable in the back box to a suitable length.
2. Remove the cable shield and drain wire up to the cable sheath. Strip the wires.
3. Connect the wires to the connector [3] according to fig. A.

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B DIP switch, 4 pole


DIP switch 1 – 3: Channel coding

If functional units are requested, the room lamps have to be assigned to different channels. For the channel number to be set please refer to the installation plan. Set the channel number according to fig. B.

DIP switch 4: Signalling

OFF = Fresh calls of the call category “Calls” are signalled with a short attention blink. This makes it possible to distinguish between fresh calls and answered calls. (OFF = factory setting)

ON = Fresh calls of call category “Calls” are signalled as continuous light according to the German standard DIN VDE 0834-1.

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LED module colours

Lamp:	77 0170 00 77 0175 00	77 0170 01 77 0175 01	77 0170 10 77 0175 10
	red: calls	blue: cardiac alarm	red: calls
	green: staff 1	blue: cardiac alarm	green: staff 1
	yellow: staff 2	white: WC call	yellow: staff 2
	no LED module	no LED module	white: WC call

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C Mounting

At the point of delivery the room lamp is assembled and must be dismantled as follows, see fig. D:

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
Wall installation on a one-gang back box [1] fixed with screws on the lateral holes of the back box *):

1. Prepare the wires and connect them to the connector [3] (see section "Connection").
2. Plug the connector [3] from the rear onto the PCB [5] (in the rear of the housing [2]).
3. Bolt the housing [2] to the back box [1] with the screws [4]. The arrow on the rear of the housing points toward the top.
4. Press the light dome with insert module [8] onto the decorative frame until it latches in audibly.

*) **NOTE!** To mount the room lamp on the vertical holes of the back box, the PCB must be removed temporarily:

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NOTE! The room lamp, 4 sections, equals to 1 RAN user to be set at the room terminal. The previous model equalled to 2 RAN users. Distinctive mark of the previous model: It has no DIP switch [7].

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- [1] * Back box
[2] Housing (glass decorative frame looks different)
[3] * Connector (70 0807 00)
[4] * Back box screws
[5] Printed circuit board (PCB) with LED modules

- [6] Fastening screw for the PCB
[7] DIP switch
[8] Light dome with insert module
* Not included in the scope of delivery, please order separately.

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EN - Installation Instructions

Room lamp, 3 sections, with doorplate, order no. 77 0171 00

Signalling of all call types and staff presences. Doorplate as label field for room designation.

Room lamp, 4 sections, with doorplate, order no. 77 0171 10

Same as 77 0171 00, but with additional display for WC call.

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

LED module colours

Lamp:	77 0171 00	77 0171 10
	red: calls	red: calls
	green: staff 1	green: staff 1
	yellow: staff 2	yellow: staff 2
	no LED module	white: WC call

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
1. Prepare the wires and connect them to the connector [3] (see section "Connection").
2. Plug the connector [3] from the rear onto the PCB [5] (in the rear of the housing [2]).
3. Bolt the housing [2] to the back box [1] with the screws [4]. The arrow on the rear of the housing points toward the top.
4. Press the light dome with insert module [8] onto the decorative frame until it latches in audibly.
5. Insert the label strip [9] and the protection cover [10].

***) NOTE!** To mount the room lamp on the vertical holes of the back box, the PCB must be removed temporarily:

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- [1] *Back box
[2] Housing
[3] *Connector (70 0807 00)
[4] *Back box screws
[5] Printed circuit board (PCB) with LED modules
[6] Fastening screw for the PCB

- [7] DIP switch
[8] Light dome with insert module
[9] *Label strip (Height x Width: 70 x 92 mm)
[10] Protection cover

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**Room lamp, 4 sections, C/A/PR/WC**

- Order no. 77 0176 00: silver design frame
- Order no. 77 0176 01: white design frame

Intended use

Room signal lamp for connection to the room bus RAN of a nurse call system of the Flamenco system family. Four light sections for visual signalling of calls and presence. Installation on walls in corridors and rooms.

LED modules	Colour	Signalling (C/A/PR/WC)	Spare part order no.
	Red	C: Call	77 0190 00
	Blue	A: Cardiac alarm	77 0190 04
	Green	PR: Presence	77 0190 02
	White	WC: WC call	77 0190 03

Signalling

Event	Colour	Fresh call
Call	Red	Standard ¹ : Flashing Option ¹ : Steady light
WC call	Red + White	Standard ¹ : Flashing Option ¹ : Steady light
Emergency call	Red	Rapid flashing
WC emergency call	Red + White	Rapid flashing
Cardiac alarm	Blue	Flickering
Presence	Green	Steady light

¹) The option can be set by DIP switch.

The following applies to all types of call: By answering a fresh call, the flashing changes to a steady light. This means that answered calls light are lit steadily.

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
CAUTION! The printed circuit board and the LED modules include electrostatic sensitive components. Avoid touching.

Defining the mounting position (→ fig. B)

- Mounting on a 1-gang back box (for solid wall: 17 0100 00, for partition wall: 17 5100 00).
- The room lamp must be easily assignable to the room.
- External sources of light must not prevent the room lamp from being identified.
- Mounting height above floor level: 150 – 250 cm.

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Connecting the connector (→ fig. C)

1. Strip the connecting cable in the back box to a suitable length.
2. Remove the cable shield and drain wire up to the cable sheath. Strip the wires. Stripping length: 7 mm.
3. Connect the wires to the connector [3] according to fig. C.

CAUTION! The unused wires must not be connected through.

DIP switch, 4-pole (→ fig. D)


DIP switch 1 – 3: Channel coding

Room lamps can be assigned to different channels to form functional units. For the channel number to be set, see the relevant installation plan. Set the channel number as shown in fig. D.

DIP switch 4: Signalling

The section “Signalling” describes a standard (factory setting) and an option. You set these with DIP switch 4:

- DIP switch 4 OFF = Standard (factory setting)
- DIP switch 4 ON = Option

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Note on the ComTerminal in the room

The ComTerminal in the room has presence buttons for staff 1 (green) and staff 2 (yellow). However, the room lamp only displays presence messages of staff 1. Therefore, the presence button for staff 2 must be disabled via the configuration menu of the ComTerminal.

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Mounting (→ fig. E)

At the point of delivery the room signal lamp is assembled and must be dismantled as follows, see fig. F:



- Lightly compress the light dome [8] from top and bottom. Then pull off the light dome [8].

Wall installation on a one-gang back box [1] fixed with screws on the lateral holes of the back box²:

1. Prepare the wires and connect them to the connector [3] (see section “Connecting the connector”).
2. Plug the connector [3] from the rear onto the PCB [5] (in the rear of the housing [2]).
3. Bolt the housing [2] to the back box [1] with the screws [4]. The arrow on the rear of the housing points toward the top.
4. Press the light dome with insert module [8] onto the design frame until it latches in audibly.

²) **NOTE!** To mount the room lamp on the vertical holes of the back box, the PCB must be removed temporarily:

1. Loosen the fastening screw [6] of the PCB and remove the PCB [5] (first bottom, then top).
2. The PCB [5] is installed in the reverse order.

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- | | |
|--------------------------------------|--|
| [1] * Back box, 1-gang | [5] Printed circuit board (PCB) with LED modules |
| [2] Housing | [6] Fastening screw for the PCB |
| [3] * Connector, 3-pole (70 0807 00) | [7] DIP switch |
| [4] * Back box screws | [8] Light dome |

* Not included in the scope of delivery, please order separately.

Technical data

Standby current consumption	20 mA
Additional current consumption	30 mA per light section
Standby power consumption	0.48 W
Additional power consumption	0.72 W per light section
Dimensions (HxWxD)	110 x 150 x 40 mm
Housing material	ABS
Light dome material	Styrene-acrylonitrile
Connecting capacity	up to 1.5 mm ²

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77 0180 10: Room lamp universal, 3 sections

77 0185 10: Room lamp universal, 3 sections, glass decor

for signalling of all call types (red), presence of staff 1 (green) as well as
additional display for WC call (white).

77 0180 00: Room lamp universal, 4 sections

77 0185 00: Room lamp universal, 4 sections, glass decor

as 77 0180 10, but additionally presence of staff 2 (yellow).

77 0182 10: Room lamp universal, 2 sections

77 0185 00: Room lamp universal, 2 sections, glass decor

for signalling presence of staff 1 (green) and telephone call (white).

77 0182 50: Room lamp universal, 1 section

77 0185 50: Room lamp universal, 1 section, glass decor


for signalling of all call types (red).

NOTE! The complete installation of the system is described in
the Technical Manual.

CAUTION! The LED modules include electrostatic sensitive com-
ponents. Avoid touching.

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- [1] *Back box, 1-gang
(solid wall: 17 0100 00,
partition wall: 17 5100 00)
- [2] *Connector, 7-pole
(70 0807 07)
- [3] Housing

- [4] *Back box screws
 - [5] PCB with LED modules
 - [6] Fastening screw for the PCB
 - [7] Light dome with insert module
- * Not included in the scope of deliv-
ery, please order separately.

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B Mounting

At the point of delivery the room lamp is assembled and must be dis-
mantled as follows, see fig. D:

- Lightly press the light dome [7] together from top and bottom. Then
pull off the light dome [7].

Wall installation on a one-gang back box [1] fixed with screws on the lat-
eral holes of the back box, see fig. B:

1. Prepare the wires and connect them to the connector [2] (see sec-
tion "Connection").
2. Plug the connector [2] from the rear onto the PCB [5] (in the rear of
the housing [3]).
3. Bolt the housing [3] to the back box [1] with the screws [4]. The ar-
row on the rear of the housing points toward the top.
4. Press the light dome with insert module [7] onto the decorative
frame.

To mount the PCB on the vertical holes of the back box, the PCB must
be removed temporarily:

1. Loosen the fastening screw [6] of the PCB and remove the PCB (first
bottom, then top).
2. The PCB is installed in the reverse order.

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







Connection

1. Strip the connection cable in the back box to a suitable length. Remove the cable shield and drain wire up to the cable sheath. Strip the wires.
2. Connect the wires to the connector [2] according to the connection diagram.
3. If required: Plug the LED module onto the PCB (on the rear of the housing) according to section „LED module colours“.



Note: The colour of an LED module can be identified by a colour spot of the same colour on the surface of the LED module.

LED module colours

	Connection point	77 0180 10 77 0185 10	77 0180 00 77 0185 00
	L1	Red: Calls	Red: Calls
	L2	Green: Staff 1	Green: Staff 1
	L3	No LED module	Yellow: Staff 2
	L4	White: WC call	White: WC call

	Connection point	77 0182 10 77 0185 00	77 0182 50 77 0185 50
	L1	No LED module	Red: Calls
	L2	Green: Staff 1	No LED module
	L3	No: LED module	No LED module
	L4	White: Phone call	No LED module

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Pull cord switch insert

Intended for connection to the room bus (RAN) of a Tunstall nurse call system.

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
A Product description

The pull cord switch serves for triggering calls or WC calls. A pull cord [7] and a red call button [9] serve for triggering calls. An LED in the call button [9] lights up faintly to locate the button in the dark (location light). The LED lights up brightly as soon as a call is triggered (reassurance light). If presence in the room is switched on when the call is triggered, an emergency call is triggered.

A safety release in the pull cord [7] opens when the pulling force exceeds a limit value of approx. 25 Nm. The safety release is easily closed again by simply pushing it together.

As an option: WC call cancel function

WC calls are cancelled in the standard case by pressing a grey cancel button at the call location. The pull-cord switch can be configured in such a way that WC calls can not only be triggered with the pull cord [7] but also cancelled if presence is switched on.

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- [1] *Back box, 1-gang
(solid wall: 17 0100 00,
partition wall: 17 5100 00)
- [2] Mounting plate
- [3] *Back box screws
- [4] *Cover frame, 1-gang

- [5] Terminal block
(plug-in terminals)
- [6] Pull cord switch insert
- [7] Pull cord (2.50 m) incl. safety release
- [8] Call handle
- [9] Call button (red)

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B Mounting

CAUTION! Pull cord switches in shower stalls must be installed at least 20 cm above the highest possible position of the shower head.

Ceiling or wall installation on one-gang back box:

1. Bolt the mounting plate [2] to the back box [1]. The four square recesses in the mounting plate are located at the top and bottom.
2. Prepare the wires and feed them through the mounting plate (see section "Connection").
3. Connect the pull cord switch insert [6] according to section "Connection".
4. Hold the cover frame [4] in front of the mounting plate on the wall.
5. Plug the pull cord switch insert [6] onto the mounting plate so that the latching hook of the pull cord switch insert fits into the square recesses of the mounting plate. The call button is located at the top.
6. The call handle [8] must be located between 10 and 20 cm above the floor. The safety release must be located above the position of the pull cord where it will be pulled to trigger a call. To maintain these characteristics of the pull cord [7], shorten the pull cord [7] at the safety release or in the call handle and re-knot it there again.

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EN - Installation Instructions

C Dismantling

1. Pull the cover frame [4] together with the pull cord switch insert [6] off the mounting plate [2].
2. For disconnecting the wires, refer to fig. E.



D Connection

1. Pull the connection cable out through the mounting plate and strip it to a suitable length. Remove the cable shield and drain wire up to the cable sheath. Strip the wires.

2. Connect the wires according to the connection diagram.

NOTE: The terminal block [5] can be pulled off temporarily.

3. If the pull cord switch insert is to be used for triggering WC calls, a wire bridge must be fitted between connection points LT1 and LT2. No wire bridge is to be set for triggering standard calls (= factory setting).

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Configuration

The following functions can be configured if required:

Channel number (1 – 6)

A channel number can be set on the pull cord switch to specify the call location in the room in more detail. The channel number can, for example, identify beds or WC cubicles.

- Factory setting = channel 0, i.e. no channel set.

WC call cancel function

The WC call cancel function is available from production week 24 in 2021. The production week can be found in the serial number on the unit: XXXXXXXX-XX-**2421**-XXXXX.

If the pull cord switch is set to "WC call" (= wire bridge between LT1 and LT2), you can configure the pull cord switch so that WC calls (and WC emergency calls) can not only be triggered but also cancelled with the pull cord.

CAUTION! The channel number also applies to the WC call cancel function. If a channel number is set, WC calls that were triggered with this pull cord switch or another call switch with the same channel number can be cancelled using this pull cord switch. WC calls from all channel numbers can be cancelled with a pull cord switch set to channel 0 (no channel).

NOTE! To comply with the German standard DIN VDE 0834, an additional cancel switch must be installed.

- Factory setting = WC call cancel function is disabled.

How to proceed to the configuration

CAUTION! When configuring, both the channel number (0 – 6) and the WC call cancel function (on or off) are always configured. It is not possible to configure only one parameter.

For configuring the pull cord switch within the first ten seconds after the power supply is switched on, start with the following procedure. (Alternatively, the terminal block [5] can be pulled off briefly and then re-plugged.)

1. Press and hold the button [9] for 5 seconds until a triple tone sounds. (During the 5 seconds, the LED must be off).

The LED in the button flashes according to the channel number set: Channel 1 = once, Channel 2 = twice and so on, Channel 0 (no channel) = LED lights up for 2 seconds.

2. Set the channel number with the button as follows: Channel 1 = press 1 time briefly, Channel 2 = press 2 times briefly, etc. To set channel 0 (no channel), do not press the button.
3. Configure the WC call cancel function with the pull cord as follows: Enable = pull 1 time (short tone is emitted). Disable = do not pull the pull cord.
4. To save both settings, press and hold the button again for 5 seconds until a tone sounds.

The LED in the button flashes according to the channel number set: Channel 1 = once, Channel 2 = twice and so on, Channel 0 (no channel) = LED lights up for 2 seconds.

Timeout: If the button is not pressed within 30 seconds, a tone sounds and the call pull switch exits configuration mode. The new settings are not saved.

Function test: After configuration, be sure to perform a function test of the configured functions.

Function test of the WC call cancel function

With the WC call cancel function, WC calls and WC emergency calls can be cancelled when presence is switched on:

1. Trigger a WC call with the pull cord switch.
2. Press a presence button in the room.
3. Hold the pull cord pulled for 5 seconds until a triple tone sounds and then release.

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Presentation of the switches

...*) The letter at the end of the order number represents the frame of the switch:

A = Frame dimensions (HxW): 91 x 91 mm.

F = Frame dimensions (HxW): 80 x 80 mm.

C = Frame dimensions (HxW): 107 x 107 mm.

Caution! The frame **C** is made from real glass!

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
Call switch, order no. 77 0211 00...*)


Two call buttons (red) for raising calls. The LEDs (red) are on with a weak light for finding the buttons in the dark (night location light). The LEDs are on with a bright light, when a call was raised (reassurance light).

Optional setting: WC calls are raised instead of calls.

Call switch/WC, order no. 77 0211 01...*)

same as 77 0211 00..., but for WC call.

Optional setting: Calls are raised instead of WC calls.

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

Staff presence switch, order no. 77 0212 00...*)

Two presence buttons for switching on and off the staff presence for two different staff groups. The LEDs are on, when the respective presence is switched on (reminder light).

- Presence button (green) with LED (green) for staff 1.
- Presence button (yellow) with LED (yellow) for staff 2.


“Call tone” function available, but disabled in the factory setting; for details refer to page 4.

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
Cancel switch/WC, order no. 77 0213 00...*)

Two cancel buttons (grey) for cancelling WC calls and WC emergency calls. The LEDs (yellow) are on, when a WC call or a WC emergency is raised that can be cancelled with these buttons (reminder light).

“Call tone” function available, but disabled in the factory setting; for details refer to page 4.

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

EN - Installation Instructions

Cardiac alarm switch, order no. 77 0214 00...*)

Two alarm buttons (blue) for raising cardiac alarms, if staff presence is switched on in the room. When the alarm switch is active, the LEDs (red) are on with a weak light for finding the buttons in the dark (night location light). The LEDs are on with a bright light, when an alarm was raised (reassurance light).

Optional setting: Cardiac alarms can also be raised, if no staff presence is switched on in the room.

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

Call switch/WC with cancel switch, 77 0217 00...*)

Call button (red) for raising WC calls. The LED (red) is on with a weak light for finding the button in the dark (night location light). The LED is on with a bright light, when a WC call was raised (reassurance light).

Cancel button (grey) for cancelling WC calls and WC emergency calls. The LED (yellow) is on, when a WC call or a WC emergency is raised that can be cancelled with this button (reminder light).

“Call tone” function available, but disabled in the factory setting; for details refer to page 4.

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
Call switch with privacy switch, 77 0218 00...*)

Call button (red) for raising calls. The LED (red) is on with a weak light for finding the button in the dark (night location light). The LED is on with a bright light, when a call was raised (reassurance light).

Privacy button (grey) for switching on and off the privacy in the room, if staff presence is switched on in the room. The LED (white) is on, if the privacy was switched on with the privacy button.

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

Staff presence combination with call tone, 77 0219 00...*)

Call button (red) for raising calls. The LED (red) is on with a weak light for finding the button in the dark (night location light). The LED is on with a bright light, when a call was raised (reassurance light).

Presence button (green) for switching on and off the staff presence of staff group 1. The LED (green) is on, when presence 1 is switched on (reminder light).

“Call tone” function enabled: The integrated beeper signals forwarded calls, if staff presence is switched on in the room. The beeper also attracts attention for announcements.

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A Mounting the pattress

Connect the 3-pole connector [4] to the room bus RAN (without speech) according to Fig. A.

Risk of malfunctions in the nurse call system!

- The cable shield and the shield wire of the RAN room bus (without speech) are not connected and must be removed in the back box up to the cable sheath.
- The unused wires must not be interconnected.


B Mounting

Wall mounting on 1-gang back box (see Fig. B):

1. Screw the support ring [7] to the back box [8] with the back box screws [5]. The openings [6] in the support ring [7] for the locking hooks [2] of the switch module [1] must be located on the top and on the bottom.
2. Lead the connected connector [4] through the support ring [7] out of the back box [8].
3. Place the switch module [1] in the frame [3]. The LEDs of the switch module [1] must be located on the left side.
4. Plug the connected connector [4] to the socket on the rear side of the switch module [1].
5. Press the switch module [1] together with the frame [3] onto the support ring [7] until it engages.

Dismantling

1. Lever the switch module [1] together with the frame [3] off the support ring [7].
2. Remove the connector [4] from the switch module [1].
3. Undo the screws [5] of the support ring [7].
4. Remove the support ring [7].

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Enabling the “Call tone” function

The „call tone” function is available for the following switches:

- Staff presence switch (77 0212 00...)
- Cancel switch/WC (77 0213 00...)
- Call switch/WC with cancel switch (77 0217 00...)
- Staff presence combination with call tone (77 0219 00...)

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This function is only enabled in the factory setting of the staff presence combination with call tone. For all other switches it is disabled.

CAUTION! Calls from other rooms must be signalled with a call forwarding tone in rooms with presence activated. If the room terminal is not able to sound call forwarding tones (e.g. ControlTerminal) a switch must sound the call forwarding tones (e.g. staff presence combination with call tone).

This is how you enable the “call tone” function

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1. Press and hold one button of the switch for approx. 10 seconds until the upper LED lights up yellow. (Sequence: After approx. 5 seconds the switch emits a beep and the upper LED lights up pink. After a further 5 seconds the switch emits a second beep and the upper LED lights up yellow. Now release the button.)

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- Call tone disabled = Upper LED flashes once in red.
- Call tone enabled = Upper LED flashes once in green.

2. By pressing one button repeatedly the setting switches between Green (for “Call tone enabled”) and Red (for “Call tone disabled”). Press the button until the desired setting.

3. To save the setting press and hold one button for approx. 5 seconds until the switch emits a beep.

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Enabling the optional setting

For the following switches you can enable the optional setting instead of the factory setting, refer to the switch description:

- Call switch (77 0211 00...)
- Call switch/WC (77 0211 01...)
- Cardiac alarm switch (77 0214 00...)

On the next page is described, how you enable the optional setting.

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This is how you enable the optional setting

Within the first 10 seconds after the power supply is switched on, you have to start with the following procedure. (Alternatively, the connector of the switch can be pulled off briefly and then re-plugged.)

1. Press and hold one button of the switch for approx. 10 seconds until the upper LED lights up white. (Sequence: After approx. 5 seconds the switch emits a beep and the upper LED lights up pink. After a further 5 seconds the switch emits a second beep and the upper LED lights up white. Now release the button.)

The upper LED signals the current setting:

- Factory setting = Upper LED flashes once in red.
 - Optional setting = Upper LED flashes once in green.
2. By pressing one button repeatedly the setting switches between Green (for optional setting) and Red (for factory setting). Press the button until the desired setting.
 3. To save the setting press and hold one button for approx. 5 seconds until the switch emits a beep.

The new setting is saved. The upper LED signals the saved setting.

NOTE! If you do not press a button within 30 seconds, the new setting is not saved.

Setting a channel number (0 – 7)

If functional units are requested, the switches have to be assigned to different channels. Refer to your installation plan to find out, whether channel numbers have to be assigned. Factory set to channel 0.

Within the first 10 seconds after the power supply is switched on, you have to start with the following procedure. (Alternatively, the connector of the switch can be pulled off briefly and then re-plugged.)

1. Press and hold one button of the switch for approx. 5 seconds until the switch emits a beep and the upper LED lights up pink.

The upper LED signals the current setting:

- Channel 0 = Upper LED flashes once in red (factory-set).
 - Channel 1 = Upper LED flashes once in green.
 - Channel 2 = Upper LED flashes twice in green.
 - and so on ...
 - Channel 7 = Upper LED flashes seven times in green.
2. Set the channel number predefined in your installation plan as follows:
 - Channel 0 = Do not press a button.
 - Channel 1 = Press one button once.
 - Channel 2 = Press one button twice.
 - and so on ...
 - Channel 7 = Press one button seven times.
 3. To save the setting press and hold one button of the switch until the switch emits a beep.

The new setting is saved. The upper LED signals the saved setting.

NOTE! If you do not press a button within 30 seconds, the new setting is not saved.

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
Pull cord call switch, order no. 77 0215 00...*)

2,50 m pull cord and call handle for raising calls or WC calls (factory setting: calls). The LEDs (red) are on with a weak light for finding the pull cord call switch in the dark (night location light). The LEDs are on with a bright light, when a call was raised (reassurance light).

The safety release opens, if the pulling force exceeds a limit value.

Pull cord call switch/WC, order no. 77 0215 01...*)

same as order no. 77 0215 00..., but factory setting: WC call.

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
R(ed) = +24V-Si

Setting a channel number (0 – 7)

If functional units are requested, the pull cord call switches have to be assigned to different channels. Refer to your installation plan to find out, whether channel numbers have to be assigned. Factory set to channel 0.

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1. Pull and keep pulled the cord for approx. 5 seconds until the switch emits a beep and the upper LED lights up pink.
 - ✓ The upper LED signals the current channel number:
 - Channel 0 = Upper LED flashes once in red (factory setting).
 - Channel 1 = Upper LED flashes once in green.
 - Channel 2 = Upper LED flashes twice in green.
 - and so on ...
 - Channel 7 = Upper LED flashes seven times in green.
2. Set the channel number predefined in your installation plan as follows:
 - Channel 0 = Do not pull the cord.
 - Channel 1 = Pull the cord once.
 - Channel 2 = Pull the cord twice.
 - and so on ...
 - Channel 7 = Pull the cord seven times.
3. To save the setting pull and keep pulled the pull cord until the switch emits a beep.
 - ✓ The new setting is saved.

Note! If you do not pull the cord within 30 seconds, the new setting is not saved.

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1. Pull and keep pulled the cord for approx. 10 seconds until the upper LED lights up white. (Sequence: After approx. 5 seconds the switch emits a beep and the upper LED lights up pink. After a further 5 seconds the switch emits a second beep and the upper LED lights up white. Now release the cord.)
 - ✓ The upper LED signals the current setting for the call type:
 - Factory setting = Upper LED flashes once in red.
 - Optional setting = Upper LED flashes once in green.
2. Set the call type as follows:
 - Factory setting = Do not pull the cord.
 - Optional setting = Pull the cord once.
3. To save the setting pull and keep pulled the pull cord until the switch emits a beep.
 - ✓ The new setting is saved.
 - Note! If you do not pull the cord within 30 seconds, the new setting is not saved.
 - ✓ The upper LED signals the set call type:
 - Factory setting = Upper LED flashes once in red.
 - Optional setting = Upper LED flashes once in green.

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1. Screw the support ring 7 to the back box 8* with the back box screws 5*. The openings 6 in the support ring 7 for the locking hooks 2 of the switch module 1 must be located on the top and on the bottom, see Fig. B.
2. Lead the connected connector 4* through the support ring 7 out of the back box 8*.
3. Place the switch module 1 in the frame 3. The LEDs of the switch module 1 must be located on the left side, see Fig. B.
4. Plug the connected connector 4* to the socket on the rear side of the switch module 1.
5. Press the switch module 1 together with the frame 3 onto the support ring 7 until it engages.
6. The call handle 11 must be located between 10 and 20 cm above the floor. Therefore shorten the pull cord 10 to the desired length and re-knot it in the call handle 11.

Safety release

The safety release 9 opens, if the pulling force exceeds a limit value. To re-establish the connection just plug the two parts of the safety release together.

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
C = Frame dimensions (HxW): 107 x 107 mm. **Caution!** The frame **C** is made from real glass!

Note! The complete installation of the system is described in the Technical Manual.

Pneumatic call switch, order no. 77 0216 00...*)

2 m air hose and red rubber bulb for raising calls or WC calls by squeezing the bulb (factory setting: Calls). The LEDs (red) are on with a weak light for finding the pneumatic call switch in the dark (night location light). The LEDs are on with a bright light, when a call was raised (reassurance light).


Pneumatic call switch/WC, order no. 77 0216 01...*)

same as order no. 77 0216 00..., but factory setting: WC call.

Connecting the connector

Connect a 3-pole connector (order no. 70 0807 00) to the room bus RAN (without speech) according to Fig. A.

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
1. Screw the support ring 7 to the back box 8* with the back box screws 5*. The openings 6 in the support ring 7 for the locking hooks 2 of the switch module 1 must be located on the top and on the bottom, see Fig. B.
2. Lead the connected connector 4* through the support ring 7 out of the back box 8*.
3. Place the switch module 1 in the frame 3. The LEDs of the switch module 1 must be located on the left side, see Fig. B.
4. Plug the connected connector 4* to the socket on the rear side of the switch module 1.
5. Press the switch module 1 together with the frame 3 onto the support ring 7 until it engages.

Warning! Call activation might be prevented. Do not kink the air hose.

Dismantling

1. Lever the switch module 1 together with the frame 3 off the support ring 7.
2. Remove the connector 4* from the switch module 1.
3. Undo the screws 5* of the support ring 7.
4. Remove the support ring 7.

B

Setting a channel number (0–7)

If functional units are requested, the pneumatic call switches have to be assigned to different channels. Refer to your installation plan to find out, whether channel numbers have to be assigned. Factory set to channel 0.

To set a channel number: Within the first ten seconds after the power supply is switched on, you have to start with the procedure described in the following (Alternatively, the connector of the switch can be pulled off briefly and then re-plugged):

1. Squeeze and keep squeezed the rubber bulb for approx. 5 seconds until the switch emits a beep and the upper LED lights up pink.
 - ✓ The upper LED signals the current channel number:
 - Channel 0 = Upper LED flashes once in red.
 - Channel 1 = Upper LED flashes once in green.
 - Channel 2 = Upper LED flashes twice in green.
 - and so on ...
 - Channel 7 = Upper LED flashes seven times in green.
2. Set the channel number predefined in your installation plan as follows:
 - Channel 0 = Do not squeeze the rubber bulb.
 - Channel 1 = Squeeze the rubber bulb once.
 - Channel 2 = Squeeze the rubber bulb twice.
 - and so on ...
 - Channel 7 = Squeeze the rubber bulb seven times.
3. To save the setting squeeze and keep squeezed the rubber bulb until the switch emits a beep.
 - ✓ The new setting is saved.

Note! If you do not squeeze the rubber bulb within 30 seconds, the new setting is not saved.

 - ✓ The upper LED signals the set channel number:
 - Channel 0 = Upper LED flashes once in red.
 - Channel 1 = Upper LED flashes once in green.
 - Channel 2 = Upper LED flashes twice in green.
 - and so on ...
 - Channel 7 = Upper LED flashes seven times in green.

Setting the call type (if needed)

You can set the call type to be initiated by the switch.

- Pneumatic call switch, order no. 77 0216 00...; factory setting: Call. Optional setting: WC call.
- Pneumatic call switch/WC, order no. 77 0216 01...; factory setting: WC call. Optional setting: Call

To set the call type: Within the first ten seconds after the power supply is switched on, you have to start with the procedure described in the following (Alternatively, the connector of the switch can be pulled off briefly and then re-plugged):

1. Squeeze and keep squeezed the rubber bulb for approx. 10 seconds until the upper LED lights up white.
(Sequence: After approx. 5 seconds the switch emits a beep and the upper LED lights up pink. After a further 5 seconds the switch emits a second beep and the upper LED lights up white. Now release the rubber bulb.)
 - ✓ The upper LED signals the current setting for the call type:
 - Factory setting = Upper LED flashes once in red.
 - Optional setting = Upper LED flashes once in green.
2. Set the call type as follows:
 - Factory setting = Do not squeeze the rubber bulb.
 - Optional setting = Squeeze the rubber bulb once.
3. To save the setting squeeze and keep squeezed the rubber bulb until the switch emits a beep.
 - ✓ The new setting is saved.

Note! If you do not squeeze the rubber bulb within 30 seconds, the new setting is not saved.

 - ✓ The upper LED signals the set call type:
 - Factory setting = Upper LED flashes once in red.
 - Optional setting = Upper LED flashes once in green.

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Cancel switch/WC with presence switch, order no. 77 0221 00... *)

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F = Frame dimensions (HxW): 80 x 80 mm

C = Frame dimensions (HxW): 107 x 107 mm.

CAUTION! The frame **C** is made from real glass!

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
“Call tone” function: The integrated beeper signals forwarded calls, if staff presence is switched on in the room. The beeper also attracts attention for announcements.

2. Technical data

Power Supply	24 V DC
Standby current consumption	8 mA
Housing material	Polycarbonate
Housing colour	Studio white, similar to RAL 9016
Ambient temperature	+5°C – +40°C

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- | | |
|--------------------------------------|---|
| [1] Switch module | [6] Openings for the locking hooks |
| [2] Locking hooks | |
| [3] Frame | [7] Support ring |
| [4] * Connector, 3-pole (70 0807 00) | [8] * Back box, 1-gang (solid wall: 17 0100 00, partition wall: 17 5100 00) |
| [5] * Back box screws | |

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- Channel 0 = Upper LED flashes once in red (factory-setting).
 - Channel 1 = Upper LED flashes once in green.
 - Channel 2 = Upper LED flashes twice in green.
 - and so on ...
 - Channel 7 = Upper LED flashes seven times in green.
2. Set the channel number predefined in your installation plan as follows:
 - Channel 0 = Do not press a button.
 - Channel 1 = Press one button once.
 - Channel 2 = Press one button twice.
 - and so on ...
 - Channel 7 = Press one button seven times.

3. To save the setting press and hold one button of the switch until the switch emits a beep.

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7. Disabling the “Call tone” function (if required)

The “call tone” function is enabled on the cancel switch/WC with presence switch: An integrated beeper signals forwarded calls, if staff presence is switched on in the room. The beeper also attracts attention for announcements.

This function can be disabled if required.

WARNING! Calls from other rooms must be signalled with a call forwarding tone in rooms with presence activated. If the room terminal is not able to sound call forwarding tones (e.g. ControlTerminal) a switch must sound the call forwarding tones (e.g. cancel switch/WC with presence switch).

This is how you disable the “call tone” function

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- Call tone enabled = Upper LED flashes once in green.
- Call tone disabled = Upper LED flashes once in red.

2. By pressing one button repeatedly the setting switches between Green (for “Call tone enabled”) and Red (for “Call tone disabled”). Press the button until the desired setting.
3. To save the setting press and hold one button for approx. 5 seconds until the switch emits a beep.

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
- [1] Touch display
- [2] Alarm button (blue)
- [3] Call button (red)
- [4] Presence button (green) for staff 1
- [5] Microphone
- [6] Presence button (yellow) for staff 2
- [7] Loudspeaker

ComTerminal^{IP}, 24V, order no. 76 0510 10ComTerminal^{IP}, PoE, order no. 76 0510 00

1. Intended use

The ComTerminal^{IP} is intended for connection to the IP network of a Flamenco^{IP+} nurse call system. The ComTerminal^{IP} controls and monitors all room functions in compliance with DIN VDE 0834. The room bus RAN is connected to the ComTerminal^{IP}.

ComTerminal^{IP}, 24V (76 0510 10) is connected to the 24 V DC power supply of the ward.

ComTerminal^{IP}, PoE (76 0510 00) is supplied by PoE+.

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
CAUTION! The printed circuit boards include electrostatic sensitive components. Avoid touching.

CAUTION Glass – fragile! The front panel of the ComTerminal^{IP} is made of glass. Handle with care!

2. Back box installation

When installing the back box consider the space required for the ComTerminal^{IP}. Measuring from the centre of the box, the available space must be at least 120 mm above and below.

3. Product description (→ fig. A)

The red call button [3] serves for triggering calls. The call button is faintly lit to locate the button in the dark. It lights up brightly as soon as a call is triggered.

The blue alarm button [2] is only active when presence is switched on. It serves to trigger cardiac alarms. The alarm button is faintly lit to locate the button in the dark. It lights up brightly as soon as a cardiac alarm is triggered.

The green presence button [4] serves for switching the presence of the staff 1 on and off. The yellow presence button [6] serves for switching the presence of the staff 2 on and off. The presence button lights up when the corresponding presence is switched on.

The touch display [1], microphone [5] and two loudspeakers [7] are used to answer forwarded calls from other rooms when presence is switched on.

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4. Mounting the pattress (→ fig. B)

Wall mounting on two-gang back box.*)


Preparation: Ideally, insert the Ethernet cable into the back box from below and use a Keystone module with RJ45 jack for the connection. The RJ45 jack must be connected to the RJ45 socket on the connection field of the ComTerminal^{IP} via the short patch cable [12]. The keystone module is not included in the scope of delivery of the ComTerminal^{IP}.

1. Insert one plug of the short patch cable [12] into the RJ45 jack of the Keystone module. Prepare the wires of the other cables (see section “7. Connecting to the connection field”).
2. Feed the short patch cable [12] and the prepared wires through outlet in the pattress [10].
3. Fix the pattress [10] using the back box screws [11] to the back box [8].
4. Connect the connection wires to the connection field in the pattress [10] according to section “7. Connecting to the connection field”.
5. Insert the second plug of the patch cable [12] into the RJ45 socket on the connection field.

***) NOTE!** If mounting on a back box is not possible, the ComTerminal^{IP} may be fixed to the wall with dowels and screws using the same holes in the pattress.

5. Mounting the front panel (→ fig. C)

1. Place front panel [14] onto the top of the pattress [10] creating a hinge [13].
2. Press the front panel downwards onto the pattress with low pressing power (glass!) until the catch [15] engages.
3. For improved security the security screw [16] can be fitted into the hole below the yellow presence button.

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- | | |
|--|--|
| [7] Two loudspeakers | [12] Short patch cable |
| [8] * Back box | [13] Hinge |
| [9] Two connection wires to the loudspeakers | [14] Front panel |
| [10] Pattress | [15] Catch |
| [11] * Four back box screws | [16] Security screw (M3, 10 x 3 mm, cross-head PZ) |

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6. Dismantling (→ fig. D)

1. If a security screw [16] is fitted into the base of the ComTerminal^{IP} remove it.

CAUTION Glass – fragile! The front panel [14] is fastened only with the catch [15].

2. Use a slotted screwdriver to push the catch [15] upwards until the front panel [14] loosens.
3. Remove the front panel [14].
4. Lay down the front panel [14] protected.

7. Connecting to the connection field (→ fig. E)

24 V power supply (only for 76 0510 10)	
Connector	Screw terminal, 4-pole
Screw	M2.5
Cable type	NYM 2x2.5 mm ²
Wire stripping length	6 mm
Connection capacity	0.05 – 3.31 mm ²


IP network for the nurse call system	
Connector	RJ45 socket
Cable type for 76 0510 00 (PoE+)	Min. CAT5e, shielded, Ø min. 0.64 mm (22 AWG)
Cable type for 76 0510 10: Min. CAT5e, shielded	
The Ethernet port is internally equipped with a separation point according to the German standard DIN VDE 834. An external separator is not required.	

Room bus (RAN)	
Connector	Screw terminal, 8-pole
Screw	M3
Cable type for room bus without speech	J-Y(St)Y 2x2x0.8 or CCS32. Optional: Min. CAT5, shielded, Ø min. 0.57 mm (23 AWG).
Cable type for room bus with speech	2x J-Y(St)Y 2x2x0.8 or CCS32. Options: 2x CAT5/CAT6, shielded, Ø min. 0.57 mm (23 AWG) or 1x CAT5/CAT6 S/FTP, Ø min. 0.57 mm (23 AWG) or 1x CAT7, Ø min. 0.57 mm (23 AWG).
Wire stripping length	6 mm
Connection capacity	0.20 – 2.50 mm ² , i.e. max. 4 wires per connection point

Important instructions for connection of the room bus!



- The cable for speech transmission, i.e. the cable containing the wires “MIC”, “LS”, and “SCH-B”, may be stripped for a max. of 30 mm only, to ensure the protective screen function. The shield wire must be connected to the connection point “SCH-B”.
- The cable shield and the shield wire of the cable containing the wires “RAN”, “+24V-Si” and “0V” are not connected and must be removed in the back box up to the cable sheath.
- The unused wires must not be interconnected.
- Unused screws must be screwed in.

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
- [17] Connection field
 [18] IP network for the nurse
 call system
 [19] Room bus RAN
 [20] 24 V power supply
 (only for 76 0510 10)
 [21] LED lights up green, when
 PoE+ voltage is applied
 (only for 76 0510 00)

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



Kabel: CCS32 (max. 4 Adern pro Klemmpunkt)

**Room bus (RAN) with speech**


Cable: CCS32 (max. 4 wires per connection point)

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8. Using the “Settings” menu

The ComTerminal^{IP} must be configured in the “Settings” menu.

8.1 Opening the “Settings” menu (→ fig. G)

In the standby state, the touch display shows the system time.


1. Triple tap the touch display in the free space.
- ✓ The icon for “Settings” (gear) appears.
2. Tap the gear icon.
- ✓ The “Settings” menu appears.

8.2 Basic operation of the menu

- To scroll through the menu items, drag with a finger down or up.
- To select a menu item, tap the menu item.
- To return to the previous menu level, tap the < icon at the top left of the touch display.

8.3 Closing the “Settings” menu

- To close the “Settings” menu, tap the < icon in the “Settings” menu.

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Menu item	Meaning
Language	Display language of this device.
Room type	Room type, e.g. Patient room.
RAN	Setting of number of room devices. Check whether the stored number of room devices is equal to the number of operational room devices. Test if room devices are ready to operate.
Network	Network configuration of this device.
Contrast	Setting of display contrast.
Key tone	Switching the key tone of this device on or off. Key sound ON is recommended.
Info	Revisions of this device.
Audio test	Only for Tunstall technicians.
Button functions	Function of the buttons of this device.
System	Enable password protection for the “Settings” menu and set the password. Enable Extended mode. Restore factory settings of this device. Restart this device.

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9. Making the mandatory settings

9.1 Setting the user language

In the factory setting, the display language is set to “Automatic”, i.e. the system language set in the nurse call system (SystemOrganizer) is used. If required, you can set a different display language for this device.

1. Tap “Language” in the “Settings” menu.
2. Tap “Automatic” so that it is set to “Off”.
3. Tap the desired language.
4. Return to the main menu by tapping <.

9.2 Selecting the room type

1. Tap “Room Type” in the “Settings” menu.
2. Tap the correct room type:
 - “Patient room”
 - “Staff room” (Calls with presence switched on are not signalled with higher priority. This setting is only recommended in staff rooms of systems with ManagementCenter.)
 - “Children’s room” (Privacy always OFF)
3. Return to the main menu by tapping <.

9.3 Setting of RAN number (0 – 30)

RAN number = Number of devices in the room (switches, connection sockets, room lamps, etc.) that are connected via RAN. (**No** plug-in connected devices, e.g. pear push switches).

1. Tap “RAN” in the “Settings” menu.
2. Tap “RAN number”.
3. Enter the RAN number in the input field.
4. Return to the main menu by tapping <.

9.4 Network

1. Tap “Network” in the “Settings” menu.
- ✓ The “Network” window appears. The current network configuration of the ComTerminal^{IP} is displayed.

For normal operation, the ComTerminal^{IP} must be set to DHCP, i.e. “DHCP” must be highlighted in white.

NOTE! In case of service it may be necessary to set the network configuration manually. To do this, tap “Manual” so that it is highlighted in white and then set the network configuration.

2. Return to the main menu by tapping <.

9.5 Extended mode

The Extended mode must be enabled if the ComTerminal^{IP} is coupled with an optional CareCom system. In all other cases, the extended mode must be disabled.

1. Tap “System” in the “Settings” menu.
2. Change to “On” by tapping on “Off” next to “Extension mode” and vice versa.
“Off” = Extension mode inactive (= factory setting).
“On” = extension mode active.
3. Return to the main menu by tapping <.

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10.1 Button functions

1. Tap “Button functions” in the “Settings” menu.
2. Tap the desired button function.
3. Return to the main menu by tapping <.

10.1.1 Presence 1 button

Green presence button on ComTerminal^{IP} for staff 1:

- “On” = Button is enabled (= factory setting).
- “Off” = Button is without function, although a green circle is visible on the ComTerminal^{IP}. Presence 1 cannot be switched on the ComTerminal^{IP}. To comply with the German standard DIN VDE 0834, a separate presence switch must be installed instead and be easily accessible.

10.1.2 Presence 2 button

Yellow presence button on ComTerminal^{IP} for staff 2:

- “On” = Button is enabled (= factory setting).
- “Off” = Button is without function, although a yellow circle is visible on the ComTerminal^{IP}. Presence 2 cannot be switched on the ComTerminal^{IP}.

10.1.3 Presence cancels WC

- “Off” = WC calls and WC emergency calls must be cancelled with a grey WC cancel button (= factory setting). The factory setting complies with the requirements of the German standard DIN VDE 0834.
- “On” = WC calls and WC emergency calls can be cancelled with a grey WC cancel button or by switching off the presence at the ComTerminal^{IP}. According to the German standard DIN VDE 0834, this setting is only permitted if the call location can be seen from the ComTerminal^{IP}, e.g. in a ward bathroom without cubicles.

10.1.5 Alarm button

NOTE! The alarm button on the ComTerminal^{IP} is only visible and active while presence is switched on.

- “Press briefly” = A cardiac alarm is triggered by briefly pressing the blue alarm button (= factory setting).
- “2 seconds” = A cardiac alarm is triggered when the blue alarm button has been pressed and hold for 2 seconds. This setting deviates from the requirements of the German standard DIN VDE 0834.
- “3 seconds” = A cardiac alarm is triggered when the blue alarm button has been pressed and hold for 3 seconds. This setting deviates from the requirements of the German standard DIN VDE 0834.
- “Off” = The blue alarm button on the ComTerminal^{IP} is disabled. No cardiac alarm can be triggered on the ComTerminal^{IP}.

11. Note on the screen saver

In standby state, the display shows the system time. If no button is pressed for 2 minutes, the display goes dark.

- To reactivate the display, tap any button or the display.

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12. Testing the room bus RAN

12.1 Testing the RAN status

Check whether the set RAN number (refer to section “9.3 Setting of RAN number” equals to the number of operational room devices.

1. Tap “RAN” in the “Settings” menu.
- ✓ The status is displayed under the heading “RAN status”:
 - “OK! Active devices = RAN number” = The set RAN number matches the number of operational devices on the room bus.
 - “ERROR! Active devices > RAN number” = More operational devices recognized on the room bus than have been set.
 - “ERROR! Active devices < RAN number” = Less operational devices recognized on the room bus than have been set.

If an ERROR has been displayed, perform the test as follows:

12.2 Performing a RAN test

Test whether all devices on the RAN have been connected correctly and are ready for operation:

1. Tap “RAN” in the “Settings” menu.
2. Tap on “RAN test”.
- ✓ The RAN test starts. During the RAN test “On” is displayed beside the word “RAN Test”.
 - The LEDs of all devices on the room bus must flash. Devices not flashing have been connected incorrectly or are defective.
 - Replace defective devices.
 - Connect devices correctly that have not been connected correctly.
3. To end the RAN test, tap “RAN test” again.
- ✓ The RAN test is ended. “Off” is displayed beside “RAN Test”.

12.2.1 Inactive devices


In case of the RAN status “ERROR! Active devices < RAN number”, you can display the inactive devices. To do this, tap “RAN status”.

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ComTerminal-E Flamenco, order no. 77 0510 50

Communication terminal with two-way speech communication. Red call button, blue alarm button, green presence button (staff 1), yellow presence button (staff 2) and 4 function keys as soft keys. Microphone and loudspeakers.

Control and monitoring of all room functions in compliance with DIN VDE 0834.

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
Glass – fragile! The operating front cover of the ComTerminal is made of glass. Handle with care!


Back box installation

When installing the back box consider the space required for the ComTerminal. Measuring from the centre of the box, the available space must be at least 120 mm above and below.

A Mounting the pattress

1. Loosen the four back box screws [3].
2. Insert the flat cable [2] through the outlet in the pattress [7] and plug it into the appropriate socket on the connection socket [5].
3. Fix the pattress [7] using the back box screws [3] to the back box [4]. CAUTION! Do not disconnect the two wire connections [6] to the loudspeakers [1].

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- | | |
|--|--|
| 1 - Two loudspeakers | 6 - Two wire connections to the loudspeakers |
| 2 - Flat cable of the pattress | 7 - Pattress |
| 3 - * Four back box screws | 8 - Flat cable of the front panel |
| 4 - * Back box | 9 - Hinge |
| 5 - * Connection socket ComTerminal (order no. 70 0491 00) | 10 - Catch |
| | 11 - Security screw |
| | 12 - Front panel |

* Not included with ComTerminal-E delivery

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
B Mounting the front panel

1. Plug the flat cable [8] of the front panel [12] into the appropriate socket on the pattress [7].
2. Place front panel [12] onto the top of the pattress [7] creating a hinge [9].
3. Press the front panel downwards onto the pattress with low pressing power (glass!) until the catch [10] engages. Make sure that the flat cable [8] is not pinched.
4. For improved security the security screw [11] can be fitted into the hole below the yellow presence button.

Finally, the ComTerminal has to be configured according to the pages "Configuration Instructions".

C Dismantling

1. If a security screw [11] is fitted into the base of the ComTerminal, remove it.
2. Use a screwdriver to push the catch [10] upwards until the front panel [12] loosens.
3. Remove the front panel [12].
4. Disconnect the flat cable [8] from the pattress [7].
5. Lay down the front panel [12] protected. Caution: Electrostatic sensitive components! Caution: Glass!
6. Disconnect the flat cable [2] from the connection socket [5].
7. Loosen the four screws [3] from the connection socket.
8. Demount the pattress [7] from the back box [4].

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- 1 - Two loudspeakers
- 2 - Flat cable of the pattress
- 3 - * Four back box screws
- 4 - * Back box
- 5 - * Connection socket ComTerminal (order no. 70 0491 00)


- 6 - Two wire connections to the loudspeakers
- 7 - Pattress
- 8 - Flat cable of the front panel
- 9 - Hinge
- 10 - Catch
- 11 - Security screw
- 12 - Front panel

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
Symbols in the configuration menu

●	Circle key	Set the marked menu item.
▲	Upwards arrow key	Scroll upwards through the list.
▼	Downwards arrow key	Scroll downwards through the list.
⏏	Home key	Cancel, without storing.

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Menu items

Language:	Selecting the user language, e.g. English.
Room type:	Selecting the room type, e.g. patient room.
RAN > RAN number:	Setting of number of room devices (= RAN number).
RAN > Test RAN:	Test if room devices are ready to operate and are correctly connected to the RAN.
RAN > Status:	Check whether the stored RAN number is equal to the number of operational room devices.
Address:	Setting of room address.
Cleaning time:	Setting the time (seconds) how long it takes to clean (wipe) the front of the ComTerminal.
Contrast:	Function use only for Tunstall technicians.
Key sound:	Switching the key sound of the ComTerminal on or off. "Key sound On" (factory setting) is recommended.
Info:	Displaying the ComTerminal's software revision.
Audio Test:	Only for Tunstall technicians.
Key Function:	Setting the function of the Alarm key and the Presence Keys on the ComTerminal.
Reset:	Function use only for Tunstall technicians. (Restart the ComTerminal).

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Selecting the room type

1. Mark "Room type" using the arrow keys. Then press the circle key.
2. Mark the desired room type using the arrow keys:
 - ♦ Patient room
 - ♦ Patient room with cancel key for WC call at ComTerminal
 - ♦ Staff room (emergency calls are signalled as normal calls)
 - ♦ Children's room (no key to switch the privacy function at ComTerminal, Privacy OFF)
 - ♦ Connection terminal (keys and display at ComTerminal non-operational)
3. Press the circle key to set the selection.

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RAN number = Number of devices in the room (switches, connection sockets, room lamps, etc.) that are connected via RAN. (**No** pear push switch).

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3. Entering of optional settings

The remaining menu items provide optional settings. Set these if required. Examples:

Key Function

Caution! Switched off keys are not available for nursing staff!

Alarm Key

With the factory settings the alarm key on the ComTerminal is operational and is initiated by pressing it.

If there is the risk, that the alarm button is initiated inadvertently, it might make sense to set a delay time (2 or 3 seconds). In that case the alarm will be initiated after the alarm key has been pressed and then hold for 2 or 3 seconds.

Caution! Setting a delay time does not comply with the German standard DIN VDE 0834. A delay time may only be set, where compliance with the DIN VDE 0834 is not required.

In well-founded exceptional cases you can switch of the alarm key function.

Note! Die Alarm Key is only active while the staff presence is switched on in the room.

Presence Key 1 / Presence Key 2

With the factory settings the presence keys on the ComTerminal are operational.

In well-founded exceptional cases you can switch of the key function for these keys.

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4. Checking the room bus RAN

Status

Check whether the set RAN number (refer to section “Setting of RAN number” of chapter “2. Entering of necessary settings”) equals to the number of operational room devices.

1. Mark “RAN” using the arrow keys. Then press the circle key.
2. Mark “Status” using the arrow keys. Then press the circle key.

Display: **OK**: The set RAN number is equal to the number of operational room devices.

Display: **Error**: The set RAN number is not equal to the number of operational room devices.

3. Press the home key to end. In case of an **error** perform a “Test RAN”.

Test RAN

Test if room devices are ready to operate and if they are correctly connected to the RAN.

Check the correct setting of the RAN number.

1. Mark “RAN” using the arrow keys. Then press the circle key.
2. Test start: Mark “Test RAN” using the arrow keys. Then press the circle key.

The LEDs of all connected room devices must flash. Non-flashing devices are wrongly connected or defective. Replace defective devices. Correct any wrong connection.

3. Compare the number of flashing room devices with the set RAN number. If the set RAN number is wrong, correct the RAN number setting after the RAN test as described in section “Setting of RAN number” of chapter “2. Entering of necessary settings”.
4. End of test: Press the home key.

5. Exit the configuration menu

When all settings are made and all tests are completed, you have to exit the configuration menu:

- Press home key several times until normal operation display appears.

NOTE! If for a period of three minutes no function key is pressed, the system will leave the configuration menu.

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
RoomTerminal Flamenco, order no. 77 0520 00

Terminal for use in patient rooms, duty rooms, and function rooms, without speech communication. Red call button, blue alarm button, green presence button (staff 1), yellow presence button (staff 2) and 4 function keys as soft keys.


Control and monitoring of all room functions in compliance with DIN VDE 0834. Connection to the group bus (OSYnet) and the room bus (RAN).

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Back box installation


When installing the back box consider the space required for the RoomTerminal. Measuring from the centre of the box, the available space must be at least 120 mm above and below.

A Mounting the pattress

1. Insert the connection cables through the outlet in the pattress [2].
2. Fix the pattress [2] using the back box screws [1] to the back box [3].

NOTE! If mounting on a back box is not possible, the RoomTerminal may be fixed to the wall with dowels and screws using the same holes in the pattress.

Next, connect the connection cables to the connection field in the pattress [2] according to page "Connections".

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- | | |
|----------------------------|-------------------------------|
| 1 - * Four back box screws | 5 - Flat cable |
| 2 - Pattress | 6 - Hinge |
| 3 - * Back box | 7 - Catch |
| 4 - Front panel | 8 - Security screw (optional) |

* Not included with RoomTerminal delivery

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
D Anschlüsse**B Mounting the front panel**

1. Plug the flat cable [5] of the front panel [4] into the appropriate socket on the connection field in the pattress [2].
2. Place front panel [4] onto the top of the pattress [2] creating a hinge [6].
3. Press the front panel downwards onto the pattress with low pressing power until the catch [7] engages. Make sure that the flat cable [5] is not pinched.
4. For improved security the security screw [8] can be fitted into the hole below the yellow presence button.

Finally, the RoomTerminal has to be configured according to the pages "Configuration Instructions".



C Dismantling

1. If a security screw [8] is fitted into the base of the RoomTerminal, remove it.
2. Use a screwdriver to push the catch [7] upwards until the front panel [4] loosens.
3. Remove the front panel [4].
4. Disconnect the flat cable [5] from the pattress [2].
5. Lay down the front panel [4] protected.
Caution: Electrostatic sensitive components!

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

Connection capacity	0.20 – 2.50 mm ²
Screwdriver blade	0.6 x 3.5 mm
Insulation strip length	6 mm

E**Spannungsversorgung**Kabel: NYM 2x2x2,5 mm²**Power supply**Cable: NYM 2x2x2.5 mm²**Zimmerbus (RAN) ohne Sprechen**

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




* Note! The NF connections are not required in the following installations: There are no ComTerminals (= terminals with speech) connected to this group bus OSYnet AND there shall not be connected any ComTerminals to this group bus OSYnet in the future.

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
1. Start configuration menu

- Press the left and right function keys simultaneously (3 sec.) until the configuration menu is shown in the display.

Symbols in the configuration menu

●	Circle key	Set the marked menu item.
▲	Upwards arrow key	Scroll upwards through the list.
▼	Downwards arrow key	Scroll downwards through the list.
⏠	Home key	Cancel, without storing.

Menu items

Language:	Selecting the user language, e.g. English.
Room type:	Selecting the room type, e.g. patient room.
RAN > RAN number:	Setting of number of room devices (= RAN number).
RAN > Test RAN:	Test if room devices are ready to operate and if they are correctly connected to the RAN.
RAN > Status:	Check whether the stored RAN number is equal to the number of operational room devices.
Address:	Setting of room address.
Contrast:	Function use only for Tunstall technicians.
Key sound:	Switching the key sound of the RoomTerminal on or off. "Key sound On" (factory setting) is recommended.
Info:	Displaying the RoomTerminal's software revision.
Reset:	Function use only for Tunstall technicians. (Restart the RoomTerminal).

2. Entering of necessary settings

Selecting the user language

1. Mark "Language" using the arrow keys. Then press the circle key.
2. Mark the desired language using the arrow keys: D = German; GB = English; etc.
3. Press the circle key to set the selection.

Selecting the room type

1. Mark "Room type" using the arrow keys. Then press the circle key.
2. Mark the desired room type using the arrow keys:
 - ♦ Patient room
 - ♦ Patient room with cancel key for WC call at RoomTerminal
 - ♦ Staff room (emergency calls are signalled as normal calls)
 - ♦ Connection terminal (keys and display at RoomTerminal non-operational)
3. Press the circle key to set the selection.

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1. Mark “RAN” using the arrow keys. Then press the circle key.
2. Mark “RAN number” using the arrow keys. Then press the circle key.
3. Mark the desired RAN number using the arrow keys.
4. Press the circle key to set the selection.

Setting of room address (0 – 110)

1. Mark “Address” using the arrow keys. Then press the circle key.
2. Mark the desired room address using the arrow keys.
3. Press the circle key to set the selection.

3. Entering of optional settings

The remaining menu items provide optional settings. Set these if required.

4. Checking the room bus RAN**Status**

Check whether the set RAN number (refer to “Setting of RAN number”) equals to the number of operational room devices.

1. Mark “RAN” using the arrow keys. Then press the circle key.
2. Mark “Status” using the arrow keys. Then press the circle key.

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ControlTerminal Flamenco, order no. 77 0550 00

Control and monitoring of all room functions in compliance with DIN VDE 0834. Direct connection to the group bus (OSYnet) and the room bus (RAN).

Optical signalling of all call types and staff presences as well as additional display for WC call.

ControlTerminal Flamenco, glass decor order no. 77 0555 00

Same as 77 0550 00, but decorative glass frame.

NOTE! Prior to mounting the ControlTerminal has to be configured using the ControlTerminal ConfigSet, order no. 77 0920 00. The software in the ConfigSet must have revision **1.04** or higher!

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Connections

1. Strip the connecting cables in the back box to a suitable length.
2. Remove the cable shield and drain wire up to the cable sheath. Strip the wires.
3. Connect the wires to the two connectors according to fig. E.

RAN-TEST (DIP switch)

Fig. D: Test if all room devices are ready to operate and if they are correctly connected to the room bus (RAN):

1. Set DIP switch “RAN-Test” to ON (left).
The LEDs of all connected room devices must flash. (Non-flashing devices are wrongly connected or defective. Replace defective devices. Correct any wrong connection.)
2. Compare the number of connected room devices with the RAN number configured with the ControlTerminal ConfigSet. (If the configured RAN number is not equal to the number of connected room devices, you have to correct the RAN number setting with the ControlTerminal ConfigSet.)
3. To end the test, set DIP switch “RAN-TEST” back to OFF (right).

Technical data

Power supply	24 V DC
Standby current consumption	42 mA
Current consumption per light section	30 mA
4-pole connector: - Connection capacity - Insulation strip length	0.20 – 2.50 mm ² 7 mm
8-pole connector: - Connection capacity - Insulation strip length	0.14 – 1.50 mm ² 7 mm
Dimensions (HxWxD)	110 x 150 x 40 mm
Degree of protection	IP 20
Ambient temperature	+5 °C – +40 °C
Relative humidity	0 % – 85 %

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Mounting


At the point of delivery the ControlTerminal is assembled and must be dismantled as follows, see fig. A:

- Lightly compress the light dome [7] from top and bottom. Then pull off the light dome [7].

Wall mounting, see fig. B:

1. Connect the connection wires according to chapter “Connections” of this product leaflet to the connectors (screw clamps) on the mounting plate [2]. Thereby insert two insertion bridges as shown.
2. Unplug the four LED modules [6] and lay them aside.
3. Fit the mounting plate [2] to the back box [1] using the four back box screws [3]. For the correct mounting direction the connector for power supply must be at the top.
4. Place the casing [4] (incl. the pre-mounted PCB) onto the mounting plate [2]. Press very carefully, until plugs and sockets of the connectors are firmly connected.
5. Fit the casing [4] (incl. the pre-mounted PCB) with the two fixing screws [5] to the mounting plate [2].
6. Plug the four LED modules onto the PCB according to the fig. C.
7. Press the light dome with the insert module [7] onto the casing [4] until it locks in place.

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
LED modules

- red: Calls
- green: Staff 1
- yellow: Staff 2
- white: WC call

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
RESET: To reset the ControlTerminal you have to set DIP switch RESET to ON for one second, then back to OFF.

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

- 1 - * Back box
- 2 - ** Mounting plate with connectors
- 3 - * Four back box fixing screws
- 4 - Casing incl. pre-mounted PCB
- 5 - ** Two fixing screws


- 6 - Four LED modules
- 7 - Light dome with insert module

Order separately:


- * Back box (incl. four screws)
- ** Mounting kit, order no. 77 0960 00

E**Spannungsversorgung**Kabel: NYM 2x2x2,5 mm²**Power supply**Cable: NYM 2x2x2.5 mm²**Zimmerbus (RAN) ohne Sprechen**

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




*** NOTE!** The NF wires exist only at the group bus OSYnet, if ComTerminals (= terminals with speech) are connected to the same group bus. These wires are not connected to the ControlTerminal. Each of them has to be connected through using the push-wire connector (order no. 00 0222 88, wire diameter 0.6 – 0.8 mm).

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ControlTerminal with doorplate Flamenco, order no. 77 0551 00

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Optical signalling of all call types and staff presences as well as additional display for WC call. Doorplate as label field for room designation.

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RAN-TEST (DIP switch)

Fig. C: Test if all room devices are ready to operate and if they are correctly connected to the room bus (RAN):

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Technical data

Power supply	24 V DC
Standby current consumption	42 mA
Current consumption per light section	30 mA
4-pole connector: - Connection capacity - Insulation strip length	0.20 – 2.50 mm ² 7 mm
8-pole connector: - Connection capacity - Insulation strip length	0.14 – 1.50 mm ² 7 mm
Dimensions (HxWxD)	190 x 150 x 40 mm
Field for the label strip (HxW)	70 x 92 mm
Degree of protection	IP 20
Ambient temperature	+5 °C – +40 °C
Relative humidity	0 % – 85 %

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Mounting



At the point of delivery the ControlTerminal with doorplate is assembled and must be dismantled as follows, see fig. **A**:

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Wall mounting, see fig. **B**:

1. Connect the connection wires according to chapter “Connections” of this product leaflet to the connectors (screw clamps) on the mounting plate [2]. Thereby insert two insertion bridges as shown.
2. Unplug the four LED modules [6] and lay them aside.
3. Fit the mounting plate [2] to the back box [1] using the four back box screws [3]. For the correct mounting direction the connector for power supply must be at the top.
4. Place the casing [4] (incl. the pre-mounted PCB) onto the mounting plate [2]. Press very carefully, until plugs and sockets of the connectors are firmly connected.
5. Fit the casing [4] (incl. the pre-mounted PCB) with the two fixing screws [5] to the mounting plate [2].
6. Plug the four LED modules onto the PCB according to the fig. **D**.
7. Press the light dome with the insert module [7] onto the casing [4] until it locks in place.

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
Risk of short circuit! When setting the DIP switch do not contact any other electronic components.

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


- 1 - * Back box
- 2 - ** Mounting plate with connectors
- 3 - * Four back box fixing screws
- 4 - Casing incl. pre-mounted PCB
- 5 - ** Two fixing screws
- 6 - Four LED modules


- 7 - Light dome with insert module
 - 8 - * Label strip
 - 9 - Protection cover
- Order separately:**
- * Back box (incl. four screws)
 - ** Mounting kit, order no. 77 0960 00

E**Spannungsversorgung**Kabel: NYM 2x2x2,5 mm²**Power supply**Cable: NYM 2x2x2.5 mm²**Zimmerbus (RAN) ohne Sprechen**

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



Connection socket with call switch, order no. 70 0171 03

Connection socket with two 8pole jacks for connection of call devices like pear push switches or medical monitoring devices. The connection for a third call device (external call switch) is made on the back side with a special connection cable.

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A Mounting

- | | |
|---|--|
| 1* Back box | 4 Cover |
| 2* Connector (70 0807 01) | 5 Fixing clamp (if screw fixing is not possible) |
| 3 Mounting frame with printed circuit board (pre-mounted) | |
- * not included in delivery.

1. Plug the connector 2 into the green socket on the printed circuit board 3.
2. If required: Plug the terminal of the connection cable (50 0308 02) for the external call device into the black socket on the printed circuit board 3.
3. Fix the mounting frame 3 using the screwpoints of the back box 1. If screw fixing is not possible, use the fixing clamp 5 to fix the mounting frame to the box. The fixing clamp cannot be used with back boxes for partition walls.
4. Fit the cover 4 over the mounting frame 3 and lock it.

B Connections

As external call devices only normally-closed contacts may be used. The external call device is connected via the special connection cable, order no. 50 0308 02 (length: 50 cm). The maximum cable length for the external call device is 2.5 m.




C Setting the bed number

The corresponding bed number can be specified by means of the coding switches 1, 2 and 3 of the printed circuit board.

*1 Attention: Diagnostic call only in combination with programmed bed number. With bed number 6 no diagnostic call can be used.


D Dismantling

1. Place the screwdriver in one of the two notches and take off the cover.
2. Remove the screws at the mounting frame 3.
3. Remove the connector 2.

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Connection socket with call switch, bedhead unit, order no. 70 0171 50

Connection socket with two 8-pole sockets for connection of call devices like pear push switches or medical monitoring devices. The connection of an additional external call device incl. location light and reassurance light is made on the back side with a special connection cable. Two outputs for light control. Mounting in a medical supply unit.

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
Electromagnetic compatibility! Due to insufficient suppression at fluorescent lamps in medical supply units faults in the call system may occur. These external interferences possibly can be avoided by installing suppressor elements (varistor circuits). Varistor circuits are available from the manufacturers of the medical supply units. Tunstall offers a proven overvoltage filter 230 V (order no. 70 0890 97).


Medical supply units may differ clearly in EMC behaviour. Even two supply units of the same type may differ, if the installation is different.

When laying call system cables in medical supply units the regulations of European standard ISO 11197 have to be followed.

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



Bed no.:

External call device:

Coding switch 4 must always be set to OFF.

*1 Attention: A diagnostic call can only be used, if a bed number is set. With bed number 6 no diagnostic call can be used.

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

Mounting


The mounting depends on the individual model of the medical supply unit. The connection socket is mounted by the manufacturer of the medical supply unit.

Rückwärtige Montage: / Mounting from backside:

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





As external call device a normally-closed or a normally-open contact may be used. The external call device is connected via the special connection cable, order no. 50 0308 02 (length: 50 cm). The maximum cable length for the external call device is 2.5 m.

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Connection socket with call switch, order no. 70 0171 60...*)

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F = Frame dimensions (HxW): 80 x 80 mm

C = Frame dimensions (HxW): 107 x 107 mm.

CAUTION! Frame C is made from real glass!

Intended use

Connection socket with two 8-pole sockets for connection of call devices like pear push switches or medical monitoring devices. The connection of an additional external call device incl. location light and reassurance light is made on the back side with a special connection cable. Two outputs for light control.

Pressing the call button, triggering a connected pear push switch or triggering the call device connected to the external call input triggers the "Call" call type. An emergency call is triggered when the presence is activated. A medical monitoring device connected to a socket via a diagnostic connection cable (70 0812 10) triggers a diagnostic call.

Connection socket with call switch, WC order no. 70 0171 62...

As 70 0171 60..., except for the following differences:

Pressing the call button, triggering a connected pear push switch or triggering the call device connected to the external call input triggers the "WC call" call type. A WC emergency call is triggered when the presence is activated.

Connecting a medical monitoring device to a socket via a diagnostic connection cable (70 0812 10) is not permitted. The connection socket is not able to trigger a diagnostic call.

Connection socket with call switch, Alarm order no. 70 0171 65...

As 70 0171 60..., except for the following differences:

A call device connected to a socket via a diagnostic connection cable (70 0812 10) triggers the call type "Cardiac alarm". For information on how to connect the connection cable to the call device, please refer to the installation instructions supplied with the connection cable. Please note that this is an exceptional application of the connection cable and that the information on the diagnostic call in the installation instructions is therefore not applicable.




Connecting a medical monitoring device to a socket via a diagnostic connection cable (70 0812 10) is not permitted. The connection socket is not able to trigger a diagnostic call.

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


A Setting coding switches**Bed no.:****External call device:**

*) ATTENTION! A diagnostic call can only be used, if a bed number is set. With bed number 6 no diagnostic call can be used.

Coding switch 4 must always be set to **OFF**.

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Mounting (→ fig. B)

1. Plug the 5-pole connector [2] into the 5-pole socket on the rear side of the PCB [4].
2. If an external call device shall be connected, plug the connector of the connection cable [3] for the external call device into the appropriate socket on the rear side of the PCB [4].
3. Screw the PCB [4] to the back box [1] with the back box screws [5].




CAUTION! Risk of wedging or breaking the red call button! Proceed carefully when pressing the central plate [8] onto the PCB [4].

4. Press the central plate [8] together with the intermediate frame [7] and the frame [6] onto the PCB [4] until it engages audibly on both sides.

Dismantling (→ fig. C)

1. Lever frame [6] together with the intermediate frame [7] and the central plate [8] off the wall with a screw driver.
2. Undo the screws [5] and remove the PCB [4] from the back box [1].
3. Unplug the connectors from the rear of the PCB [4].

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
Technical data

Standby current consumption	21 mA
Central plate colour	Studio white, similar to RAL 9016
Central plate material	Polycarbonate

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Connection socket combi, order no. 70 0424 00, 70 0425 00

Intended use

A connection socket with two different sockets for connection of patient units and call devices. Designed for the following functions: transmission of calls, control of two lights, TV transmission, transmission of entertainment programmes, blinds control. The use of the supported functions depends on the connected patient unit or call device. Suitable for wall mounting.

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Notes on the connections (→ fig. D, E)**Handling of the push wire connectors**

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Shield wire

For connecting the shield wire, use the screw driver to hold the spring contact open, while you are inserting the shield wire. Afterwards check for a firm wire contact by pulling at the wire.

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TV


The available TV functions depend on the TV set in use. Before the final installation, technicians shall check with the ordering party for the project to clarify the details regarding the connection and the functions.

Technical data


Power Supply	24 V DC
Standby current consumption	30 mA
Dimensions (HxWxD)	90 x 181 x 45 mm
Housing colour 70 0424 00	RAL 9018 papyrus white
Housing material 70 0424 00	ABS
Housing colour 70 0425 00	similar to RAL 9016 studio white
Housing material 70 0425 00	Polycarbonate
Degree of protection	IP20
Ambient temperature	+5°C – +40°C
Relative humidity	0 % – 85 %

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









**Connection of light relays dimming function)**

z.B. Digitally dimmable ballast OSRAM QUICKTRONIC Intelligent with DALI interface

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Connection socket combi, TVL order no. 70 0424 50, 70 0425 50

Intended use

A connection socket with two different sockets for connection of patient units and call devices. Designed for the following functions: transmission of calls, control of two lights, TV transmission, transmission of entertainment programmes, blinds control. The use of the supported functions depends on the connected patient unit or call device. TV connection e.g. via IR TV control module universal Audio (70 0360 11). Suitable for wall mounting.

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Mounting (→ fig. B)

1. Fit the socket [1] into the back box.
2. Screw the mounting frame [2] together with the socket [1] onto the back box.
3. Connect both plugs [3] with the sockets (plug A into socket A, plug B into socket B).
4. Push the cover [4] onto the mounting frame [2].

Setting DIP switches (→ fig. C)

Setting of bed numbers 1 to 6: Use DIP switches "P3", "P2", "P1" as shown in the table of figure C.


Example: Bed number 1 = "P3" OFF, "P2" OFF, "P1" ON

* Note: When using bed number 6 no diagnostic call can be used.

Setting DIP switch "TV" for use Blinds control:

- DIP switch "TV" OFF: No blinds control connected.
- DIP switch "TV" ON: Blinds control.

DIP switch "PG" must always be set to OFF.

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Connection socket combi bedhead unit, order no. 70 0434 00, 70 0435 00

Intended use

A connection socket with two different sockets for connection of patient units and call devices. Designed for the following functions: transmission of calls, control of two lights, TV transmission, transmission of entertainment programmes, blinds control. The use of the supported functions depends on the connected patient unit or call device. Suitable for mounting in medical supply units.

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Electromagnetic compatibility

Due to insufficient suppression at fluorescent lamps in medical supply units faults in the nurse call system may occur. These external interferences possibly can be avoided by installing suppressor elements (varistor circuits). Varistor circuits are commercially available from electric/electronic suppliers. Tunstall offers a proven overvoltage filter 230 V (order no. 70 0890 97).

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When laying nurse call system cables in medical supply units the regulations of European standard ISO 11197 have to be followed.

Mounting

After the settings are made the connection socket is mounted into the medical supply unit. The mounting depends on the individual model of the medical supply unit.

1. Open the housing of the medical supply unit.
2. Put in the connection socket.
3. Put on the flat band cable to the medical supply unit.
4. If there is reading light available put on also the connection to the light relay.
5. Close the housing of the medical supply unit.

Technical data

Power Supply	24 V DC
Standby current consumption	30 mA
Dimensions (HxWxD)	81 x 131 x 37 mm
Housing colour 70 0434 00	RAL 9018 papyrus white
Housing material 70 0434 00	ABS
Housing colour 70 0435 00	similar to RAL 9016 studio white
Housing material 70 0435 00	Polycarbonate
Ambient temperature	+5°C – +40°C
Relative humidity	0 % – 85 %

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Setting the connection socket (→ fig. A)

Before mounting the connection socket please adjust as follows:

1. Setting of bed number**Bed no. Jumper placed:**

Bed 1	P1
Bed 2	P2
Bed 3	P1, P2 (= factory setting)
Bed 4	P3
Bed 5	P1, P3
Bed 6 *)	P2, P3 *) When using bed number 6 no diagnostic call can be used.

2. Light option

Jumper open:	Separate outputs LT1 and LT2 for the light relays (light switches are separate).
Jumper placed:	Outputs LT1 and LT2 are connected (light switches are parallel). In this case connect only one relay.

3. Use of blinds control

Jumper 1 placed:	No blinds control connected.
Jumper 2 placed:	Blinds control.


4. TV sound channel TV1 – TV4

For this settings, always place 2 jumpers. Where star-type distribution of TV sound is provided, the sound channel is set in reference to the bed number.

Bed no. Jumper placed:

Bed 1	TV1
Bed 2	TV2
Bed 3	TV3
Bed 4	TV4

Where TV sound is parallel, select sound channel 1, i.e. TV1.

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Notes on the connections (→ fig. B, C)

Room bus (RAN)

- All cables for speech transmission (LS, MIC, SCH B) may be only stripped to a maximum of 30 mm to maintain the shielding function.
- Insulate the drain wire (SCH B) to prevent short circuiting.

Light control


The switching outputs LT1 for reading light and LT2 for room light provide an output voltage of 24 V DC and an output current of max. 200 mA, sourced from the ComTerminal power supply. The potential refers to the 0 V connections of the nurse call system. The outputs prepare a switching pulse (24 V DC, max. 200 mA) for as long as the button on the patient unit is held down (ePat[®]lite, pear push switch, or patient handset).

Wire links in the connection field

- All connection points "+24V-Si" in the connection field must be linked, e.g. by wire links.
- All connection points "0V" in the connection field must be linked, e.g. by wire links.

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
Kabel: 2x J-Y(ST)Y 2x2x0,8 (max. 2 Adern pro Klemmpunkt)

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
Kabel: CCS32 (max. 2 Adern pro Klemmpunkt)



Room bus (RAN) with speech

Cable: 2x J-Y(ST)Y 2x2x0.8 (max. 2 wires per connection point)


Room bus (RAN) with speech

Cable: CCS32 (max. 2 wires per connection point)

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






Kabel: J-Y(St)Y 2x2x0,8

**Connection of light relay for room light (light not dimmable)**

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

**Connection of light relay for room light (incl. dimming function)**

z.B. Digitally dimmable ballast OSRAM QUICKTRONIC
Intelligent with DALI interface

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Connection socket combi bedhead unit, TVL order no. 70 0434 50, 70 0435 50

A connection socket with two different sockets for connection of patient units and call devices. Designed for the following functions: transmission of calls, control of two lights, TV transmission, transmission of entertainment programmes, blinds control. The use of the supported functions depends on the connected patient unit or call device. TV connection e.g. via IR TV control module universal Audio (70 0360 11). Suitable for mounting in medical supply units.

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Electromagnetic compatibility

Due to insufficient suppression at fluorescent lamps in medical supply units faults in the nurse call system may occur. These external interferences possibly can be avoided by installing suppressor elements (varistor circuits). Varistor circuits are commercially available from electric/electronic suppliers. Tunstall offers a proven over-voltage filter 230 V (70 0890 97).

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Technical data

Power Supply	24 V DC
Standby current consumption	30 mA
Dimensions (HxWxD)	81 x 131 x 37 mm
Housing colour 70 0434 50	RAL 9018 papyrus white
Housing material 70 0434 50	ABS
Housing colour 70 0435 50	similar to RAL 9016 studio white
Housing material 70 0435 50	Polycarbonate
Ambient temperature	+5°C – +40°C
Relative humidity	0 % – 85 %

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Connection socket combi bedhead unit universal, order no. 70 0448 00

Combined connection socket with two different sockets for connection of patient handsets and other call devices. Designed for transmission of calls and for routing of control signals for external actuators.

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Suitable for mounting in medical supply units.

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Electromagnetic compatibility! Due to insufficient suppression at fluorescent lamps in medical supply units faults in the call system may occur. These external interferences possibly can be avoided by installing suppressor elements (varistor circuits). Varistor circuits are commercially available from electric/electronic suppliers. Tunstall offers a proven over-voltage filter 230 V (order no. 70 0890 97). Medical supply units may differ clearly in EMC behaviour. Even two supply units of the same type may differ, if the installation is different.

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*) When using bed number 6 no diagnostic call can be used.

2. Light option

Jumper open: Separate outputs LT1 and LT2 for the light control / light relays (light switches are separate).

Jumper placed: Outputs LT1 and LT2 are connected (light switches are parallel). In this case connect only one light control / light relay.

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

1. Open the housing of the medical supply unit.
2. Put in the connection socket.
3. Put on the flat band cable to the medical supply unit.
4. If there is readinglight available put on also the connection to the light control/light relays.
5. Close the housing of the medical supply unit.

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


Connection light control for reading light




24 V signal for the length of time the light switch remains pressed. The output may be loaded with 60 mA. Non-dimmable and dimmable light same as for room light, refer to next page.

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Connections to connection field GB

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Connection socket combi universal, order no. 70 0449 00

Combined connection socket with two different sockets for connection of patient handsets and other call devices. Designed for transmission of calls and for routing of control signals for external actuators. Suitable for wall mounting.

Use only in combination with patient handset universal, order no. 74 0748 00, and RAN interface universal, order no. 70 0848 00.

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Attention! The printed circuit board includes electrostatic sensitive components. Avoid touching.

A Mounting

1. Fit the socket 1 into the back box.
2. Screw the mounting frame 2 together with the socket onto the back box.
3. Connect both plugs 3 with the sockets (A in A, B in B).
4. Push the cover 4 onto the mounting frame 2.

B Connection

- Useable wire gauges: 0.55 - 1 mm Ø
- Shielded wires for LS, MIC: Strip to max. 30 mm!
- Insulate the guide wire (SCH) to prevent short circuiting.
- Connect only one wire per terminal.

Handling of terminals with solderless connection technique:

1. Strip the wire (8 mm).
2. Insert the stripped wire into the terminal until it stops.
3. For detaching the wire insert a screw driver with a tip of 2.5 x 0.4 mm or 3.5 x 0.5 mm into the opening and pull out the wire.

Note! For connecting the shield wire, use the screw driver to hold the spring contact open, while you are inserting the shield wire. Afterwards check for a firm wire contact by pulling at the wire.

C Setting bed number

Setting of bed numbers 1 to 6: Use switches P3, P2, P1 as shown in the table.

Example: Bed number 1 = P3 OFF, P2 OFF, P1 ON

PG and TV switches: Do not change.

*1 Attention: With bed number 6 diagnostic call is not possible.

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Schirmleitung
shield wire


C

	P3	P2	P1
1	off	off	on
2	off	on	off
3	off	on	on
4	on	off	off
5	on	off	on
6	on	on	off

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

Connection socket ComTerminal, order no. 70 0491 00

Connection socket for connecting a ComTerminal in Flamenco system to the group bus (OSYnet) or a ComTerminal in EccoLine or NewLine system or a RAN multiplexer to the ward bus (CCS bus). Suitable for 2-gang back box.

Note! The complete installation of the system is described in the installation manual.

A Placing the back box

Attention! When installing the back box please notice the required space for the ComTerminal. From centre of the back box to the top and to the bottom, a space of at least 150 mm should be allowed.

B Connections

- | | |
|--------------------------|---|
| 1 Socket for ComTerminal | 3 Connection to power supply |
| 2 Fuse 1A/Fast acting | 4 Connections for RAN, group bus and entertainment (Flamenco: no entertainment) |
| order no. 00 0130 24 | |

Important notes

- Usable cable types: Mains: 2,5 mm², other wires: 0,55 - 1 mm Ø
- All cables for music or speech transmission (shielded cables) shall be stripped to a maximum of 30 mm to maintain the shielding function.
- Insulate the shield wire (SCH) to prevent short circuiting.
- Insert only one wire per connection terminal.

C Connection technique

Handling of terminals with solderless connection technique:

- Strip the wire to 8 mm.
- Insert the stripped wire into the terminal until it stops. To make insertion easier (e.g. for CAT5 cable) put a small screw driver into the opening (see diagram) while inserting the wire.
- For detaching the wire put a small screw driver into the opening and pull out the wire.

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
Connections GB

Power	Power supply Cable: NYM 2x2.5 mm ² +24 V from previous room 0 V or to next room	
	Room bus (RAN) with speech (max. 2 connections) Cable: 2x IY(ST)Y 2x2x0.8 	
	Room bus (RAN) with speech (max. 2 connections) Cable: CCS32-Kabel 	
	Room bus (RAN) without speech (max. 4 connections) Cable: IY(ST)Y 2x2x0.8 	
Ward bus / Group bus	Group bus (OSYnet) (Flamenco only) Cable: CAT5 cable 	
	Group bus (OSYnet) (Flamenco only) Cable: IY(ST)Y 4x2x0.8 	
	Ward bus CCS (EccoLine and NewLine only) Cable: CCS 32 cable 	
Entertainm.	Entertainment, selected channel (EccoLine and NewLine only) Cable: 2xIYY per channel or similar cables ELA a to programme control unit, ELA b order no. 08 5130 00	

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






* push-wire connector, 5 pole (00 0210 21)
wire cross section: 0.5 - 2.5 mm²

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Note for ComStation^T (77 0606 20) and RoomTerminal Flamenco, desktop installation (77 0521 00)! The connections NF 1a, NF 1b, NF 2a and NF 2b are not required in the following installations: There shall neither be connected a ComStation^{CT} (with speech, 77 0606 00) nor a ComTerminal Flamenco, desktop installation (77 0511 00) to this connection socket in the future.


B ComStation^{BUS-C}, Best.-Nr./order no. 77 0605 50
ComStation^{BUS-SE}, Best.-Nr./order no. 78 0605 50

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

* push-wire connector, 4 pole (00 0222 88)
wire diameter: 0.6 - 0.8 mm

3 0V
1 Beep
1Y(St)Y 2x2x0,8
Summer (Best.-Nr. 00 0420 28), bei Bedarf
Buzzer (order no. 00 0420 28), if required

3 0V
2 AW
1Y(St)Y 2x2x0,8
Zimmerleuchte Universal,
Best.-Nr. 77 0182 10
Room lamp Universal,
order no. 77 0182 10

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
**ComStation^{IP}, order no. 76 0605 50****1. Intended use**

The ComStation^{IP} is intended for plug-in connection to the IP network of a Flamenco^{IP+} nurse call system.

The ComStation^{IP} is a ward console with the following functions:

- Answering and processing of calls
- Announcements
- Switching ward couplings on and off
- Switching shifts on and off
- Optional free or discreet speech

The ComStation^{IP} is based on an IP desk telephone from manufacturer GRANDSTREAM with adapted Android operating system. The nurse call functionality is pre-installed as an application. All other functions are inactive (kiosk mode).

2. Important documents**Flamenco system family, Technical Manual (00 8813 97)**

The complete installation of the system is described in the technical manual. The manual can be downloaded from the www.tunstall.de website.

App help

The ComStation^{IP} app has an integrated help, which describes the operation of the device in detail.

User guide (00 8813 85)

A detailed user guide can be downloaded as PDF from the www.tunstall.de website.

Quick reference guide (00 8813 87)

A quick reference guide is included with the ComStation^{IP} and can be downloaded from www.tunstall.de.

Documentation by GRANDSTREAM

The ComStation^{IP} is supplied with documents from the device manufacturer GRANDSTREAM. Keep these. These documents are not intended for nursing staff.

3. Scope of supply (→ fig. A)

- [1] Main case
- [2] Handset
- [3] Handset cord
- [4] Ethernet cable
- [5] 12 V power adaptor
- [6] Wall mount
- [7] Screen cleaning cloth
- [8] Installation instructions (this document)
- [9] Quick reference guide
- [10] Documents by GRANDSTREAM

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4. Misuse

The ComStation^{IP} is based on an IP desk telephone from manufacturer GRANDSTREAM with adapted Android operating system. The nurse call functionality is pre-installed as an application. All other functions are inactive (kiosk mode).

Take note:

- The IP desk telephone may only be used as ComStation^{IP}. No other applications may be installed or used.
- The IP desk telephone must not be connected to the Internet. Only updates from Tunstall may be installed on the IP desk telephone.
- The Android operating system must not be reset to factory settings. This would result in the loss of changes to the operating system that can only be made by Tunstall.

5. Power supply

5.1 Power supply options

The device can optionally be supplied with power as follows:

- with the supplied 12 V power adaptor* (input: 100 – 240 V AC, 50 – 60 Hz; output 12 V DC, 1.5 A).
- in a LAN with PoE/PoE+ (PoE 802.3af Class 3, PoE+ 802.3at, Class 4).

*) Only use the 12 V DC power adaptor included in the scope of delivery of the device. Use of another power adaptor may damage the unit.

5.2 Power source for safety purposes

The German standard DIN VDE 0834 stipulates that in the event of a failure of the general power supply, the nurse call system must be supplied from a power source for safety purposes. This means that the 230 AC socket respectively the PoE/PoE+ network switch or injector must be connected to an UPS/auxiliary power supply.

6. Location for the device

- Minimum distance of 20 cm between the device and the body of the user.
- The area in front of the microphone (front right) should be kept clear. The optimum distance for hands-free talking is 50 cm.
- Do not place the device in a room with high dust generation, as this can considerably reduce the service life of the device.
- No direct sunlight or other heat radiation, as this can damage the electronic components and the plastic housing.
- Not in environments where steam is generated.

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7. Connecting the device (→ fig. B)

7.1 Ports

- [1] Power jack for 12 V DC power adaptor (only required if the device is not powered via the Ethernet cable (PoE))
- [2] LAN port
- [3] RJ9 handset connector port

7.2 Proceed as follows

A **prerequisite** for using the ComStation^{IP} is that the nurse call system has been configured for the use using the SystemOrganizer configuration software.

1. Connect the handset cord to the handset and to the port for handset cord [3].
2. Connect the Ethernet cable to the LAN port [2] and to the RJ45 connection box in the room. The RJ45 connection box must be connected to the nurse call system according to the Technical Manual.

CAUTION! Make sure that the latch of the RJ plugs engages so that the plugs cannot slip out.

3. If the device is not powered by PoE or PoE+ (PoE network switch), connect the 12 V DC power adaptor to the power jack [1] and to a 230 V socket.

CAUTION! DO NOT switch off the device during system start-up or when updating the firmware. Doing so may damage the firmware.

- ✓ The device is starting up. The ComStation^{IP} app login screen appears automatically.
- 4. Tap “LOG IN”.
- ✓ You will be logged in. The Home screen of the ComStation^{IP} app appears. The ComStation^{IP} is operating.

CAUTION! You have to secure all plugs (Ethernet cable, 12 V power adaptor) against unintentional removal. It is recommended to label the RJ45 connection box and the 230 V socket with a warning: “DANGER! Do not pull out the plug! Device is part of the nurse call system”.

For information on how to operate the ComStation^{IP} or change default settings, please refer to the App Help. To open the help, tap on the question mark in the home screen.

Functional checks must be carried out as described in the technical manual, before the ComStation^{IP} can be used.

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8. Desktop installation / wall mounting (→ fig. C, D)

8.1 Desktop installation


The device has an integrated stand [2]. To use it, fold it out of the back of the device and adjust it to the desired angle. Ensure that the device is in a secure position.

8.2 Wall mounting

1. Attach the wall mount to the two mounting slots [1] on the back of the device.
2. Attach the device to the wall via the two key holes in the wall mount.
3. Pull out the tab from the handset cradle [3].
4. Rotate the tab and plug it back into the slot with the extension up [4].

8.2 Wall mounting

1. Attach the wall mount to the two mounting slots [1] on the back of the device.
2. Attach the device to the wall via the two key holes in the wall mount.
3. Pull out the tab from the handset cradle [3].
4. Rotate the tab and plug it back into the slot with the extension up [4].

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- [1] Four mounting slots for fixing the wall mount
- [2] Built-in stand for desktop installation
- [3] Handset cradle
- [4] Tab in the position for wall mounting (tab with extension up)
- [5] Tab in the position for desktop installation (tab with extension down)

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9. Technical data

Dimensions (HxWxD)	84 x 252 x 211 mm
Weight	1080 g
Power consumption	– via PoE: approx. 11 W – via Power adaptor: approx. 9 W
Ambient temperature:	– Transport, storage – Operation
	-10 °C to +60 °C 0 °C to +40 °C
Relative humidity	10 % – 90 % (non condensing)

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Connection socket ComStation^{PC}, order no. 77 0452 60 A

Connection socket ComStation^{PC}, order no. 77 0452 60 F


Connection socket ComStation^{PC}, order no. 74 0452 60 A

The connection sockets differ only in their external design. All connection sockets are designed for the connection of a ComStation^{PC}, order no. 77 0602 10.


- Preferably made for wall mounting.
- Suitable for 2-gang back box (back box not included with delivery).
- 15 pole socket as data socket, prepared for 15 pole plug with mechanical locking mechanism.
- 25 pole SUB D socket, V.24.
- 4 plug-in terminals for connection of the 24 V supply voltage, wire cross section up to 2.5 mm².
- Connection as plug-in terminals for solid wire of ø 0.55 - 1 mm.
- Contact current: max. 3 A.

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








Note! Tighten the screws of the plugs
to protect the plugs against unintentional disconnection.

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ComStation^{CT} Flamenco, order no. 77 0606 00

Terminal with two-way speech communication, designed for desktop installation in the nurse station, incl. red Call Key, blue Alarm Key, green Presence Key (staff 1), yellow Presence Key (staff 2), 4 Function Keys as soft keys, microphone and loudspeakers.


NOTE! The complete installation of the nurse call system is described in the Technical Manual.

Glass – fragile! The operating front cover of the ComStation^{CT} is made of glass. Handle with care!

Installation

1. Connect the connection plug [2] of the ComStation^{CT} Flamenco [1] to the connection socket ComStation [3].
2. Tighten the two screws of the connection plug [2], to protect the connection plug against unintentional disconnection.

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- 1 - ComStation^{CT} Flamenco 3 - * Connection socket ComStation
2 - Connection plug

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

* Not included with ComStation^{CT} delivery.

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Symbols in the configuration menu

●	Circle key	Set the marked menu item.
▲	Upwards arrow key	Scroll upwards through the list.
▼	Downwards arrow key	Scroll downwards through the list.
⌂	Home key	Cancel, without storing.

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Menu items

Language:	Selecting the user language.
RAN > RAN number:	Setting of number of room devices (= RAN number).
RAN > Test RAN:	Test if room devices are ready to operate and if they are correctly connected to the RAN.
RAN > Status:	Check whether the stored RAN number is equal to the number of operational room devices.
Address:	Setting of room address.
Cleaning time:	Setting the time (seconds) how long it takes to wipe the front of the ComStation ^{CT} .
Contrast:	Function use only for Tunstall technicians.
Key sound:	Switching the key sound on or off. "Key sound ON" (factory setting) is recommended.
Fault tone:	Setting the fault tone volume.
Info:	Displaying the ComStation's software revision.
Audio Test:	Function use only for Tunstall technicians.
Key Function:	Setting the function of Alarm Key, Presence Key 1 and Presence Key 2 on the ComStation ^{CT} .
Reset:	Function use only for Tunstall technicians. (Restart the ComStation ^{CT}).

2. Entering of necessary settings

Selecting the user language

1. Mark "Language" using the arrow keys. Then press the circle key.
2. Mark the desired language using the arrow keys: D = German; GB = English, etc.
3. Press the circle key to set the selection.

Setting of RAN number (0 – 30) (Room Area Network)

RAN number = Number of devices in the room (switches, connection sockets, room lamps, etc.) that are connected via RAN. (**No** pear push switch, **no** patient handset).

1. Mark "RAN" using the arrow keys. Then press the circle key.
2. Mark "RAN number" using the arrow keys. Then press the circle key.
3. Mark the desired RAN number using the arrow keys.
4. Press the circle key to set the selection.

Setting of room address (0 – 110)

1. Mark "Address" using the arrow keys. Then press the circle key.
2. Mark the desired room address using the arrow keys.
3. Press the circle key to set the selection.

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Setting the key function

CAUTION! Switched off keys are not available for nursing staff!

Alarm Key

With the factory settings the alarm key on the ComStation^{CT} is operational and is initiated by pressing it.

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CAUTION! Setting a delay time does not comply with the German standard DIN VDE 0834. A delay time may only be set, where compliance with the DIN VDE 0834 is not required.

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NOTE! Die Alarm Key is only active while the staff presence is switched on in the room.

Presence Key 1 / Presence Key 2

With the factory settings the presence keys are operational. In exceptional cases it may be useful to switch of the key function for these keys.

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ComStation^T Flamenco, order no. 77 0606 20


Terminal without speech communication, designed for desk-top installation in the nurse station, incl. red Call Key, blue Alarm Key, green Presence Key (staff 1), yellow Presence Key (staff 2), and 4 Function Keys as soft keys.

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Installation

1. Connect the connection plug [2] of the ComStation^T Flamenco [1] to the connection socket ComStation [3].
2. Tighten the two screws of the connection plug [2], to protect the connection plug against unintentional disconnection.

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
Funktionstasten

- 1 - ComStation^T Flamenco 3 - * Connection socket ComStation
2 - Connection plug

* Not included with ComStation^T delivery

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Menu items

Language:	Selecting the user language.
RAN > RAN number:	Setting of number of room devices (= RAN number).
RAN > Test RAN:	Test if room devices are ready to operate and if they are correctly connected to the RAN.
RAN > Status:	Check whether the stored RAN number is equal to the number of operational room devices.
Address:	Setting of room address.
Contrast:	Function use only for Tunstall technicians.
Key sound:	Switching the key sound on or off. "Key sound ON" (factory setting) is recommended.
Fault tone:	Setting the fault tone volume.
Info:	Displaying the ComStation's software revision.
Reset:	Function use only for Tunstall technicians. (Restart the ComStation ^T).

2. Entering of necessary settings

Selecting the user language

1. Mark "Language" using the arrow keys. Then press the circle key.
2. Mark the desired language using the arrow keys: D = German; GB = English etc.
3. Press the circle key to set the selection.

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3. Entering of optional settings

The remaining menu items provide optional settings. Set these if required. Example:

Fault tone volume soft, medium, loud or off

WARNING! The tone that attracts attention for a fault may only be turned off, if it is made sure, that the attention is attracted in another way.

1. Mark "Fault Tone" using the arrow keys. Then press the circle key.
2. Mark "Soft", "Medium", "Loud" or "Off" using the arrow keys.
3. Press the circle key to set the selection.

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4. Checking the room bus RAN

Status

Check whether the set RAN number (refer to section “Setting of RAN number”) equals to the number of operational room devices.

1. Mark “RAN” using the arrow keys. Then press the circle key.
2. Mark “Status” using the arrow keys. Then press the circle key.

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
3. Compare the number of flashing room devices with the set RAN number. If the set RAN number is wrong, correct the RAN number setting after the RAN test as described in section “Setting of RAN number”.
4. End of test: Press the home key.

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
Note ! Tighten the screws of the plugs
to protect the plugs against unintentional disconnection.

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Telephone interface relay, order no. 11 5350 00

for connecting analogues telephones to a nurse call system. A telephone call will be displayed as a call in the nurse call system.

- Call voltage: 32 - 80 V AC
- Power supply: 5 - 40 V
- ON/OFF switch (ON = EIN, OFF = AUS)
- Output: N/C, N/O contact
- Bridging of call intervals can be set
- Power-on time can be set.

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
A Mounting

1. Power down the equipment.
2. Place the telephone interface relay on the top hat rail and click it into place.
3. Connect as shown in the diagram below.
4. Restart the equipment.
5. Switch on the telephone interface relay by pressing down the ON/OFF switch.

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B Potentiometer Setting

When the telephone interface relay is operational, the switching contact is closed as soon as the telephone signal voltage is applied. It depends on the setting of the two potentiometers "bridging of call intervals" and "power-on time" how long the switching contact remains closed.

Setting bridging of call intervals ("Rufpausenüberbrückung")

The potentiometer for the bridging of call intervals sets the period the switching contact remains closed longer than the telephone interface relay is actuated. The setting range is linear and lies between 0 and 12 seconds. If the selected time is longer than the actuation intervals by the telephone call voltage the switching contact remains closed without interruption until the actuation stops (end of call). If the selected time is shorter than the actuation interval the switching time of the contact results from the signal period plus the selected bridging time.

Setting power-on time ("Einschaltdauer")

The potentiometer for power-on time setting allows to limit the time the switching contact remains closed. The setting range is not linear and lies between 0.25 and about 12 seconds. If the potentiometer is set for permanent ("Dauer") the switching contact is not limited and remains closed up to the end of actuation.

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Diagnostic connection cable, order no. 70 0812 10


2 m connection cable for connecting a medical electrical device to the nurse call system via a connection socket, order no. 70 xxxx xx, in order to trigger a diagnostic call in the nurse call system for alarms of the medical electrical equipment.

WARNING! The transmission of alarm conditions of medical electrical equipment (e.g. monitoring devices) to the nurse call system serves only as supporting, additional information. This is a distributed information system.

- Medical staff must **not** rely on the nurse call system displays.
- Medical staff must not leave the hearing and visual range of the medical monitoring equipment.
- Due diligence for the operation of such medical electrical equipment remains unaffected in case of the connection to the nurse call system.

Prepare the diagnostic call

1. Plug the diagnostic connection cable into the connection socket 70 xxxx xx (socket for pear push switch).

Diagnostic call is raised

If a defined threshold value is exceeded on the medical electrical device, a diagnostic call is triggered in the nurse call system:

- The room lamp is flashing red. A diagnostic call is displayed on the call handling console.
- In case the diagnostic call is not answered in a pre-set time, it is forwarded to all rooms where the presence button is activated.
- If the presence button is activated at the call location, not a fresh call but only an answered call is raised, i.e. the call will not be forwarded.

Cancel the diagnostic call

1. Enter the call location and activate the presence button.
2. Reset the medical electrical device.
3. De-activate the presence button.

The call is cancelled. The call is no longer displayed.

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RAN interface universal, order no. 70 0848 00

Interface for output of actuator signals in connection with patient handset, order no. 74 0747 00 (as of revision D1) and one of the following connection sockets:
70 0424 00 (as of revision H1), 70 0424 50 (as of revision C1),
70 0434 00 (as of revision I1), 70 0434 50 (as of revision B1).

The revision no. is written on the product label behind the order number.

The RAN interface universal provides switching outputs for controlling actuators, such as a blinds control system.

Note! This RAN user has not to be registered at the room terminal (e.g. ComTerminal). If a RAN interface universal is installed, do not change the RAN number at the room terminal.

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Caution! The printed circuit board includes electrostatic sensitive components. Avoid touching.


C Setting bed number (if required)

In the factory setting the connected actuators (e.g. blinds control system) can be activated from all beds. If the actuators shall be controlled by one bed only, the bed number has to be set using the jumpers on the backside of the PCB.


1. Unscrew one casing cover.
2. Pull out the PCB.
3. Set jumpers P1, P2, P3 according to figure C.

D Mounting

Mounting on 35 mm supporting rail by using the snap-on element on the backside of the RAN interface.

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Interface for controlling the TV set with infrared signals. Integrated audio amplifier. Supported TV sets from Philips, Samsung, LG and Grundig see installation instructions for PCM-U "P-Labor Medical Card - Universal".

Scope of delivery:

- Control module
- IR emitter
- Velcro strip
- 3 cable ties, 1 adhesive socket
- Audio cable, 1.50 m
- Network patch cord, 1 m

Configuration


You have to configure the IR TV control module universal as described in the installation instructions for PCM-U "P-Labor Medical Card - Universal". To do this, you need the IR TV installation kit, order no. 77 0360 40; one kit per nurse call system.

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
* Scope of delivery of IR TV installation kit (77 0360 40).

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
1. Mounting

1. Each control module[1] has a built-in RAN interface. In an installation with more than one control module per room, however, the RAN interface may only be active on bed 1. This means, open the housings of the control modules of bed 2, 3, and 4 and disable the RAN interface by disconnecting the connector (see fig. D).
2. Connect the cables to the control module[1] as shown in fig. E.
3. Attach the control module to the back of the TV set using the Velcro strip, see fig. A.
4. Affix the IR emitter [2] to the front of the IR receiver of the TV set.
5. Connect the audio cable [4] to the headphone jack of the TV set.
6. Connect the network patch cord to the RJ45 connection socket.
7. Configure the IR TV control module according to the installation instructions for PMC-U "P-Labor Medical Card - Universal".

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
- 8 - Audio cable connection to headphone jack of the TV set
- 9 - Potentiometer for setting the IR emitting power, left: more, right: less.
- 10 - Service connection for temporarily connection of a TV loader or PC
- 11 - IR emitter connection (left: signal, right: ground)

- 12 - Network patch cord connection to RJ45 connection socket
- 13 - Power LED lights up when voltage is applied.
- 14 - Mode LED
- LED is on: TV set is switched on.
- LED is off: TV set is switched off.
- Further signals: Refer to installation guide PMC-U.

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
TV1 TV2 TV3 TV4
 Bett 1*)
 Bett 2
 Bett 3
 Bett 4


Bed number and TV tone channel

Set the bed number and
the TV tone channel by
jumpers on the
connection socket
(70 0435 00
or 70 0434 00):

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TV audio amplifier, order no. 77 0365 00


Connection socket combi (70 0425 00 or 70 0424 00) and Connection socket combi, bedhead unit (70 0435 00 or 70 0434 00) allow in connection with TV audio amplifier (77 0365 00) the TV tone transmission to Patient handsets (74 0747 00). Note! The TV is not controlled using the Patient handset.

Connection between TV audio amplifier and connection socket via RJ45 connection socket (00 0280 40 oder 00 0280 39).

Up to 4 TV sets, i.e. 4 TV audio amplifiers, can be installed per room. The TV tone of a TV set can be transmitted to 4 Patient handsets, i.e. 4 beds.

A galvanic separation point according to DIN VDE 0834-1 and EN 60601-1 is integrated in the TV audio amplifier.

- The TV set used must have a headphone jack.
- A second 230 V socket is required per TV set for the TV audio amplifier.
- Dimensions (HxWxD): 27 x 86 x 70 mm

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- | | |
|--|---|
| 1 - TV audio amplifier | 9 - Power LED |
| 2 - Audio cable (0.5 m) | 10 - 12 V power supply with DC plug [5] connection |
| 3 - Network patch cord (1.5 m) | 11 - Audio cable connection to headphone jack of the TV set [2] |
| 4 - Power cable with Euro plug (2 m) | 12 - Volume control (clockwise = Increase volume) |
| 5 - 12 V power supply with DC plug (0.5 m) | 13 - Network patch cord [3] connection |
| 6 - 2 Velcro strips | |
| 7 - 2 cable ties | |
| 8 - 1 adhesive socket | |

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1. Connect the cables to TV audio amplifier [1] as shown in fig. B.
2. Attach the TV audio amplifier [1] and the 12 V power supply [5] to the back of the TV set using the Velcro strips [6].
3. Connect the audio cable [2] to the headphone jack of the TV set.
4. Connect the power cable [4] to the 12 V power supply [5] and to a 230 V socket.
5. Connect the network patch cord [3] to the RJ45 connection socket.

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

TV1 TV2 TV3 TV4
Bett 1*)
Bett 2
Bett 3
Bett 4


Bed number and TV tone channel

Set the bed number and the TV tone channel by jumpers on the connection socket (70 0435 00 or 70 0434 00):

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RAN interface, order no. 77 0840 00

The RAN interface is intended for the connection of a third-party device to the room bus (RAN). Available applications:

- An external call device raises the call type "call".
- An external call device raises the call type "cardiac alarm".
- An external call device raises the call type "WC call".
- An external presence switch for staff presence 1.
- A telephone call raises the call type "telephone call".
- Patient TV sets are controlled using a patient handset.

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
CAUTION! During installation the power supply should be switched off to ensure zero-potential status. Otherwise, the RAN interface could be damaged.

A Mounting


Fix the RAN interface [1] with the adhesive pads [5] to the top-hat rail clip [4]. Then clip it to the top-hat rail. See Fig. A.

A

- 1 - RAN-Schnittstelle
- 2 - Anschlussklemme
- 3 - DIP-Schalter
- 4 - Hutschieneclip
- 5 - Klebepads

- 1 - RAN interface
- 2 - Connector
- 3 - DIP switch
- 4 - Top-hat rail clip
- 5 - Adhesive pads

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ON

= OFF

RAN interface for external call device, Order no. 77 0881 00

Intended Use

The RAN interface is used to connect the call function of a patient operating device from a third-party manufacturer to the RAN room bus of a nurse call system from the Flamenco system family.

Use of a patient operating device only after approval by Tunstall GmbH.

Dimensions (HxWxD): 13 x 51 x 95 mm (without fastening clips)

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
CAUTION! Only carry out the installation in a de-energised state. Otherwise the RAN interface could be damaged.

Mounting

Clip the RAN interface with the two fastening clips onto a top-hat rail, e.g. in the medical supply unit.

- [1] Connection terminal for RAN room bus
- [2] DIP switch
- [3] Connection terminal for external system
- [4] Two fastening clips for mounting on a top-hat rail

Setting the DIP switch


DIP switches 1–3:

Set the bed number:

ON

Bed 1

Bed 2


Bed 3

Bed 4

Bed 5

Bed 6

DIP switch 4:


Contact type of the connected call button:

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


Radio receiver-T, order no. Z 00 8202 31

Plug-in connection to:

Call switch with connection socket: 19 0704 00, 10 0704 00

Call switch with 2 connect. sockets: 19 0704 20, 10 0704 20

Radio receiver-T, order no. Z 00 8202 32

Plug-in connection to:

Connection socket L200: 73 0455 00 (rev. 4), 73 0400 00

Connection socket L200 bedhead unit: 73 045 00 (rev. 4)

Radio receiver-T, order no. Z 00 8202 33

Plug-in connection to:

Connection socket combi: 70 0425 00, 70 0424 00

Connection socket combi, TVL: 70 0425 50, 70 0424 50

Connection socket combi universal: 70 0449 00

Conn. socket combi, bedhead unit: 70 0435 00, 70 0434 00

C. socket combi, bedhead unit, TVL: 70 0434 50, 70 0435 50

C. socket combi, bedhead unit, universal: 70 0448 00

Connection socket with call switch: 70 0171 60A, 70 0171 60C, 70 0171 60F, 70 0171 03, 70 0171 00


Connection socket with call switch, bedhead unit: 70 0171 50

Connection socket call devices: 70 0400 00

Radio receiver-T, order no. Z 00 8202 36

Plug-in connection to:

Call switch insert with connector: 29 0704 00...

The radio receiver-T is able to receive signals from the following radio transmitters:

- Radio trigger including call button: MyAmie (P68007/02, P68007/04), see page 3.
- Radio trigger incl call button and fall detection: iVi™ (P68005/47)
- Universal sensor (61005/30)
- Wireless step-on sensor mat (Z00800301, Z00800302)
- Radio sensor mat 869 MHz (Z00800201)
- Large surface radio pneumatic switch (75071100)
- Radio smoke detector (68005/96)

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The radio receiver-T is plug connected to a connection socket of the nurse call system in the socket for pear push switch. Activating the radio transmitter will raise the call in the nurse call system that would be raised by a pear push switch in the same socket.

NOTE! The radio transmission is not monitored. That's why the radio transmitters may only be used as additional call devices in the nurse call system.

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Check radio coverage on the location of use

The radio range between the radio receiver-T and the radio transmitters depends on the structural conditions of the building; the range is up to 30 m. Before handing over a handheld radio transmitter (MyAmie, iVi) to the user, you must check the radio coverage at the place of use. To do this, walk around the planned place of use and trigger a call over and over again. Inform the user of the radio coverage.

Documents for VarioRec6 / radio receiver-T

The radio receiver-T is an OEM product of Lehmann Electronic GmbH. In addition to these installation instructions, you must read and observe the following documents from Lehmann Electronic:

- Brief overview VarioRec6 radio receiver (supplied with the radio receiver-T), document no. LE235
- Operating instructions VarioRec6, document no. LE243

The product designation VarioRec6 is a collective term of Lehmann Electronic for various radio receivers. The radio receiver-T is only one of these. Not all topics in the VarioRec6 documents apply to the radio receiver-T. Therefore, please observe the following notes:

Available functions

The following functions mentioned in the documents for VarioRec6 are available on the radio receiver-T with the Tunstall radio transmitters:

- Teaching in radio transmitters
- Teaching out radio transmitters
- Teaching out all taught-in radio transmitters
- Acknowledge fault message
- Set master mode
- Care mode
- Set operating display active / inactive
- Reset all parameters to factory settings
- Daily message monitoring
- Vital monitoring

Functions not available

The following functions mentioned in the documents for VarioRec6 are **not** available on the radio receiver-T with the Tunstall radio transmitters:

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When used with Tunstall radio transmitters, the radio receiver-T displays the following error and fault messages mentioned in the documents for VarioRec6:

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Information about low battery level

The radio transmitters contain batteries. When the battery level of a radio transmitter is low, the LED on the radio receiver-T will flash in red to inform the nursing staff, that the battery or the radio transmitter must be changed. Please refer to the documents for the VarioRec6.

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Note: If the radio transmitters are used with a social alarm unit a call is raised to a monitoring centre in case of a low battery level. This call is **not** raised when the radio transmitters are used with the radio receiver-T!

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- LED is on (for approx. 3 seconds) after the call button has been pressed: Battery is okay. A call is raised.
- LED is flashing after the call button has been pressed: The battery level is low. A call is raised. The LED on the radio receiver-T is flashing red to indicate, that the battery is low and the MyAmie must be changed.
- LED remains dark after the call button has been pressed: The battery is empty or the MyAmie is defective! No call will be raised! The MyAmie must be changed immediately.

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Radio receiver-T UP, order no. Z 00 8202 35

The radio receiver-T UP is suitable for connection to a room terminal in a Flamenco^{IP}, Flamenco, CONCENTO^{CARE}, or CONCENTO^{PLUS} nurse call system. It receives signals from the following radio transmitters:

- Radio trigger including call button: MyAmie (P68007/02, P68007/04)
- Radio trigger including call button and fall detection: iViTM (P68005/47)
- Universal sensor (61005/30)
- Wireless step-on sensor mat (Z00800301, Z00800302)
- Radio sensor mat 869 MHz (Z00800201)
- Large surface radio pneumatic switch (75071100)
- Radio smoke detector (68005/96)

The use of further radio transmitters is possible. For this you have to contact Tunstall GmbH.

Activating the radio trigger will raise a call in the nurse call system. The type of call triggered depends on how the radio receiver is connected.

NOTE! The radio transmission is not monitored. That's why the radio transmitters may only be used as an additional call device in the nurse call system.

The radio receiver-T UP is intended for flush-mounting. A frame with 55 mm inner dimension is required.


Potential-free connection as normally closed or normally open contact.

NOTE! The complete installation of the nurse call system is described in the respective Technical Manual.

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
1. Mounting (see fig. B)

1. Strip the nurse call system's connection cable in the back box [1] to a suitable length. Remove the cable shield and drain wire up to the cable sheath. Strip the wires.
2. Unscrew the two socket head screws [9] using an Allen key (H2 type) until the radio receiver insert [8] detaches from the mounting plate [3].
3. Remove the radio receiver insert [8] from the mounting plate [3].
4. Use the back box screws [2] to screw the mounting plate [3] onto the back box [1].
5. Insert the connection wires [7] of the radio receiver through the intermediate frame [6] (if an intermediate frame is used), the frame [4], and the mounting plate [3] into the back box [1].
6. Connect the connection wires [7] according to figure:
Fig. C: Flamenco^{IP} or Flamenco system: Connection via RAN interface 77 0840 00.
Fig. D: CONCENTO^{PLUS} system, Room without Control Module ZLB: Connection via 3 pole connector, e.g. WAGO 221.
Fig. E: CONCENTO^{CARE} or CONCENTO^{PLUS} system, Room with control module ZLB: Connection via RAN interface 19 0840 00.
7. Insert the radio receiver insert [8] into the frame [4] and the intermediate frame [6] and screw it onto the spacing bolt with core thread [5] using the two socket head screws [9].
8. Program radio receiver-T UP as described in the documents for VarioRec6 U, see section 5 in these installation instructions.

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- [1] * Back box, 1-gang (solid
wall: 17 0100 00, parti-
tion wall: 17 5100 00)
- [2] * Two back box screws
- [3] Mounting plate
- [4] * Frame, e.g. 77 0210 53
- [5] Two spacing bolts with
core thread

- [6] * Intermediate frame,
e.g. 77 0210 56; not re-
quired for all frame
types.
- [7] Five connection wires
- [8] Radio receiver insert
- [9] Two socket head screws
(H2)
- * Not included with radio re-
ceiver delivery

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4. Anschlüsse


2. Check radio coverage on the location of use

The radio range between the radio receiver-T UP and the radio transmitters depends on the structural conditions of the building; the range is up to 30 m. Before handing over a hand-held radio transmitter (MyAmie, iVi) to the user, you must check the radio coverage at the place of use. To do this, walk around the planned place of use and trigger a call over and over again. Inform the user of the radio coverage.

3. Dismantling (see fig. B)

1. Unscrew the two socket head screws [9] until the radio receiver insert [8] detaches from the mounting plate [3].
2. Disconnect the connection wires [7].
3. Remove the radio receiver insert [8] together with the frame [4] and the intermediate frame [6] (if present).
4. Unscrew the screws [2] on the mounting plate [3] until you can remove the mounting plate [3].

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

5. Documents for VarioRec6 U / radio receiver-T UP

The radio receiver-T UP is an OEM product of Lehmann Electronic GmbH. In addition to these installation instructions, you must read and observe the following documents from Lehmann Electronic:

- Brief overview VarioRec6 U radio receiver (supplied with the radio receiver-T UP), document no. LE264
- Operating instructions VarioRec6 U, document no. LE265

The product designation VarioRec6 U is a collective term of Lehmann Electronic for various radio receivers. The radio receiver-T UP is only one of these. Not all topics in the VarioRec6 documents apply to the radio receiver-T UP. Therefore, please observe the following notes:

Available functions

The following functions mentioned in the documents for VarioRec6 U are available on the radio receiver-T UP with the Tunstall radio transmitters:

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Functions not available

The following functions mentioned in the documents for VarioRec6 U are **not** available on the radio receiver-T UP with the Tunstall radio transmitters:

- Radio transmitter for special functions, e.g. radio cancel switch or radio light switch

Error/failure messages during operation

When used with Tunstall radio transmitters, the radio receiver-T UP displays the following error and fault messages mentioned in the documents for VarioRec6 U:

- Transmitter battery low
- System fault
- Reception blockade
- Daily report missing
- Vitality message is missing

6. Information about low battery level

The radio transmitters contain batteries. When the battery level of a radio transmitter is low, the LED on the radio receiver-T UP will flash in red to inform the nursing staff, that the battery or the radio transmitter must be changed. Please refer to the documents for the VarioRec6 U.

ATTENTION! The nursing staff must identify and acknowledge a low battery level by watching the LED on the radio receiver-T and/or the signals on the radio transmitter. Radio transmitters with a low battery level must be replaced immediately.

Note: If the radio transmitters are used with a social alarm unit a call is raised to a monitoring centre in case of a low battery level. This call is **not** raised when the radio transmitters are used with the radio receiver-T UP.

LED on the MyAmie, order no. P68007/02, P68007/04

The MyAmie is used to raise calls by pressing the call button. The LED on the MyAmie confirms, that the button is pressed, and indicates the battery level at the same time:

- LED is on (for approx. 3 seconds) after the call button has been pressed: Battery is okay. A call is raised.
- LED is flashing after the call button has been pressed: The battery level is low. A call is raised. The LED on the radio receiver-T UP is flashing red to indicate, that the battery is low and the MyAmie must be changed.
- LED remains dark after the call button has been pressed: The battery is empty or the MyAmie is defective! No call will be raised! The MyAmie must be changed immediately.

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Loudspeaker with announcement interface, order no. 05 0024 02

Wall mounted speaker in a full metal design suitable for announcements of the nurse call system. Designed for installation in corridors and large rooms.

Adjustable amplifier for announcements from the nurse call system. Transformer for entertainment transmission.

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
CAUTION! The printed circuit board includes electrostatic sensitive components. Avoid touching.

Mounting

1. Put the connection cables through the cut-out of the rear panel.
2. Fix the rear panel to the wall with the attached screws.
3. Connect the cables according to figure B.
4. Put the housing of the speaker on the rear panel on the wall.
5. Fix the loudspeaker at the right and left side with the attached screws.

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

Technical data

Power	10/6/3/1.5 W
Frequency range	170 – 18500 Hz
SPL:	
- 1 W / 1 m, peak	98,3 dB
- 1 W / 4 m, peak	86,3 dB
- P _{max} /4 m, peak	94,0 dB
Standards	EN 54-24; building material class B2 acc. to DIN 4102-2
Supply voltage	24 V DC
Standby current consumpti.	80 mA
Dimensions (HxWxD)	193 x 253 x 82 mm
Weight	2.2 kg
Housing material	Metal
IP rating	IP21
Ambient temperature	+5 °C – +40 °C

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Door entry speaker 2, order no. 77 0351 00

1. Intended use

Door entry speaker 2 is intended for connection to OSYlink-Door entry speaker 2 (77 0801 10) in a nurse call system of the Flamenco system family. Power is supplied by the 24 V DC supply of the nurse call system. A door opener (for specification refer to page 3) can be connected to the door entry speaker and then be controlled via the nurse call system.

Limitations on use! The unit must not be used for life support systems or applications within nuclear facilities.

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
CAUTION! The printed circuit board includes electrostatic sensitive components. Avoid touching.

2. Product description (fig. A)

The call button [2] is used to initiate a door call in the nurse call system. Via microphone and loudspeaker, it is possible to speak with the nursing staff at the console. The nursing staff at the console can release the door opening mechanism. The door opener is active as long as the button at the console is pressed.

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- | | |
|------------------------------|---|
| [1] Front panel | [9] Terminal block |
| [2] Call button | [10] 3 washers (for screws with head Ø 7 mm, not included with scope of delivery) |
| [4] Lower shell | [11] Sealing block |
| [5] 2 locking screws | [12] Strain relief clamp |
| [6] 3 mounting holes | [13] Allen key |
| [7] Shaft for cable outlet | |
| [8] Connection circuit board | |

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3. Wall mounting

Wall mounting on an even surface. Uneven surfaces must be adjusted beforehand.

1. Drill three holes in the wall, see fig. C for the position of the mounting holes [6].
2. Open the housing by carefully levering off the lower shell [4], see fig. B.
3. Pull off the wired terminal block [9] on the back of the front panel [1] downwards, see fig. D. Put the front panel aside.
4. Lead the connection cable(s) out through the shaft [7] in the lower shell [4], see fig. D.
5. Fasten the lower shell [4] to the wall with screws (not included in the scope of delivery, head Ø 7 mm) and the washers [10] at the three mounting holes [6]. The washers must be inserted centred with the insulated side facing the lower shell [4]. Tighten the screws evenly and not too tightly using a hand screwdriver to avoid distortion of the lower shell.
6. The sealing block [11] has two lead-throughs for connection cables. By turning and pulling the protruding part, open the corresponding lead-through and guide the connection cable through the opening.
7. Insert the sealing block [11] into the opening provided on the lower shell [4] and then screw the strain relief clamp [12] tight with the screw provided, see fig. D.
8. Strip 6 mm of insulation from the wires of the connection cables and insert them into the connection terminals as shown in fig. D.
9. Push the wired terminal block [9] upwards again on the back of the front panel [1], see fig. D.
10. Install the cables in the lower shell pointing downwards, that they cannot be trapped when closing the shell.
11. It is essential to ensure that the seal is correctly placed in the lower shell [4] so that no liquids can ingress into the housing.
12. Before closing the housing, check that the alignment is even e.g. by applying the glass plate.
13. To close the housing, first insert the upper part of the front panel [1] into the catches in the lower shell [4] and then fold in the lower part of the front panel [1]. Do not use any force when doing this. The housing must close without resistance. If this is not the case, reinsert it and make sure that it is correctly seated in the upper catches and the seal.
14. Lock the housing on the underside by screwing in the two locking screws [5] with the Allen key [13]. Spring contacts! Tighten the screws sufficiently to ensure the device functions correctly.
15. Carry out a function test: Initiate door call by pressing the call button [2], answer the door call at a call handling console, check speech connection, check door opening mechanism.

4. Cleaning note

Liquids penetrating the housing can damage or destroy the device. Only clean the device with a dry cloth.

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5. Dismantling

1. Hold the front panel [1] with one hand.
2. With the other hand, unscrew the two locking screws [5] on the underside of the device a little with the Allen key until you can lever off the front panel [1], see fig. B.
3. Carefully lever off the front panel [1], see fig. B.
4. Pull off the wired terminal block [9] on the back of the front panel [1] downwards, see fig. D.
5. Disconnect the connected wires as follows: To release the push-in spring connection, insert a small screwdriver into the opening above the clamping point and then pull off the wire.
6. Loosen the three fastening screws and remove the lower shell [4] from the wall.

6. Connections (fig. D)

Connection to OSYlink-Door entry speaker 2	
Connector	Push-in spring connection
Conductor diameter	0.8 mm
Cable type	2-wire, 0.8 mm
Stripping length	6 mm
Max. cable length	10 m

Connection to door opener	
Connector	Push-in spring connection
Conductor diameter	0.8 mm
Cable type	2-wire, 0.8 mm
Stripping length	6 mm
Max. cable length	10 m

Door opener

The switching contact for the door opener provides an output voltage of 12 V DC and a current of max. 200 mA. The switching pulse lasts as long as the door opener button is pressed at the call handling console.


Door opener requirements

- Operating current type, i.e. the door opener unlocks when there is active electric current.
- Rated operating voltage: 8 – 24 V DC
- Max. current consumption: 200 mA
- DC resistance: min. 42 Ohm

Other door openers

If a door opener is to be used that does not meet the above requirements, an appropriate relay must be connected in between. The door opener must be supplied from a separate power supply. Requirements for the relay:

- Rated control voltage: 12 V DC
- Control current: max. 200 mA

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7. Technical data

Power supply from	from OSYlink-Door entry speaker 2
Current consumption of OSYlink-Door entry speaker 2 + Door entry speaker 2 + Door opener:	
- Standby consumption	80 mA
- During speech connection	130 mA
- Max. current consumption (= door opener actuation)	400 mA
Contact load capacity door opener	max. 200 mA, min. 42 Ohm
Separation between door opener and nurse call system	2 x MOPP according to EN 60601-1
Dimensions (HxWxD)	220 x 138 x 28 mm
Weight	approx. 750 g
Lower shell material	Plastic
Front panel material	6 mm thick safety glass
Degree of protection	IP44
Ambient temperature	
- Operation:	-20 °C – +50 °C
- Storage:	-20 °C – +70 °C
Relative humidity	non-condensing

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OSYlink-Door Entry Speaker, order no. 77 0801 00

Interface for connecting a door entry speaker, order no. 77 0350 00, to the group bus (OSYnet).

- 2-wire-connection to the door entry speaker (including data and speech)
- Connection for door opener transformer 12 V AC
- Input for local call cancelling in combination with a separate switch (N/O contact to 24 V)
- Max. current consumption: 380 mA

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A2: Wall mounting

1. Push out the three fixing clips 2 until the holes for the screws 4 appear.
2. Screw the OSYlink 1 to the wall.

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Examples: Address 1 is set by using the code switch 1. Address 4 is set by code switch 3. Address 24 is selected through addition, coding switch 4 and 5 ($8 + 16 = 24$).

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
OSYlink-Door entry speaker 2, order no. 77 0801 10

1. Intended use

Interface for the connection of a door entry speaker 2 (77 0351 00) to an OSYnet group bus in a nurse call system of the Flamenco system family. A door opener can be connected to the door entry speaker 2 and then be controlled via the nurse call system.

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2. Mounting (see fig. A)






Optionally mounting on top-hat rail or wall mounting.

A1: Mounting on top-hat rail (35 mm)

Click the OSYlink [1] onto the top-hat rail [3].

A2: Wall mounting

- Push out the three fixing clips [2] until the holes for the screws [4] appear.
- Screw the OSYlink [1] to the wall.

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- [1] OSYlink-Door entry speaker 2 [4] * 3 screws
[2] 3 fixing clips [5] * ESD protection cover
[3] * Top-hat rail (00 0276 53)

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3. Setting the address (1 – 110) (see fig. C)

For setting the address for the OSYlink on the OSYnet group bus use the code switches 1 to 7 with the values 1, 2, 4, 8, 16, 32, 64. Set the appropriate switch to ON.

Examples: Address 1 is set by using the coding switch 1, address 4 by the coding switch 3.

Address 24 is set through addition: coding switch 4 and 5 ($8+16=24$).

After the address is changed during operation a reset has to be carried out afterwards. Otherwise the new setting will not be activated.

Hardware reset

To reset the OSYlink you have to set code switch 10 to ON for one second, then back to OFF.

4. ESD protection cover (fig. D)

After connecting the cables (see fig. E), put on the ESD protective cover (order no. 00 0276 53).

The ESD protective cover serves as protection against electrostatic discharges and as protection against unintentional contact.

5. Dismantling from top-hat rail

Put a screw driver into the visible hole of the black fixing clip. Then pull the clip out upwards until the OSYlink [1] detaches from the top-hat rail [3].


6. Connections (fig. D)

OSYnet connection	
Connector	Plug-in screw-type connector
Cable type	J-Y(St)Y 4x2x0.8 or CAT7 (22 AWG), $\varnothing = 0.64$ mm
Stripping length	7 mm

Ipwr connection	
Connector	Plug-in screw-type connector
Conductor cross section	2,5 mm ²
Cable type	NYM 2x2.5 mm ²
Stripping length	7 mm

Connection to Door entry speaker 2	
Connector	Plug-in screw-type connector
Conductor diameter	0.8 mm
Cable type	2-wire, 0.8 mm
Stripping length	7 mm
Max. cable length	10 m

Connection to cancel switch (option)	
Connector	Plug-in screw-type connector
Conductor diameter	0.8 mm
Cable type	J-Y(St)Y 2x2x0.8
Stripping length	7 mm
Max. cable length	10 m

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7. Technical data

Power supply	24 V DC
Current consumption of OSYlink-Door entry speaker 2 + Door entry speaker 2 + Door opener:	
- Standby consumption	80 mA
- During speech connection	130 mA
- Max. current consumption (= door opener actuation)	400 mA
Dimensions (HxWxD)	90 x 106 x 58 mm
Housing material	Polycarbonate
Degree of protection	IP20
Ambient temperature	+5 °C – +50 °C
Relative humidity	0 % – 85 % (non-condensing)

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OSYlink-Group lamp, order no. 77 0802 00

Interface for connecting group lamps to the group bus (OSYnet). Usable as collective or direction display. Room allocation and group forming possible.

- 4 x 2 potential-free outputs for connecting 4 group lamps, with two light sections each
- Power supply selectable (jumper) internal from the call system or external
- Display compliant with DIN VDE 0834, centrally configurable
- Central, flexible configuration of flash-signal sequences and allocation of user addresses
- Max. current consumption: 26 mA + 60 mA per output

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Examples: Address 1 is set by using the code switch 1. Address 4 by the code switch 3. Address 24 is selected through addition, coding switch 4 and 5 ($8+16=24$).

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OSYlink-Universal, order no. 77 0803 00

Interface for connecting external systems and/or technical installations to the group bus (OSYnet).

Standby current consumption (all outputs OFF): 40 mA.

Max. current consumption (all outputs ON): 100 mA.

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A Mounting

Optionally mounting on supporting rail or wall mounting.

A1: Mounting on supporting rail (35 mm)

Click the OSYlink 1 onto the supporting rail 3.

A2: Wall mounting

1. Push out the three fixing clips 2 until the holes for the screws 4 appear.

2. Screw the OSYlink 1 to the wall.

B Setting of address (1 - 110)

For setting the address on the printed circuit board please use the code switches 1 to 7 with the values 1, 2, 4, 8, 16, 32, 64. Set the appropriate switch to ON.

Examples: Address 1 is set by using the code switch 1. Address 4 by the code switch 3.

Address 24 is selected through addition: coding switch 4 and 5 ($8+16=24$).

Do not change code switches 8 and 9! Factory set: OFF.

Reset

To reset the OSYlink you have to set code switch 10 to ON for one second, then back to OFF. A reset has to be made, if a code switch setting has been changed. Otherwise the new setting will not be activated.

C ESD protection cover (order no. 00 0276 54)


After you have connected the cables (see reverse side), you have to put the ESD protection cover on the OSYlink. (ESD = Electro Static Discharge)

D Dismantling from supporting rail


Put a screw driver into the visible hole of the black fixing clip. Then pull the clip upwards, until the OSYlink 1 comes loose from the supporting rail 3.

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



IN 1	Monitored input - Call	Normally open contact with parallel 4K7 resistor
IN 2	Monitored input - Emergency call	Normally open contact with parallel 4K7 resistor
IN 3	Monitored input - Alarm call	Normally open contact with parallel 4K7 resistor
IN 4	Monitored input - Call	Normally open contact with parallel 4K7 resistor
Function for IN 1 - IN 4: Contact closed = Call message Fresh call Contact opened = Off message Call cancelled Acknowledgement possible, fresh / answered		

IN 5	Input - Collective announcement (to all wards)	Normally open contact
IN 6	Input - Collective announcement (to all staff)	Normally open contact
Function for IN 5 - IN 6: Contact closed = Announcement started Contact opened = Announcement closed		

IN 7	Input - Raise call	Normally open contact
IN 8	Input - Cancel call	Normally open contact
Function for IN 7: Contact closed = Raise call Contact opened = no function		
Function for IN 8: Contact closed = Cancel call from IN7 Contact opened = no function		

OUT 1	Output can be configured. Factory setting: Collective output Calls	Solid state output 24V
OUT 2	Output can be configured. Factory setting: Collective output Emergency calls	Solid state output 24V
OUT 3	Output can be configured. Factory setting: Collective output Alarm calls	Solid state output 24V
OUT 4	Output can be configured	Solid state output 24V
OUT 5	Output can be configured	Solid state output 24V
OUT 6	Output can be configured	Solid state output 24V and/or location light feature for inputs IN 1, IN 2, IN 3, IN 4, IN 7.

OUT 7	Out. can be configured	Change-over contact potential-free, internal 24V can be used
OUT 8	Out. can be configured	Change-over contact potential-free, internal 24V can be used
The outputs are configured using the SystemOrganizer software. In case of ward coupling the outputs signal across wards.		

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OSYlink-Universal BMA, order no. 77 0803 01

Interface for connecting external systems and/or technical installations to the group bus (OSYnet). One input (IN 4) for activating the call type „Fire Alarm“, call category Alarm. The fire alarm message is active as long as the contact at the input IN 4 is closed. The acoustic signalling of the message can be acknowledged. The message is automatically cancelled, when the contact input falls back to its starting position.

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Examples: Address 1 is set by using the code switch 1. Address 4 by the code switch 3. Address 24 is selected through addition: coding switch 4 and 5 ($8 + 16 = 24$).

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IN 4 Monitored input - Fire alarm message
Normally open contact with parallel 4K7 resistor

Function for IN 4: Contact closed = Call message Fresh call
Contact opened = Off message Call cancelled
Acknowledgement possible, fresh / answered

IN 1 Monitored input - Call Normally open contact with parallel 4K7 resistor
IN 2 Monitored input - Emergency call Normally open contact with parallel 4K7 resistor
IN 3 Monitored input - Alarm Normally open contact with parallel 4K7 resistor

Function for IN 1 - IN 3: Contact closed = Call message Fresh call
Contact opened = Off message Call cancelled
Acknowledgement possible, fresh / answered

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IN 7 Input - Raise call Normally open contact
IN 8 Input - Cancel call Normally open contact

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OSYlink-Universal SZ, order no. 77 0803 02

Interface for connecting external systems and/or technical installations to the group bus (OSYnet). One input (IN 4) for activating a ward coupling (first ward coupling defined with the SystemOrganizer).

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Attention! The „ward coupling“ feature can only be used, if there is a ComStation^{BUS-C}, a ComStation^{BUS}, a ComStation^{PC}, or a ComStation^{CT} on the same ward as the OSYlink-Universal SZ (organisational structure in the SystemOrganizer).

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Examples: Address 1 is set by using the code switch 1. Address 4 by the code switch 3. Address 24 is selected through addition: coding switch 4 and 5 ($8+16=24$).

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IN 4 Monitored input - Ward coupling (first ward coupling set up in the SystemOrganizer)
Normally open contact with parallel 4K7 resistor

Function for IN 4: Contact closed = Ward coupling activated
Contact opened = Ward coupling deactivated

IN 1 Monitored input - Call Normally open contact with parallel 4K7 resistor
IN 2 Monitored input - Emergency call Normally open contact with parallel 4K7 resistor
IN 3 Monitored input - Alarm Normally open contact with parallel 4K7 resistor
Function for IN 1 - IN 3: Contact closed = Call message Fresh call
Contact opened = Off message Call cancelled
Acknowledgement possible, fresh / answered

IN 5 Input - Collective announcement (to all wards) Normally open contact
IN 6 Input - Collective announcement (to all staff) Normally open contact

Function for IN 5 - IN 6: Contact closed = Announcement started
Contact opened = Announcement closed

IN 7 Input - Raise call Normally open contact

IN 8 Input - Cancel call Normally open contact

Function for IN 7: Contact closed = Raise call
Contact opened = no function
Function for IN 8: Contact closed = Cancel call from IN7
Contact opened = no function

OUT 1 Output can be configured. Solid state output 24V
Factory setting: Collective output Calls

OUT 2 Output can be configured. Solid state output 24V
Factory setting: Collective output Emergency calls

OUT 3 Output can be configured. Solid state output 24V
Factory setting: Collective output Alarms

OUT 4 Output can be configured Solid state output 24V

OUT 5 Output can be configured Solid state output 24V

OUT 6 Output can be configured Solid state output 24V and/or location light feature for inputs IN 1, IN 2, IN 3, IN 4, IN 7.

OUT 7 Out. can be configured Change-over contact potential-free, internal 24V can be used

OUT 8 Out. can be configured Change-over contact potential-free, internal 24V can be used

The outputs are configured using the SystemOrganizer software.

In case of ward coupling the outputs signal across wards.

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Ward coupling in case of an event in the nurse call system

You can use OSYlink-Universal SZ to activate a ward coupling, if defined events happen in the nurse call system.

1. Connection

Connect input IN 4 to output OUT 8 as shown in the connection plan below.

The jumpers R_IN 4 and 24V_OUT 8 must not be set.

2. Setting up of ward coupling


Set up in the SystemOrganizer the ward coupling to be activated. For that you have to select the wards that shall be coupled. You have to set up this ward coupling before all other ward couplings. It will be on top of the ward couplings' list in the SystemOrganizer.

3. Configuring OSYlink-Universal SZ

In the SystemOrganizer you configure, which event type in the nurse call system shall activate the ward coupling. For this you configure output OUT 8 of the OSYlink-Universal SZ.

Proceed as follows in the configuration window:

1. „Outputs“: Select „OUT 8“.
2. „Call Types and Categories“: Select in case of which call type(s) and call categorie(s) the ward coupling shall be activated.
3. „Assigned Locations“: Select from the events of which rooms the ward coupling shall be activated.
4. „Blinking Rhythm activation“: Select „Inactive“.
5. „Fresh Filter“: Select any option.
6. „Show the inputs IN1-4, IN7 on the outputs“: This setting does not affect the ward coupling.

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Jumpers

R_IN4: Jumper must not be set.
24V_OUT8: Jumper must not be set.

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
OSYlink-Announcement, order no. 77 0804 00

Intended use

Interface for connecting a maximum of 5 loudspeakers with announcement interface, order no. 05 0024 01, 05 0024 02, or 05 0024 03 or a maximum of four 1-channel-audio-amplifiers 100 V / 25 W, order no. 00 0647 13, to the group bus (OSYnet).

- 1 audio output for driving the active loudspeakers or the audio amplifiers
- 1 potential free digital output for control of signalling the activity of the active loudspeakers
- Standby current consumption: 32 mA
- Max. current consumption: 60 mA

NOTE! The connection of loudspeaker with announcement interface, order no. 05 0024 02, is described on a separate installation instructions leaflet (00 8803 99).

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Mounting (→ fig. A)

Optionally mounting on supporting rail or wall mounting:

Mounting on supporting rail (35 mm) (→ fig. A1)

Click the OSYlink [1] onto the supporting rail [3].

Wall mounting (→ fig. A2)

1. Push out the three fixing clips [2] until the holes for the screws [4] appear.
2. Screw the OSYlink [1] to the wall.

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Setting of address (1 – 110) (→ fig. B)

For setting the address on the printed circuit board please use the DIP switches 1 to 7 with the values 1, 2, 4, 8, 16, 32, 64. Set the appropriate DIP switch to ON.

Examples

- Address 1 is set by using the DIP switch 1.
- Address 4 by the DIP switch 3.
- Address 24 is selected through addition: coding switch 4 and 5 ($8 + 16 = 24$).

Do not change DIP switches 8 and 9! Factory set: OFF.

Reset

To reset the OSYlink you have to set DIP switch 10 to ON for one second, then back to OFF.


A reset has to be made, if a DIP switch setting has been changed. Otherwise the new setting will not be activated.

ESD protection cover (00 0276 53) (→ fig. C)

After you have connected the cables (see reverse side), you have to put the ESD protection cover, order no. 00 0276 53 on the OSYlink (ESD = Electro Static Discharge).


Dismantling from supporting rail

Put a screw driver into the visible hole of the black fixing clip. Then pull the clip upwards, until the OSYlink [1] comes loose from the supporting rail [3].

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Connection of loudspeaker with announcement interface (order no. 05 0024 03 or 05 0024 01)

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
Anschluss von 1-Kanal-ELA-Verstärker 100 V / 25 W (Best.-Nr. 00 0647 13)

Connection of 1-channel-power-amplifier 100 V / 25 W (order no. 00 0647 13)

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OSYlink AS-CCS, order no. 77 0870 00

Interface for connecting the room terminals of a CCS ward (CCS bus) to the group bus (OSYnet) of a Flamenco nurse call system. Current consumption: 90 mA.
For connecting a half-duplex speech system (CCS 1080 W) additionally a speech amplifier, order no. 74 8000 00, with the power supply unit, order no. 21 8000 00, is needed.

OSYlink AS-VKS, order no. 77 0871 00

Interface for connecting the room terminals of a VKS ward (CCS bus) to the group bus (OSYnet) of a Flamenco^{SE} call system. Current consumption: 90 mA.
Additionally a speech amplifier, order no. 74 8000 00, with the power supply unit, order no. 21 8000 00, is needed.

OSYlink AS-L200, order no. 77 0872 00

Interface for connecting the room terminals of a L200 ward (ward bus WCB0) to the group bus (OSYnet) of a Flamenco nurse call system. Current consumption: 90 mA.

Note! The OSYlink module and connected room terminals must be configured using the SystemOrganizer of the Flamenco system.

Caution! Electrostatic sensitive components

The printed circuit board includes electrostatic sensitive components. Therefore avoid touching.

A Mounting

Optionally mounting on supporting rail or wall mounting in distribution boxes or technical rooms.

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A2: Wall mounting

- Push out the three fixing clips 2 until the holes for the screws 4 appear.
- Screw the OSYlink 1 to the wall.

B Setting the 10 pole coding switch

Setting the address (1 – 110)

For setting the address for the OSYlink on the group bus OSYnet use the code switches 1 to 7 with the values 1, 2, 4, 8, 16, 32, 64. Set the appropriate switch to ON.

Examples: Address 1 is set by using the code switch 1. Address 4 by the code switch 3. Address 24 is selected through addition: coding switch 4 and 5 (8 + 16 = 24).

Setting the connected system (according to figure B)


For setting the connected system use code switch 8 „W/G“ (W = simplex speech, G = duplex speech) and code switch 9 „1080/2000“.

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C System Status LEDs

OSYlink AS-CCS or OSYlink AS-VKS:


Flashing slowly (1s / 1s) = Connection to OSYnet is ok.
Steady light = No connection to OSYnet.

Flashing slowly (1s / 1s) = Speech line is operational.
Flashing quickly = Active speech connection or announcement.

Flashing slowly (1s / 1s) = CCS-Bus is ok.

Flashing quickly = Minimum one fault on the CCS bus.
Steady light = OSYlink AS is not operational! It is waiting for configuration data from the OSY-ControlCenter

OSYlink AS-L200:

Flashing slowly (1s / 1s) = Connection to OSYnet ok.
Steady light = No connection to OSYnet.

Flashing slowly (1s / 1s) = Ward bus WCB0 is ok.
Flashing quickly = Minimum one fault on the ward bus WCB0.
Steady light = OSYlink AS is not operational! It is waiting for configuration data from the OSY-ControlCenter

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IP-SystemManager, order no. 76 2100 00

The control of a Flamenco^{IP+} or Flamenco^{IP} nurse call system is performed by decentralised LAN networked IP-SystemManagers.

All IP-SystemManagers are equipped with identical hardware. However, each IP-SystemManager has a defined function. For this function the required software and a project specific database were factory installed. That means, each IP-SystemManager must be installed at that location, for which it has been prepared at the factory.

For installation location and IP address refer to the device label.

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CAUTION! The printed circuit board behind the service cover includes electrostatic sensitive components. Avoid touching.

A Product description

- | | |
|-----------------------|---|
| 1 - Housing | * Not included in the scope of delivery, please order separately. |
| 2 - Service cover | |
| 3 - * DIN rail | |
| 4 - 2 DIN rail clamps | |

B Mounting

Installation in dry rooms only. Always easily accessible to authorised persons (access space at least 60 cm wide).


The IP-SystemManager must be protected against unintentional contact, e.g. by installing it in a 19" wall mounting cabinet, 7HU, GD (order no. 77 3519 07).

In case of installation in switch cabinets or similar, dissipated heat shall be discharged by means of forced ventilation, if necessary.

Mounting on DIN rail according to fig. B. e.g. with 19" mounting set (order no. 76 0900 01).

C Dismantling

Dismantling from the DIN rail according to fig. C.

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Patchkabel, mindestens Kat. 5e.

Connections

★ The connections marked with a star are only enabled, if the IP-SystemManager was prepared accordingly at the factory, refer to the project specific documents.

Connection "Power"

Connection to the 24 V power supply ring of the ward or directly to a 24 V power supply unit. Insert insertion bridges!

Conductor cross-section: 2.5 mm²

Connector, 4-pole: spare part order no. 00 0211 33.

Two insertion bridges: spare part order no. 00 0220 52.

Connection "UPS Input"

Input for the signalling contact of a connected uninterruptible power supply (UPS), which signals buffer operation. This signalling contact is available at:

- UPS Control unit 10A (77 3411 00): "Buffering" contact.
- DC-UPC with capacitor storage 24V, 6kWs (77 3413 00): "Buffering" contact.
- Power supply unit UPS (77 3400 00): "Battery active" contact.
- Power supply unit UPS 60 (77 3400 10): "Battery active" contact.

Conductor cross-section: 0.14 mm² – 1.5 mm²

Connector, 2-pole: spare part order no. 00 0211 47.

★Connection "Fault Relay"

The signalling of the output is in parallel with the "Fault" LED (below "Fault Relay"), for details refer to the chapter "LED displays".

Fault message relay output.

Change over contact, potential free.

Maximum contact load: 2 A, 24 V

Conductor cross-section: 0.14 mm² – 1.5 mm²

Connector, 3-pole: spare part order no. 00 0211 45.

★Connections "OSYnet 1" and "NF 1"

Connection of a group bus OSYnet. For the maximum number of bus users that can be connected, please refer to the Planning manual.

CAUTION! The IP-SystemManager must not be installed in the middle of the group bus line. It must be connected at one end of the group bus line, as it includes an integrated bus terminating resistor.

In nurse call systems without speech communication, connection "NF 1" is not used.

For permitted cable types, refer to the Technical manual.

Connector, 3-pole, for OSYnet: spare part order no.

00 0211 45; Connector, 4-pole, for NF: Spare part order no. 00 0211 36.

★Connections "OSYnet 2" and "NF 2"

Not used.

Connection "Ethernet"

Connection to the Ethernet of the nurse call system. This connects the IP-SystemManager with the other IP-SystemManagers and the IP based devices of the nurse call system.

NOTE! The Ethernet port is internally equipped with a separation point according to EN 60601-1 (2 x MOPP). An external separation device is not required.

Ethernet, 10/100 Mbit LAN connection, RJ45 socket.

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***Connection "RS232 A"**

ESPA 4.4.4 protocol.

For example, connection of a radio paging system, DECT system, fire alarm system or medical electrical system.

NOTE! If room terminals are connected to the OSYnet port and an RS232 port is used a separator device according to EN 60601-1 (2 x MOPP) must be connected between the connected device and the RS232 port, e.g. the interface isolator RS232 (order no. 76 5010 00).

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
* The connections marked with a star are only enabled, if the IP-SystemManager was prepared accordingly at the factory, refer to the project specific documents.

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EN Connecting

1. Strip the connection cables to a suitable length. Remove the cable shield and drain wire up to the cable sheath. Strip off the wires 7 mm.

CAUTION! If the IP-SystemManager is operated without insertion bridges at the power connections, high operating currents from the system flow through the IP-SystemManager. This may damage the device. The two insertion bridges **must** be inserted according to fig. E!

2. Connect the wires to the connectors (pluggable screw-type connectors) according to fig. D.
3. Insert the patch cable's RJ45 plug into the Ethernet socket.

LED displays**Green LED "Power"**

LED lights up permanently. = Power supply available at the IP-SystemManager.

Green LED "Status"

LED flashes once a second. = IP-SystemManager is ok.

LED lights up permanently or is off. = Fault! IP-SystemManager is not operational.

Red LED "Fault" below "Fault Relay"

LED lights up permanently. = There is at least one fault.

Signalling in parallel to "Fault Relay" output.

This is a collective indicator. The assigned fault locations depend on the operating mode (see device label) of the IP-SystemManager:

- Operating mode "System" or "System + Local"

Collective indicator for faults in the complete nurse call system. Indicated fault types:

- An OSYnet bus user is defective or disconnected.
- A RAN bus user is defective or disconnected.
- The power supply on a ward is interrupted. This fault is recognized indirectly by a failure in OSYnet bus user polling.
- The power supply at an IP-SystemManager is interrupted. This fault is recognized indirectly by a failure in OSYnet bus user polling.
- An OSYnet bus line is interrupted.
- LAN connection to an IP-SystemManager is interrupted.
- An IP-SystemManager is not operational.
- The power supply unit is in buffer mode. Indication in parallel to the red LED "Fault" below "UPS Input".

- Operating mode "Local"

Collective indicator for faults in the ward connected to this IP-SystemManager. Indicated fault types:

- An OSYnet bus user is defective or disconnected.
- A RAN bus user is defective or disconnected.
- The power supply unit is in buffer mode. Indication in parallel to the red LED "Fault" below "UPS Input".

NOTE! The red LED "Fault" (below "Fault Relay") and the "Fault Relay" connection of IP-SystemManagers with operating mode "BMA/MED", "Voice Gateway" and "Voice Gateway Master" are inactive.

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EN - Installation Instructions

Red LED "Fault" below "UPS Input"

LED lights up permanently. = The uninterruptible power supply operates in buffer mode.

OSYnet 1: Green LED "Rx", yellow LED "Tx"

LEDs are flashing. = Data traffic on the group bus, connected to "OSYnet 1".

Yellow LED lights up permanently. = Group bus OSYnet 1 has not yet been configured with the SystemOrganizer.

OSYnet 2: Green LED "Rx", yellow LED "Tx"

Yellow LED lights up permanently. = There is no active use licence for OSYnet 2.

Yellow LED "NF 1"

LED lights up permanently. = There is a speech connection or an announcement in the group connected to OSYnet 1.

Yellow LED "NF 2"

Not used.

Ethernet, green LED "Link", yellow LED "Speed"

The LEDs are signalling "Link" and "Speed" of the Ethernet connection.

RS232 A, green LED "Rx", yellow LED "Tx"

LEDs are flashing. = Data traffic between the IP-SystemManager and the device connected to „RS232 A“.

RS232 B, green LED "Rx", yellow LED "Tx"

As for "RS232 A", but for RS232 B.

F Service cover


The service cover [2] is removed by levering off. After use re-attach the service cover.

5 - 6 grommets for fixing the service cover

6 - Li battery, 3 V, CR2032

7 - Reset button

8 - DIP switch for addressing: Do not change the factory setting!

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Replace battery

Behind the service cover [2] there is a lithium button cell [6] (3 V, CR2032) to buffer the internal clock in case of power failure.

In normal operation of the nurse call system, i.e. the IP-SystemManager is supplied with power, the battery should be replaced every 2 years. If the IP-SystemManager is not supplied with power, e.g. when the device is stored, an earlier replacement may be necessary due to the battery consumption. Replacement battery: High-quality Li battery, 3 V, CR2032

Behaviour of the system when the battery is empty

An empty battery leads to the following situation: After a failure of the supply voltage of the IP-SystemManager, the internal clock restarts at 00:00. An incorrect time is displayed on the display devices of the nurse call system (ComStation, corridor display, ComTerminals, RoomTerminals). The time is correct again when a connected TimeServer sends the time. If no TimeServer is connected, the time must be reset manually via the SystemOrganizer.

Battery disposal

The used battery must be returned to a collection point for used batteries or to Tunstall.


Technical data

Dimensions (HxWxD):	132 x 216 x 48 mm
Housing material:	PS UL94-HB
Service cover material:	ABS
Weight:	415 g
IP protection class:	IP 20
Nominal voltage:	24 V DC
Permitted voltage range:	20 – 26 V DC
Ambient temperature:	+5 °C – +40 °C
Relative humidity:	0% – 85% (non condensing)
Standby current consumption:	120 mA

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RAN-Multiplexer Flamenco, Best.-Nr. 77 0270 00

RAN multiplexer Flamenco, order no. 77 0270 00

Multiplexer controls the nurse call functions of up to 5 rooms. Connection to the group bus (OSYnet), the room bus (RAN), and the 24 V DC power supply of the nurse call system via connection socket ComTerminal.

The connection socket ComTerminal (order no. 70 0491 00) is not included in the scope of delivery of the RAN multiplexer.

NOTE! The basic installation of the nurse call system is described in the Technical Manual.

1. Product characteristics

- Max. 5 rooms can be connected
- Max. 50 room users (switches, connection sockets, room lamps, etc.) can be connected, distributed as needed among the 5 rooms and the collective display (as of software revision 1.01)
- A room lamp can be used as a collective display (display of messages from all connected rooms) (as of software revision 1.01)
- Total max. 100 m room bus cable (RAN)
- A failure of the collective display is indicated in the nurse call system as a fault of room 1.

2. Technical data

Dimensions (HxWxD):	205 x 110 x 34 mm
Power supply:	24 V DC
Standby current consumption without connected RAN users:	38 mA
Housing material:	ABS
Ambient temperature:	+5 °C – +40 °C
Relative humidity:	0% – 85% (non condensing)

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3. Preparatory work for renovations

NOTE! This chapter 3 does not apply to a new installation of a RAN multiplexer Flamenco.

The RAN multiplexer Flamenco, order no. 77 0270 00, makes it possible to switch over a room bus (RAN), which was connected to a CCS bus via a RAN multiplexer, order no. 70 0270 00, to a group bus (OSYnet).

The installed connection socket ComTerminal, via which the RAN multiplexer (70...) was connected, can still be used for the new RAN multiplexer Flamenco.

To prepare for the update you need to:

- write down the address of the old RAN multiplexer (70...)
- count the number of connected RAN users in each room.

3.1 Dismantling the old RAN multiplexer (70...) > Fig. A

1. To remove the lower cover [6] of the RAN multiplexer [5], press a screw driver into a lateral recess and lever out the cover [6].
2. Write down the RAN multiplexer's address, coded by the 10-pole DIP switch [7], as described in the chapter 3.2 and shown in Fig. B.
3. Remove the RAN multiplexer [5] from the mounting plate [3] by releasing it at the snap lug and simultaneously pushing it upwards.
4. Disconnect the flat cable from the connection socket [2].
5. Unscrew the four fixing screws [4] until the mounting plate [3] comes loose. Remove the mounting plate [3]. The screws [4] remain in the connection socket [2].

Do not dismantle the connection socket [2]. It will be reused with the new RAN multiplexer Flamenco.


3.2 Address of the old RAN multiplexer (70...) > Fig. B

The RAN multiplexer's address is coded by DIP switches 1 to 6 of the 10-pole DIP switch [7].


DIP switch 1 codes address 1. DIP switch 2 codes address 2, DIP switch 3 codes address 4 etc, see Fig. B. The appropriate DIP switches are set to „ON“. Other addresses are coded by addition, e.g. DIP switch 4 + DIP switch 5 = Address 24.

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4. Connecting the connection socket > Fig. C

1. Connect the wires of the connection cables and the RAN expander [9] to the connection socket [2] according to Fig. C.
2. Insert fine wire fuse 2 A/F 5x20 [8] into the fuse holder of the connection socket [2].

ATTENTION, fuse! The connection socket (70 0491 00) is delivered with a fuse 1 A/F. For operation with the RAN multiplexer this fuse must be replaced by a fine wire fuse 2 A/F 5x20 (00 0130 36).


ATTENTION, terminating resistor! In the last device on the group bus (OSYnet) an 120 ohms resistor (00 0040 76) must be fitted between the connection points DH and DL.

Handling of terminals with solderless connection technique

Prerequisite: The wires must be stripped to 8 mm.

- Connect a wire: Insert the stripped wire into the terminal until it stops. To make insertion easier put a small screw driver into the opening (see Fig. C) while inserting the wire.
- Disconnect a wire: For loosen the wire put a small screw driver into the opening (see Fig. C) and pull out the wire.

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

5. Mounting > Fig. D

1. Connect the cables to the ComTerminal connection as described in chapter 4.
2. Insert fine wire fuse 2 A/F 5x20 [8] into the fuse holder of the connection socket [2] as described in chapter 4.
3. Loosen the four back box screws [4].
4. Insert the flat cable [10] through the outlet in the housing rear panel [11] and plug it into the appropriate socket on the connection socket [2].
5. Fix the housing rear panel [11] using the back box screws [4] to the back box [1].
6. Plug the flat cable [13] of the housing front panel [12] into the appropriate socket on the housing rear panel [11].
7. Place the housing front panel [12] onto the top of the housing rear panel [11] creating a hinge. Now press the housing front panel [12] downwards onto the rear panel [11] with low pressing power until the catch [14] engages. Make sure that the flat cable [13] is not pinched.
8. For improved security a securing screw [15] can be fitted into the existing hole.

Finally, the RAN multiplexer Flamenco has to be configured, see „Configuration Instructions“ as of page 6.

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- | | |
|--|--|
| 1 - * Back box | 10 - Flat cable of the housing rear panel |
| 2 - * Connection socket ComTerminal (70 0491 00) | 11 - Housing rear panel |
| 4 - * Back box screws | 12 - Housing front panel |
| 8 - Fuse 2 A/F (00 0130 36) | 13 - Flat cable of the housing front panel |
| 9 - RAN expander (70 0908 10) | 14 - Catch |
| | 15 - Securing screw |

* Not included in the scope of delivery of the RAN multiplexer Flamenco.

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6. Room addresses and channel numbers

A maximum of 5 rooms can be connected to one RAN multiplexer. The RAN multiplexer converts the channel numbers of its connected room devices (call switches, connection sockets etc.) to room addresses, see Fig. E.

The room address of the first connected room must be set at the RAN multiplexer. The four following addresses, i.e. the addresses of the other rooms, will be assigned automatically.

The room addresses are always assigned in steps of five, even if only three rooms are connected to the RAN multiplexer. The unused room addresses remain blank.

When renovating, set the same values for room addresses and channel numbers that were set in the old installation.

7. Set channel numbers at room devices

Room devices = Switches, connection sockets, room lamps etc. **No** pear push switch!

- Set the channel number at each room device.

Every device in a room must have the same channel number.

The five rooms connected to a RAN multiplexer are numbered from 1 to 5, i.e. channel number 1 to 5.

NOTE! At connection sockets the channel number equals to the bed number.


Refer to the installation instructions delivered with the room device to learn how to set its channel number or bed number.

Set channel number for a collective display

It is possible to use a room lamp as a collective display, i.e. it signals the calls and presences of all 5 rooms connected to the RAN multiplexer Flamenco.

- To set a room lamp as a collective display, the channel number 0 must be set.

NOTE! Channel 0 must not be set on any other devices (switches, connection sockets, etc.).

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
Symbols in the configuration menu

✓	Checkmark	Set the marked menu item.
▲	Upwards arrow	Scroll upwards through the list.
▼	Downwards arrow	Scroll downwards through the list.
⌂	Home	Cancel, without storing.

8.2 Make settings

8.2.1 Setting of configuration menu language

1. Mark „Language“ using the arrow keys. Then press the check-mark.
2. Mark the desired language using the arrow keys:
D = German; GB = English etc.
3. Press the check-mark to set the marked language.

8.2.2 Setting of room address

Here you set the same address that was set at old RAN multiplexer (70...) at this location, see chapter 3.2. For general explanation of the room address, see chapter 6.

1. Mark „Address“ using the arrow keys. Then press the check-mark.
2. Mark the desired room address using the arrow keys.
3. Press the check-mark to set the marked address.

8.2.3 Setting of RAN numbers (Room Area Network)

Here you set number of room devices for each connected room. Room devices = Switches, connection sockets, room lamps etc. **No** pear push switch!

1. Mark „RAN“ using the arrow keys. Then press the check-mark.
2. Mark „RAN number“ using the arrow keys. Then press the check-mark.
3. Mark „Room 1“ using the arrow keys. Then press the check-mark.
4. Mark the RAN number for room 1 using the arrow keys.
5. Press the check-mark to set the marked RAN number for room 1.
6. Carry out step 3. to 5. for all connected rooms.

NOTE! Usually, one room lamp corresponds to one RAN user. In the former Tunstall product range, there were room lamps that unusually corresponded to 2 RAN users. This is noted in the documentation for the respective room lamps.

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8.2.4 Setting RAN number for collective display

If a room lamp is set to channel number 0, i.e. it is used as a collective display, you must set the RAN number for the collective display:

1. Mark „RAN“ using the arrow keys. Then press the check-mark.
2. Mark „RAN number“ using the arrow keys. Then press the check-mark.
3. Mark „Collective display“ using the arrow keys. Then press the check-mark.
4. Mark the correct value using the arrow keys:
0 = No collective display installed.
1 = Collective display installed. The room lamp corresponds to one RAN user.
2 = Collective display installed. The room lamp corresponds to two RAN users.

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5. Press the check-mark to set the marked value.

8.3 Displaying the software revision

If you want to know the RAN multiplexer's software revision:

1. Mark „Info“ using the arrow keys. Then press the check-mark.

The RAN multiplexer's software revision is displayed.

2. Press the home key to return to the menu.


8.4 Checking the room bus RAN

8.4.1 Status

Check whether the set RAN numbers for room devices and collective display equal to the numbers of operational devices, see “8.2.3 Setting of RAN numbers” and “8.2.4 Setting RAN number for collective display”.

1. Mark „RAN“ using the arrow keys. Then press the check-mark.
2. Mark „Status“ using the arrow keys. Then press the check-mark.

NOTE! The Test can take some time.

- (a) Sum of all operational room devices
- (b) Sum of all set room devices
- (c) C1: Collective display (1 RAN user) operational
C2: Collective display (2 RAN users) operational
- (d) Operational room devices per room
Room 1 | Room 2 | Room 3 | Room 4
- (e) „OK“ : All set RAN numbers (room devices, collective display) match to the numbers of operational devices.

„Error“: **Not** all set RAN numbers (room devices, collective display) match to the numbers of operational devices! Fix error.

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- Press the home key to end. In case of an **error** perform a „Test RAN“.

8.4.2 Test RAN

Test if all devices (room devices, collective display) are operational and if they are correctly connected to the RAN.

1. Mark „RAN“ using the arrow keys. Then press the check-mark.
2. Test start: Mark „Test RAN“ using the arrow keys. Then press the check-mark.

HINWEIS! The Test can take some time.

The LEDs of all connected devices must flash. Non-flashing devices are wrongly connected or defective. Replace defective devices. Correct any wrong connection.

3. Compare the number of flashing room devices with the set RAN number for room devices. If the set RAN number is wrong, correct the RAN number setting after the RAN test as described on the previous page.
4. End of test: Press the home key.

8.5 Exit the configuration menu

When all settings are made and all tests are completed, exit the configuration menu:

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9. Configuring the nurse call system (SystemOrganizer)

The complete nurse call system must be configured using the PC software „SystemOrganizer“. How to configure a nurse call system with the SystemOrganizer is described in the online help for the SystemOrganizer.

For configuring the RAN multiplexer Flamenco within the SystemOrganizer you need to know the following additional information:

For each RAN multiplexer Flamenco 5 locations must be created. For each of these 5 locations the following properties must be set:

- Location Type: Patient Room
- Installed Device: RoomTerminal
- The 5 locations must be in the same ward and have 5 consecutive addresses.

NOTE! Even if less than 5 rooms are connected to the RAN multiplexer Flamenco, 5 locations must be created in the SystemOrganizer.

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

Installation location requirements

- Limited access location for technical equipment
- Dry location (max. relative humidity 75% at approx. 18 °C)
- Ambient temperatures between 0° and 40 °C
- The OSY-ControlCenter shall be readily accessible (access aisle 60 cm wide).
- Consider sufficient free space for cable entries.



Ventilation

The heat transfer must not be blocked. The selected location must ensure sufficient air circulation above and below the OSY-ControlCenter. Therefore, check for minimum distance of 50 mm (3 U in a rack) to walls, ceiling and other equipment. The venting ports must not be covered up. Where the OSY-ControlCenter is installed in a 19-inch rack, forced-air ventilation may be required.

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

DANGER! Potential equalisation required.

Power supply unit 5A, DIN rail (77 3410 50)
+ UPS control unit 10A (77 3411 00)
+ Battery module (77 3412 00)

0V

+24V

Ferrite

(order no. 00 0080 34)


2 x 2.5 mm²

0V

+24V

10 cm

① Connection for potential equalisation on the threaded bolt in the housing:

Internal potential equalisation connection

② Fault message output


- Potential-free output.
- Collective output.
- Contact load: max. 100 mA / 24 V.
- In case of a fault contact opens at the screw clamps.

Fault message output switches in case of:

- Disconnected or defective OSYnet bus user.
- Disconnected or defective RAN bus user.
- Interrupted power supply on the ward.
- Interrupted power supply at the OSY-ControlCenter.
- Interrupted OSYnet bus line.
- OSY-ControlCenter does not boot-up.

③ Assignment for the 9 pole Sub-D connector:

- RxD = 2
- Tx/D = 3
- GND = 5


Important notes

- Check all cables for continuity, short circuits and ground faults.
- When the data cables are correctly connected, the values are: voltage between DH and DG and voltage between DL and DG approx. 2.5 V.
- **ATTENTION!** The connection plan on the following page is a standard connection plan. Please follow the deviations in your project specific connection plan.

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Connection of network switch

Depending on its components the OSY-ControlCenter may include a maximum of 6 RJ45 sockets. Each socket has to be connected to an arbitrary RJ45 socket in the network switch. But the following clear arrangement is recommended.

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Battery replacement

WARNING! Improperly performed work on the control components of the OSY-ControlCenter can lead to failures and malfunctions of the nurse call system. Therefore, work on the control components may only be carried out by Tunstall engineers.

The nurse call system is controlled by embedded PCs mounted on plug-in cards. For each embedded PC, there is a lithium button cell (3 V, CR2032) on the plug-in card to buffer the internal clock in case of a power failure.

All batteries of the OSY-ControlCenter should be replaced every 2 years. Have the batteries replaced by a Tunstall engineer!

Replacement battery: High-quality Li battery, 3 V, CR2032

Behaviour of the system when a battery is empty

The empty battery of an embedded PC leads to the following situation: After a failure of the supply voltage of the OSY-ControlCenter, the internal clock for the embedded PC restarts at 00:00. An incorrect time is displayed on the display devices of the ward (ComStation, corridor display, ComTerminals, RoomTerminals). The time is correct again when a connected TimeServer sends the time. If no TimeServer is connected, the time must be reset manually via the SystemOrganizer.

Disposal instructions


Used batteries must not be disposed of in household waste. Used batteries must be returned to a collection point for used batteries or to Tunstall.

The symbol of the crossed-out waste bin indicates that this electronic device must not be disposed of in household waste. It must be returned for reuse to a collection point for waste electrical equipment or to Tunstall. If the old device contains personal data, you as the end user are responsible for deleting this data before returning it. If possible, remove the batteries from the old device before returning it for disposal.


The control LEDs help Tunstall engineers to analyse faults.

- ① Linie / Line
- ② Durchsageverstärker / Announcement amplifier
- ③ Klingelsignal / Ringing signal
- ④ Schleifenspannung / Loop current

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LED 1 is ON: Speech connection between the telephone line and nurse call system.
LED 2: is always OFF.
LED 3 flickering: Ringing signal on the telephone line.
LED 4 is ON: Telephone off-hook.

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OSYnet-SafetyController, order no. 77 2100 00

Safety controller for direct connection to the group bus (OSYnet). Usable in parallel with OSY-ControlCenter, per physical group. Intended use as redundancy for maintaining the system operation in case of OSY-ControlCenter fail.

- Controlling the base light call functions of one physical group (DIN VDE 0834)
- Control of room lamps and displays
- Autonomous and automatic switching between normal and safety operation
- 1 potential free digital output (change-over contact) for signalling the operation mode
- 1 potential free digital output (change-over contact) for fault signalling
- Digital output: internal: 24 VDC/60 mA, external: 24 VAC/1 A
- Max. current consumption: 230 mA

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A Mounting

Mounting on supporting rail or wall mounting in the electrical room of the ward.

A1: Mounting on supporting rail (35 mm)

Click the OSYnet-SafetyController 1 onto the supporting rail 3.

A2: Wall mounting

1. Push out the three fixing clips 2 until the holes for the screws 4 appear.
2. Screw the OSYnet-SafetyController 1 to the wall.

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At the OSYnet-SafetyController you have to set **address 126**. If you install a second OSYnet-SafetyController in the same physical group, you have to set **address 125** at the second OSYnet-SafetyController.

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Reset

To reset the OSYnet-SafetyController you have to set code switch 10 to ON for one second, then back to OFF. A reset has to be made, if a code switch setting has been changed. Otherwise the new setting will not be activated.

B Control LED (red)


Light OFF:	OSYnet-SafetyController without power supply or defect
Flashing slowly: (1 s ON - 1 s OFF)	OSY-ControlCenter controls the group. OSYnet-SafetyController ready-to-operate (standby)
Flashing quickly: (200 ms ON - 200 ms OFF)	OSYnet-SafetyController controls the group.
Light ON for more than 10 s:	OSYnet-SafetyController not ready to operate.


C ESD protection cover (order no. 00 0276 53)

After you have connected the cables (see reverse side), you have to put the ESD protection cover on the OSYnet-SafetyController. (ESD = Electro Static Discharge)

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Put a screw driver into the visible hole of the black fixing clip. Then pull the clip upwards, until the OSYnet-SafetyController 1 comes loose from the supporting rail 3.

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Output for signalling the operation mode

OUT 1 closed & OUT 2 open = Normal operation

OUT 1 open & OUT 2 closed = Safety operation (SafetyController controls the group)

Jumper OUT 1/2 set = Internal power supply (24V DC, maximum 60 mA for triggering a relay)

Jumper OUT 1/2 not set = External power supply (24V AC or DC, maximum 1 A)

Output for fault signalling

Faults of group bus users

OUT 3 closed & OUT 4 open = no fault

OUT 3 open & OUT 4 closed = Fault!

Jumper OUT 3/4 set = Internal power supply (24V DC, maximum 60 mA for triggering a relay)

Jumper OUT 3/4 not set = External power supply (24V AC or DC, maximum 1 A)

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
Power supply unit 10A, DIN rail, Order no. 77 3410 00

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
WARNING! Risk of fire:

The power supply unit is suitable for Flamenco, Flamenco^{SE}, EccoLine, NewLine, CONCENTO and VKS 3000 systems.

The power supply unit must not be used in legacy systems CCS (nurse call system) and VKS 2000 (prison system). In these systems the cables may burn in case of a short-circuit.

If this power supply needs to be used in CCS or in VKS 2000 system, each room terminal has to be equipped with a fuse as shown below:

Connections
in the terminal

① Plug-in screw terminal with fuseholder (order no. 00 0224 81), fuse 3.15 A, slow blow (order no. 00 0130 23)

② Push-wire connector for junction boxes (order no. 00 0210 21)

NOTE! VKS 2000 terminals are identified by the type no. 22.96xx or 22.98xx on the printed circuit board.

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Product Description

The 77 3410 00 is a DIN-rail mountable single-phase-input power supply for industrial use and for use in medical equipment. It provides a floating, stabilized and galvanically separated SELV/PELV output voltage.

The specialty of this power supply is the 2MOPP (two means of patient protection) safety approval for medical use and the fulfilment of the required EMC tests for professional healthcare facility and home healthcare environments.

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

Intended Use

This device is designed for installation in an enclosure and is intended for use in medical equipment.


Do not use this device in HF surgical equipment or similar equipment. In systems where failure or failure of this equipment may cause serious injury or endanger human life, this equipment must be used only if technical and / or organizational measures are taken to correct the effects of the malfunction.

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Installation Instructions

Install device in an enclosure providing protection against electrical, mechanical and fire hazards.

Install the device onto a DIN-rail according to EN 60715 with the input terminals on the bottom of the device. Other mounting orientations require a reduction in output current.

Use shielded or unshielded cables, twisted or untwisted cables of any length. Make sure that the wiring is correct by following all local and national codes. Use appropriate copper cables that are designed for a minimum operating temperature of 90°C for ambient temperatures between +45°C and +60°C and 105°C for ambient temperatures up to +70°C. Ensure that all strands of a stranded wire enter the terminal connection.

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The enclosure of the device provides a degree of protection of IP20. The enclosure does not provide protection against spilled liquids.

The isolation of the device is designed to withstand impulse voltages of overvoltage category III according to IEC 60664-1.

The device is designed as "Class of Protection I" equipment according to IEC 61140. Do not use without a proper PE (Protective Earth) connection.

The device is suitable to be supplied from TN, TT or IT mains networks. The continuous voltage between the input terminal and the PE potential must not exceed 300Vac.

The input can also be powered from batteries or similar DC sources. The continuous voltage between the supply voltage and the PE/ground potential must not exceed 375Vdc.

A disconnecting means shall be provided for the input of the device.

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The device is designed for altitudes up to 3000m (9840ft).

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The device is designed, tested and approved for branch circuits up to 32A (IEC) and 30A (UL) without additional protection device. If an external fuse is utilized, do not use circuit breakers smaller than 6A B- or C-Characteristic to avoid a nuisance tripping of the circuit breaker.

The maximum surrounding air temperature is +70°C (+158°F). The operational temperature is the same as the ambient or surrounding air temperature and is defined 2cm below the device.

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Portable RF communication equipment should be used no closer than 30cm (12 inch) to the power supply or the power supply connection cables.

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

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Functional Description

The output is electronically protected against no-load, overload and short circuit and can supply any kind of loads, including unlimited inductive and capacitive loads. If capacitors with a capacitance >1.5F are connected, the unit might charge the capacitor in an intermittent mode.

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The output voltage can be adjusted with a small flat-blade screwdriver on the front of the unit.

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The device is equipped with an over-temperature protection. In case of a high temperature, the output shuts down and starts automatically again after cooling off.

At heavy overloads (when output voltage falls below 13V), the device delivers continuous output current for 2s. After this, the output is switched off for 18s before a new start attempt is automatically performed. This cycle is repeated as long as the overload exists.

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In case of an internal defect, a redundant circuit limits the maximum output voltage to 32V. The output shuts down and automatically attempts to restart.

Technical data

All values are typical figures specified at 230Vac 50Hz input voltage, 24V 10A output load, 25°C ambient temperature and after a 5 minutes run-in time unless otherwise noted.

Output voltage	DC 24V	Nominal
Adjustment range	24 – 28Vdc	Factory setting 24.1V
Output current continuous	12.0 – 10.3A	Below +45°C ambient
	10.0 – 8.6A	At +60°C ambient
	7.5 – 6.5A	At +70°C ambient
	Derate linearly between +45°C and +70°C	
Input voltage AC	AC 100 – 240V	-15% / +10%
Mains frequency	50 – 60Hz	±6%
Input current AC	2.15 / 1.13A	At 120 / 230Vac
Power factor	0.99 / 0.97	At 120 / 230Vac
Input voltage DC	DC 110 – 150V	±20%
Input current DC	2.35A	At 110Vdc
Input inrush current	6 / 9A pk	At 40°C, 120 / 230Vac
Efficiency	93.6 / 95.2%	At 120 / 230Vac
Power losses	16.4 / 12.1W	At 120 / 230Vac
Hold-up time	37ms	
Temperature range	-25 to +70°C	
Max. wire size (lit wire)	4mm ²	Power terminals
Wire size AWG	AWG 20-10	Power terminals
Max. wire diameter	2.8mm	Power terminals
Wire stripping length	10mm / 0.4inch	Power terminals
Max. wire size (lit wire)	1.5mm ²	DC-OK terminals
Wire size AWG	AWG 24–16	DC-OK terminals
Max. wire diameter	1.6 mm	DC-OK terminals
Wire stripping length	7mm / 0.28inch	DC-OK terminals
Size (wxhxd)	39x124x117mm	Without DIN-rail
Weight	620g / 1.37lb	

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
Power supply unit 5A, DIN rail, Order no. 77 3410 50

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
WARNING! Risk of fire:

The power supply unit is suitable for Flamenco, Flamenco^{SE}, EccoLine, NewLine, CONCENTO, and VKS 3000 systems.

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Connections
in the terminal

- ① Plug-in screw terminal with fuseholder (order no. 00 0224 81), fuse 3.15 A, slow blow (order no. 00 0130 23)
- ② Push-wire connector for junction boxes (order no. 00 0210 21)

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Product Description

The 77 3410 50 is a DIN-rail mountable single-phase-input power supply for industrial use and for use in medical equipment. It provides a floating, stabilized and galvanically separated SELV/PELV output voltage. The specialty of this power supply is the 2MOPP (two means of patient protection) safety approval for medical use and the fulfilment of the required EMC tests for professional healthcare facility and home healthcare environments.

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


Do not use this device in HF surgical equipment or similar apertures. In systems where failure or failure of this equipment may cause serious injury or endanger human life, this equipment must be used only if technical and / or organizational measures are taken to correct the effects of the malfunction.

If this device is used in a manner outside of its specification, the protection provided by the device may be impaired.

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
Use shielded or unshielded cables, twisted or untwisted cables of any length. Make sure that the wiring is correct by following all local and national codes. Use appropriate copper cables that are designed for a minimum operating temperature of 90°C. Use ferrules for wires on the input terminals.

Ensure that all strands of a stranded wire enter the terminal connection.

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The enclosure of the device provides a degree of protection of IP20. The housing does not provide protection against spilled liquids.

The device is designed for overvoltage category II zones.

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Please note: This device incorporates a single fuse in the L-input. For medical use an additional external fuse in the N-input may be necessary, depending on the type of medical equipment.

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Functional Description

The output is electronically protected against no-load, overload and short circuit and can supply any kind of loads, including inductive and capacitive loads. If capacitors with a capacitance >0.3F are connected, the unit might charge the capacitor in an intermittent mode (Hiccup behavior).

Do not apply return voltages from the load to the output higher than 35V.

The output voltage can be adjusted with a small flat-blade screwdriver on the front.

The green DC OK LED reports an output voltage above 90% of the adjusted voltage of a running device.

The DC OK relay monitors the output voltage and the contact is closed when the DC OK LED is on. Contact ratings: 60Vdc 0.3A, 30Vdc 1A, 30Vac 0.5A for resistive loads.

The device is equipped with an over-temperature protection. In case of a high temperature, the output shuts down and starts automatically after cooling off.

At heavy overloads (when output voltage falls below 13V), the device delivers continuous output current for 1s. After this, the output is switched off for approx. 9s before a new start attempt is automatically performed. This cycle is repeated as long as the overload exists.

Do not connect outputs of devices in a parallel connection for higher output currents.

Do not connect outputs of devices in a series connection for higher output voltages.

In case of an internal defect, a redundant circuit limits the maximum output voltage to 32V. The output switches off and performs three re-start attempts. If the failure continues, the output shuts down. Cycle input power to reset.

Technical data

All values are typical figures specified at 230Vac 50Hz input voltage, 24V 5A output load, 25°C ambient temperature and after a 5 minutes run-in time unless otherwise noted.

Output voltage	DC 24V	Nominal
Adjustment range	24 – 28Vdc	Factory setting 24.1V
Output current	5.0 – 4.3A	
Input voltage AC	AC 100 – 240V	-15% / +10%
Mains frequency	50 – 60Hz	±6%
Input current AC	1.09 / 0.6A	At 120 / 230Vac
Power factor	0.98 / 0.91	At 120 / 230Vac
Input inrush current	5 / 6A pk	At 120 / 230Vac, 40°C, cold start
Efficiency	93.6 / 94.3%	At 120 / 230Vac
Power losses	8.2 / 7.3W	At 120 / 230Vac
Hold-up time	35 / 35ms	At 120 / 230Vac
Temperature range	-25 to +60°C	
Max. wire size (litz wire)	4mm ²	Power terminals
Wire size AWG	AWG 20-10	Power terminals
Max. wire diameter	2.8mm	Power terminals
Wire stripping length	10mm / 0.4inch	Power terminals
Max. wire size (litz wire)	1.5mm ²	DC-OK terminals
Wire size AWG	AWG 24–16	DC-OK terminals
Max. wire diameter	1.6 mm	DC-OK terminals
Wire stripping length	7mm / 0.28inch	DC-OK terminals
Size (wxhxd)	32x124x102mm	Without DIN-rail
Weight	440g / 0.97lb	

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
Read this first!

Before operating this device, please read this manual thoroughly and retain this manual for future reference! This device may only be installed and put into operation by qualified personnel. If damage or malfunction should occur during operation, immediately turn power off and send device to the factory for inspection. The device does not contain serviceable parts. The information presented in this document is believed to be accurate and reliable and may change without notice. For any clarifications the English translation will be used.

WARNING! Risk of electrical shock, fire, personal injury, or death:

- Turn power off before working on the device. Protect against inadvertent re-powering.
- Do not open, modify or repair the device.
- Use caution to prevent any foreign objects from entering the housing.
- Do not use in wet locations or in areas where moisture or condensation can be expected.
- Do not touch during power-on and immediately after power-off. Hot surfaces may cause burns.


Product Description

The 77341100 is an uninterruptible power supply controller (DC-UPS), which is used in combination with a 24V power supply and one external 12V-battery to bridge power failures. When the power supply provides sufficient voltages, the DC-UPS controller charges the battery. When the power supply voltage fails, the energy stored in the battery is re-leased to the DC bus in a regulated process.

Intended Use

This device is designed for installation in an enclosure and is intended for commercial use, such as in industrial control, process control, monitoring and measurement equipment or the like. In systems where failure or failure of this equipment may cause serious injury or endanger human life, this equipment must be used only if technical and / or organizational measures are taken to correct the effects of the malfunction.

Functional Diagram

Important note for nurse call systems

The German standard DIN VDE 0834 for nurse call systems requires compliance with EN 60601-1 regarding electrical safety. This means: The 24V power supply unit must be earthed. All other devices of the modular power supply, e.g. the UPS control unit 10 A (77341100), must be mounted insulated against protective earth (PE). Recommendation for insulated mounting: Install the UPS control unit on a DIN rail raiser HAGER HYA036 (order no. 76090004, width 106mm). In addition, a distance of 1cm must be maintained from the insulated mounted device to all earthed surfaces and equipment.

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Installation Instructions

Turn power off and disconnect the battery fuse before working on the device.

Install the device in an enclosure providing protection against electrical, mechanical and fire hazards.

Install the device onto a DIN-rail according to EN 60715 with the input terminals on the top of the unit.

The input can be powered from a regulated power supply or a similar DC source. The voltage between the input and ground must not exceed 60Vdc continuously. (WARNING! For nurse call systems, the protective earth (PE) must not be connected to the UPS control unit.) The ripple voltage in the low frequency range between 50Hz and 10kHz must be negligible when used in marine applications.

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Use an appropriately sized 24V power supply, which can deliver the additional internal current consumption, required to charge the battery.

Use only VRLA lead acid batteries with a capacity between 3.9Ah and 40Ah.

Check for correct input and battery polarity. The device will not operate when the voltage is reversed.

Make sure that the wiring is correct by following all local and national codes. Use appropriate copper cables that are designed for a minimum operating temperature of 60°C for ambient temperatures up to +45°C, 75°C for ambient temperatures up to +60°C and 90°C for ambient temperatures up to +70°C. Ensure that all strands of a stranded wire enter the terminal connection.

Do not use wires smaller than 2.5mm² (or AWG 12) and not longer than 2x1.5m between the battery and the DC-UPS controller. Longer or smaller gauge wires can change performance of the system.

Use a 30A fuse (ATOF® 287 030 from Littelfuse or an UL listed fuse with the same characteristics) in the battery circuit. The battery fuse protects the wires between the battery and the DC-UPS and shall be located close to the battery.

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A disconnecting means shall be provided for the input and the battery input of the device.

The device is designed for convection cooling and does not require an external fan. Do not obstruct airflow and do not cover ventilation grid! Keep the following minimum installation clearances: 40mm on top, 20mm on the bottom, 5mm left and right side. Increase the 5mm to 15mm in case the adjacent device is a heat source. When the device is permanently loaded with less than 50%, the 5mm can be reduced to zero.

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Green Status LED

This LED is on solid, when the battery is charged more than 85%. During battery charging, the LED is flashing with a slow frequency and rapid flashing is an indication of a buffer event.

Yellow Diagnosis LED

This LED reports an overload situation, an expired buffer time (flashes for 15 minutes), an activated inhibit input or when the battery should be replaced.

Red Check Wiring LED

This LED indicates a failure in the installation such as too low input or battery voltage, incorrect wiring, blown or missing battery fuse.

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Buffer-time Limiter

The buffer time during a buffer event can be set to 10s, 30s, 1min, 3min, 10min or infinity with the dial on the front. This can save battery energy for faster recharge and expand the lifetime of the battery. When dial is set to infinity, buffering continues until stopped by the battery deep-discharge protection.

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EN - Installation Instructions

End-of-charge Voltage

The end-of-charge-voltage shall be set according to the expected temperature in which the battery is located. The dial on the front allows a setting between +10 and +40°C. If in doubt, set the unit to 40°C.

Inhibit Input

The inhibit input disables or stops buffering. In power supply mode, a continuous input of 10-30V is required. During buffering, a short input of 10-30V for at least 250ms is required to stop buffering. The current of the inhibit input is limited to 6mA by the DC-UPS control unit.

Ready, Buffering and Replace Battery Relay Contacts

The Ready contact is closed when the battery is 85% charged and the unit is ready to buffer.

The Buffering contact is closed when the unit is in buffer mode.

The Replace Battery contact is closed when a battery replacement is necessary.


Contact ratings: 60Vdc 0.3A, 30Vdc 1A, 30Vac 0.5A for resistive loads.

Technical data

All values are typical figures specified at 24Vdc input voltage, 10A output current in power supply mode at 25°C ambient, no charging and after a 5 minutes run-in time unless otherwise noted.

Input voltage	DC 24V	-20%/+25%
Minimum input voltage to start charging and to enable battery mode	22.8Vdc	
Transfer voltage to switch into battery mode	22.5Vdc	
Internal current consumption	1.1A	Includes charging current
Output in power supply mode:		
Input to output voltage loss	230mV	At 10A output load
	330mV	At 15A output load
Output current	15A	Below +60°C
	11.3A	At +70°C
	Derate linearly between +60°C and +70°C	
Output in battery mode:		
Output voltage	22.45V	At no output load
	22.25V	At 10A output load
Output current	10A	Below +60°C
	7.5A	At +70°C
	Derate linearly between +60°C and +70°C	
	15A	Short term, up to 5s
Power losses	5.5W	At 10A output load
Temperature range	-25°C to +70°C	
Max. wire size (litz wire)	4mm ²	For power terminals
Wire size AWG	AWG 20-10	For power terminals
Max. wire diameter	2.8mm	For power terminals
Wire stripping length	10mm/ 0.4inch	For power terminals
Max. wire size (litz wire)	1.5mm ²	For signal terminals
Wire size AWG	AWG 22-14	For signal terminals
Max. wire diameter	1.5mm	For signal terminals
Wire stripping length	6mm/ 0.25inch	For signal terminals
Tightening torque	0.4Nm/ 3.5lb.inch	For signal terminals
Size (wxhxd)	49x124x117 mm	Without DIN-rail
Weight	530g / 1.17lb	

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
Product Description

The battery module (77 3412 00) utilizes one non-spillable maintenance-free VRLA lead-acid battery, which is assembled in a package with all wiring and fusing included.

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
Installation Instructions

This device does not contain serviceable parts, except the fuse or replacement battery. If malfunction should occur during installation or operation, first check the fuse and battery for proper functionality. If the unit is damaged or if the malfunction continues, immediately turn power off and send device to Tunstall for inspection with a description of the fault.

Do not install the battery fuse until the installation work is complete.

Replacement fuses must be of the same type or an UL listed fuse with the same nominal value and characteristics.

Do not store batteries or battery modules at temperatures higher than 25°C or performance and lifetime may change.

Before initial use, check the last charging which is marked on the battery. If it is older than 9 month, replace the battery.

Do not store batteries longer than 9 months without charging or performance and lifetime may change. Every 9 months or whenever the open circuit voltage (OCV) falls below 12.6V, the battery module should be charged. Apply a charging voltage of 13.8V (at 25°C) for 72 hours. Batteries with an OCV below 12V should be disposed.

Do not short-circuit the battery poles or power wires and always check correct battery polarity.

Mount the device onto a flat surface. Do not install the battery module upside down.

Install the device in a cool area inside an enclosure providing protection against electrical, mechanical and fire hazards. Ensure sufficient ventilation by following the requirements of the EN 50272-2.

The device is designed for pollution degree 2 areas in controlled environments. No condensation or frost allowed.

The enclosure of the device provides a degree of protection of IP00.

Do not use batteries other than the original batteries supplied by the manufacturer to avoid conflicts with agency approvals.

The maximum surrounding air temperature is +60°C (+140°F). The operational temperature is the same as the ambient or surrounding air temperature and is defined 2cm below the device. Charging is allowed up to +50°C.

Frequently charging below +5°C reduces the lifetime of the battery.

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Instructions for replacement and disposal of batteries

Batteries need to be replaced periodically. The time interval depends on the ambient temperature in which the battery is used and on the number and depth of discharge cycles. In a typical application it is recommended to replace the batteries at least every 5 years. Disconnect battery fuse before replacing the batteries. Do not over-tighten screws of the terminal bolts of the battery. Recommended tightening torque is 3Nm (26 lb.inch).

Always dispose batteries through a recycling organization. Optionally, batteries can be returned to the supplier or the battery manufacturer. Batteries must be completely discharged before recycling.

CAUTION: The unit contains wet cell batteries. Ventilation of the installation room in accordance with the rules of Canadian Electrical Code, part 1 is required if used in Canada.

Installation Instructions for Hazardous Location Areas

The device is suitable for use in Class I Division 2 Groups A, B, C, D locations and for use in Group II Category 3 (Zone 2) environments.

Hazardous Location Classification: ATEX: EPS 15 ATEX 1 025 X, II 3G EX ec nC IIC T4 Gc

WARNING EXPLOSION HAZARDS!

Substitution of components may impair suitability for this environment. Do not disconnect the unit or remove or insert the fuse or battery unless power has been switched off or the area is known to be non-hazardous.

A suitable enclosure must be provided for the end product which has a minimum protection of IP23 and fulfils the requirements of the EN 60079-0. Sufficient ventilation must be ensured in the final installation.

Do not use batteries other than the original batteries supplied by the manufacturer or batteries, which are approved according to EN/ IEC 60079-7.

Technical data

All values are typical figures unless otherwise noted.

Voltage	DC 12V	
Nominal battery capacity	26Ah	
Battery current	Max. 30A	Discharge current
	Max. 5.5A	Charging current
Recommended charging voltage	13.9V	At +10°C
	13.75V	At +20°C
	13.6V	At +30°C
	13.45V	At +40°C
Self-discharge rate	3% / month	At +20°C
Battery fuse	30A	ATOF® 287 030 from Littelfuse.
Temperature range	-10°C to +50°C	For charging
	-15°C to +60°C	For discharging
Storage temperature range	-20°C to +50°C	
Max. wire size (litz wire)	4mm ²	Battery module terminal
Wire size AWG	AWG 12-10	Battery module terminal
Max. wire diameter	3.5mm	Battery module terminal
Wire stripping length	7mm / 0.28inch	Battery module terminal
Tightening torque	0.5Nm / 4.3lb.inch	Battery module terminal
Size (wxhxd)	214x179x153 mm	Without DIN-rail
Weight	10.1kg / 22.3lb	With battery
	700g / 1.54lb	Without battery

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
DC-UPS with capacitor storage 24V, 6kWs Order no. 77 3413 00

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Product Description


The 77 3413 00 is an uninterruptible power supply controller (DC-UPS) with included storage capacitors, which is used in combination with a 24V power supply to bridge power failures. When the power supply provides sufficient voltages, the DC-UPS controller charges the capacitors. When the power supply voltage fails, the energy stored in the capacitors is released to the DC bus in a regulated process.

Intended Use

This device is designed for installation in an enclosure and is intended for commercial use, such as in industrial control, process control, monitoring and measurement equipment or the like.

In systems where failure or failure of this equipment may cause serious injury or endanger human life, this equipment must be used only if technical and / or organizational measures are taken to correct the effects of the malfunction.

Functional Diagram

*) WARNUNG! For nurse call systems, the protective earth (PE) must not be connected to the DC-UPS. The DC-UPS must be mounted insulated from the protective earth (PE).

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The German standard DIN VDE 0834 for nurse call systems requires compliance with EN 60601-1 regarding electrical safety. This means: The 24V power supply unit must be earthed. All other devices of the modular power supply, e.g. the DC-UPS (77341300), must be mounted insulated from the protective earth (PE). Recommendation for insulated mounting: Install the DC-UPS on two DIN rail raisers HAGER HYA036 (order no. 76090004, width 106mm each). In addition, a distance of 1cm must be maintained from the insulated mounted device to all earthed surfaces and equipment.

Installation Instructions

Install the device in an enclosure providing protection against electrical, mechanical and fire hazards. Install the device onto a DIN-rail according to EN 60715 with the input terminals on the top of the unit. The input can be powered from a regulated power supply or a similar DC source. The voltage between the input and ground must not exceed 60Vdc continuously. (WARNING! For nurse call systems, the protective earth (PE) must not be connected to the DC-UPS.) The ripple voltage in the low frequency range between 50Hz and 10kHz must be negligible when used in marine applications.

The input must be powered from a PELV or SELV source or an "Isolated Secondary Circuit" in order to maintain a SELV or PELV output.

Use an appropriately sized 24V power supply, which can deliver the additional internal current consumption, required to charge the capacitors.

Check for correct input polarity. The device will not operate when the voltage is reversed.

Make sure that the wiring is correct by following all local and national codes. Use appropriate copper cables that are designed for a minimum operating temperature of 60°C for ambient temperatures up to +40°C and 80°C for ambient temperatures up to +60°C. Ensure that all strands of a stranded wire enter the terminal connection.

The device is designed for pollution degree 2 areas in controlled environments. No condensation or frost is allowed.

The device is designed as "Class of Protection III" equipment according to IEC 61140.

The enclosure of the device provides a degree of protection of IP20.

A disconnecting means shall be provided for the output of the device.

The device is designed for convection cooling and does not require an external fan. Do not obstruct airflow and do not cover ventilation grid!


Keep the following minimum installation clearances: 40mm on top, 20mm on the bottom, 5mm left and right side. Increase the 5mm to 15mm in case the adjacent device is a heat source. When the device is permanently loaded with less than 50%, the 5mm can be reduced to zero.

The device is designed for altitudes up to 6000m (19685ft).

The maximum surrounding air temperature is +60°C (+140°F). The operational temperature is the same as the ambient or surrounding air temperature and is defined 2cm below the device.

The device is designed to operate in areas between 5% and 95% relative humidity.

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WARNING EXPLOSION HAZARDS!

Substitution of components may impair suitability for this environment. Do not disconnect the device or change unit settings unless power has been switched off or the area is known to be non-hazardous.

A suitable enclosure must be provided for the end product which has a minimum protection of IP54 and fulfils the requirements of the EN 60079-0.

Functional Description

The output is electronically protected against no-load, overload and short circuit and can supply any kind of loads, including unlimited inductive loads and capacitive loads.

Do not apply return voltages from the load to the output terminals higher than 35V.

Green Status LED

This LED is on solid, when the capacitor is charged. During charging, the LED is flashing with a slow frequency and rapid flashing is an indication of a buffer event.

Yellow Diagnosis LED

This LED reports an overload situation, a high temperature, an expired buffer time or an activated inhibit input.

Red Warning LED

This LED indicates a failure in the input voltage or an activated PC-mode.

LED-Blinkmuster**Grüne Status-LED:**

Bereit

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Gelbe Diagnose-LED:

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Hohe Temperatur

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Inhibit aktiv

Rote Warnungs-LED:

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LED Signal Pattern**Green LED "Status":**

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EN - Installation Instructions

Ready and Buffering Relay Contacts

The Ready contact is closed when the capacitor is charged and the unit is ready to buffer.

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PC-Mode Input

The PC-mode always turns the output off for at least 5s after a buffer event lasting longer than 1s, independent of whether the 24V may have recovered during this time. This function ensures that the PC gets a restart signal. To enable a safe shut-down of the system, the forced turn off of the output is delayed with a constant time of 70s. To activate the PC-mode, connect the two pins marked with „PC-mode“ together on the signal connector together.

Technical data

All values are typical figures specified at 24Vdc input voltage, 15A output current in power supply mode at 25°C ambient, no charging and after a 5 minutes run-in time unless otherwise noted.

Input voltage	DC 24V	-20%/+25%
Minimum input voltage to start charging and to enable buffer mode	22.8Vdc	
Transfer voltage to switch into buffer mode	22.5Vdc	
Internal current consumption	1.1A	Includes charging current
Output in power supply mode:		
Input to output voltage loss	230mV	At 10A output load
	330mV	At 15A output load
Output current	15A	Continuous
Output in buffer mode:		
Output voltage	22.45V	At no output load
	22.25V	At 10A output load
Output current	15A	Continuous
Buffer time	340s	At 0.5A output current
	39s	At 5A output current
	16.5s	At 10A output current
	9s	At 15A output current
Charging time	16 minutes	For initial charging
	4 minutes 30s	For recharging shortly after a 5A buffer event
Power losses	4.6W	At 10A output load
Temperature range	-40°C to +60°C	
Max. wire size (litz wire)	4mm ²	For power terminals
Wire size AWG	AWG 20-10	For power terminals
Max. wire diameter	2.8mm	For power terminals
Wire stripping length	10mm/ 0.4inch	For power terminals
Max. wire size (litz wire)	1.5mm ²	For signal terminals
Wire size AWG	AWG 24-16	For signal terminals
Max. wire diameter	1.5mm	For signal terminals
Wire stripping length	8mm/ 0.3inch	For signal terminals
Size (wxhxd)	126x124x117 mm	Without DIN-rail
Weight	1150g / 2.54lb	

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EMC supplementary filter 24V DC, DIN rail, Order no. 77 3417 00

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- Turn power off before working on the device. Protect against inadvertent re-powering.
- Do not open, modify or repair the device.
- Use caution to prevent any foreign objects from entering the housing.
- Do not use in wet locations or in areas where moisture or condensation can be expected.

Intended Use

This EMC supplementary filter is designed as supplement to the power supply unit 10A, DIN rail (77 3410 00) for optimum limitation of conducted interference emissions to meet EN 61000-6-3.

Installation in direct vicinity to the power supply unit.

Installation Instructions

Install the device directly next to the power supply unit in the same enclosure providing protection against electrical, mechanical and fire hazards. Keep the minimum installation clearance of 15 mm to the power supply unit.


Install the device onto a DIN-rail according to EN 60715 with the input terminals on the bottom of the device.

Use shielded or unshielded cables, twisted or untwisted cables. Make sure that the wiring is correct by following all local and national codes. Use appropriate copper cables that are designed for a minimum operating temperature of 90°C for ambient temperatures between +45°C and +60°C and 105°C for ambient temperatures up to +70°C. Ensure that all strands of a stranded wire enter the terminal connection.

The housing of the device provides a degree of protection of IP20. The housing does not provide protection against spilled liquids.

The maximum ambient air temperature is +70°C. The device is designed to operate in areas between 5% and 95% relative humidity.

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

EN - Installation Instructions

19" Mounting Set, order no. 76 0900 01


Intended for the installation of up to two IP-SystemManagers or the installation of OSYlink modules in a 19" system cabinet.

Technical data

Dimensions (HxWxD)	132.5 x 482.6 x 180 mm
Type of mounting	19", 3 U
Top-hat rail	TS 35/7.5, length: 425 mm

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
- 1 - Frame
- 2 - Top-hat rail
- 3 - Front cover
- 4 - Four magnets for fixing the front cover
- 5 - Four screws (with washer and nut) for fixing the frame into the 19" system cabinet

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Mounting set for 10" wall mounting cabinet, order no. 76 0900 02

Mounting set for installation of equipment into 10" wall mounting cabinet. Sheet steel frame. Pre-mounted support rail for fixing the devices.

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- 1 - Frame
- 2 - 4 mounting holes for fixing the frame into the cabinet
- 3 - Supporting rail, 203 mm
- 4 - 4 pan-head screws, M6x16
- 5 - 4 cage nuts, M6

Mounting

- Fix the frame [1] at the four mounting holes [2] in the 10" cabinet using the pan-head screws [4] and the cage nuts [5].
- Clip the devices onto the supporting rail [3].

Scope of supply

- Frame
- Support rail, 203 mm, pre-mounted
- 4 pan head screws, M6x16
- 4 cage nuts, M6
- 4 cable ties, 180 mm

Technical data

Dimensions (HxWxD)	132.5 x 260 x 180 mm
Type of mounting	10", 3 U
Supporting rail	TS 35/7.5, length: 203 mm

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

C Kontroll-LEDs

OSYnet-Y-RepeaterOpto, order no. 77 4000 00

Data repeater for electrical isolation of the group bus (OSYnet). Intended for branching and/or isolation of the group bus OSYnet into 2 autonomous sections, Y-branching. One repeater can be used per physical group.

- Direct connection to the group bus (OSYnet) and its branch sections.
- Mounting on top-hat rail (35 mm), several units can be mounted side-by-side.
- Dimensions (HxWxD): 85 x 25 x 83 mm
- Power supply: 24 V DC (8 ... 30 V DC)
- Current consumption (at 20 °C): approx. 60 mA

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A Connection of OSYnet sections

For screw connection of the OSYnet sections you need three OSYnet connecting plugs, order no. 77 0950 00.

Connection point CAN Net 1 is electrically isolated from the other parts. That's why the OSYnet section, that should be electrically isolated from the others, has to be connected to connection point CAN Net 1.



For the pin assignment refer to the figure A.

B Connection of power supply


Connect the power supply cable to the connection clamp on the top **or** the bottom of the OSYnet-Y-RepeaterOpto (screw connection).

Attention! The power supply may not be connected through. I.e. it is not allowed to use one connection clamp as 24 V input and the other connection clamp as 24 V output e.g. for supply of further devices.

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



OSYnet-Gateway, order no. 77 4001 00


Active gateway for electrical isolation of the group bus (OSYnet). Intended for extension of the maximum cable length to additional 700 m. One gateway can be used per physical group.

- Direct connection to the group bus (OSYnet) and its branch sections
- Mounting on top-hat rail (35 mm), several units can be mounted side-by-side
- Dimensions (HxWxD): 105 x 37 x 35 mm
- Power supply: 24 V DC
- Max. current consumption (at 20 °C): 50 mA

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2) Note! The connection points DH and DL of OSYnet 1 as well as DH and DL of OSYnet 2 have to be connected by a terminating resistor (120 Ohm, order no. 00 0040 76).

B Control LEDs

- ☐ Power — LED on, if OSYnet-Gateway is supplied with power.
- ☐ OSYnet 1 Active — LED on for data traffic OSYnet section 1
- ☐ OSYnet 2 Active — LED on for data traffic OSYnet section 2

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© 2024 Tunstall GmbH