
- 76 2200 00 OSYnet-SafetyController 2
- 78 2200 00 OSYnet-SafetyController 2 SE

Installationsanleitung • Installation Instructions • Notice d'installation

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# Table of Contents

1. Intended use .....	9
2. Safety instructions .....	9
3. Product description .....	9
4. Mounting/Dismantling .....	11
5. 10-pole DIP switch .....	12
6. ESD protective cover .....	12
7. Connecting the cables .....	13
8. Control LED (red) .....	15
9. Technical data .....	15

## 1. Intended use

The devices are safety controllers for connection to an OSYnet group bus. They are used to maintain the basic system functions for an OSYnet group if the control unit (IP-SystemManager or OSY-ControlCenter) fails. According to the German standard DIN VDE 0834, this redundancy is required if the control unit is not installed locally on the ward or is supplied with a separate power supply.

OSYnet-SafetyController 2 (76 2200 00) is intended for connection to the OSYnet group bus of a Flamenco<sup>IP</sup> or Flamenco system (hospital).

OSYnet-SafetyController 2 SE (78 2200 00) is intended for connection to the OSYnet group bus of a Flamenco<sup>SE</sup> system (prison).

## 2. Safety instructions

### NOTICE

The complete installation of the system is described in the technical manual.

### CAUTION

The printed circuit board includes electrostatic sensitive components. Avoid touching.


## 3. Product description

If the communication between the OSYnet group bus and the control unit (IP-SystemManager or OSY-ControlCenter) fails, the OSYnet-SafetyController 2 (SE) automatically takes over the control of the devices connected to the OSYnet group bus. The OSYnet-SafetyController operation does not offer the full range of functions as the normal operation, i.e. operation with the IP-SystemManager or OSY-ControlCenter.

**Limitations during OSYnet-SafetyController operation:**

- Display and handling of calls only within the own OSYnet group.
- No connection to other OSYnet groups.
- No connection to devices connected via IP (corridor display<sup>IP</sup>, ComStation<sup>IP</sup> etc.).
- No connection to ManagementCenter.
- Display of call types, call locations and volume of the call tone as configured for the OSYnet-SafetyController 2 (SE) in the SystemOrganizer.
- Only one call is displayed at a time on ComStations. This is the oldest call with the highest priority (call category).
- No announcements, no calling from ComStations.
- No operation of group lamps or direction lamps.
- No indication of faults.
- Flamenco, Flamenco<sup>IP</sup> only: No shift operation, no ward coupling.
- No call recording.
- No operation of OSYlink-Universal.
- No operation fo ComStation<sup>BUS</sup> (77 0605 00).

## 4. Mounting/Dismantling

### 4.1. Mounting on top-hat rail

- Click the device onto the top-hat rail.

## 4.2. Wall mounting

1. Push out the three fixing clips until the holes for the screws appear.
2. Screw the device to the wall.

## 4.3. Dismantling from top-hat rail

Put a screw driver into the visible hole of the black fixing clip. Then pull the clip out upwards until the device detaches from the top-hat rail.

# 5. 10-pole DIP switch

## 5.1. Setting the address (126, 125)

For setting the address for the device on the OSYnet group bus use the code switches with the values 1, 2, 4, 8, 16, 32, 64. Set the appropriate switch to ON. Other addresses are set through addition.

You need to set address 126. If you install a second OSYnet SafetyController 2 (SE) on the same group bus, you must set address 125 on it.

If the address is changed during operation, a reset must be carried out afterwards. Otherwise the new setting will not be applied.

## 5.2. Reset


To carry out a hardware reset, set the coding switch 10 (Reset) to ON for one second and then back to OFF.

# 6. ESD protective cover


The ESD protective cover (00 0276 53) serves as protection against electrostatic discharges and as protection against contact.

- After connecting the cables (see Section 7: “Connecting the cables” (page 13)), put on the ESD protective cover.

## 7. Connecting the cables

Connection terminals: Pluggable screw terminals, stripping length: 7 mm. Power supply connection: 1 wire per terminal; Other connections: 2 wires per terminal.

<b>OSYnet group bus</b> Cable: CAT7 (22 AWG) $\varnothing = 0.64 \text{ mm}$ or CAT5, shielded (23 AWG) $\varnothing \text{ min. } 0.57 \text{ mm}$	DL — white		from previous and to next OSYnet user
	DH — blue		
	DG — green		
	DG — white		
	NF 1a — brown		
	NF 1b — white		
	NF 2a — orange		
	NF 2b — white		
<b>OSYnet group bus</b> Cable: J-Y(St)Y 4x2x0.8	DL — red		from previous and to next OSYnet user
	DH — blue		
	DG — green		
	DG — white		
	NF 1a — brown		
	NF 1b — white		
	NF 2a — yellow		
	NF 2b — white		

**Table 1. Output for signalling the operation mode**

Characteristic	Meaning	Setting
Function	The output signals whether the group is controlled by the IP-System-Manager/OSY-ControlCenter (= normal operation) or by the OSYnet-SafetyController (= OSYnet-SafetyController operation).	
Type of output	Potential-free digital output, changeover contact.	
Contacts	Signalling: Normal operation.	OUT 1 closed, OUT 2 open.
	Signalling: OSYnet-SafetyController operation.	OUT 1 open, OUT 2 closed.
Jumper	Power supplied by the nurse call system (24 V DC, max. 60 mA for controlling a relay)	Jumper OUT 1/2 inserted.
	External power supply (24 V AC or DC, max. 1A)	Jumper OUT 1/2 not inserted.

**Table 2. Output for fault signalling**

Characteristic	Meaning	Setting
Function	The output signals whether a user on the OSYnet group bus has reported a fault.	
Type of output	Potential-free digital output, changeover contact.	
Contacts	Signalling: No fault.	OUT 3 closed, OUT 4 open.
	Signalling: Fault!	OUT 3 open, OUT 4 closed.

Characteristic	Meaning	Setting
Jumper	Power supplied by the nurse call system (24 V DC, max. 60 mA for controlling a relay)	Jumper OUT 3/4 inserted.
	External power supply (24 V AC or DC, max. 1A)	Jumper OUT 3/4 not inserted.

## 8. Control LED (red)

**Table 3. Signaling of the control LED (red)**

Signaling	Means
Does not light up.	OSYnet SafetyController 2 (SE) without power supply or defective.
Flashing slowly (1 s ON - 1 s OFF).	Normal operation = IP-SystemManager or OSYnet-SafetyController is controlling the group. OSYnet-SafetyController 2 (SE) ready for operation.
Flashing quickly (200 ms ON - 200 ms OFF).	OSYnet-SafetyController operation = OSYnet-SafetyController 2 (SE) is controlling the group.
Continuous light longer than 10 s.	OSYnet-SafetyController 2 (SE) not ready for operation.

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© 2023 Tunstall GmbH
Tunstall GmbH • Orkotten 66 • 48291 Telgte • Germany • [www.tunstall.de](http://www.tunstall.de)