Installationsanleitung • Installation Instructions

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Table of Contents

1. Intended Use	9
2. Safety notes	9
3. Product description	9
4. Mounting	10
4.1. Mounting the pattress	11
4.2. Mounting the front panel	11
5. Connecting	12
6. Dismantling	14
7. Technical data	14

1. Intended Use

Table 1. Products

76 0550 00	ConnectionTerminal IP, PoE with VoIP and bedside speech
76 0550 10	ConnectionTerminal IP, 24V DC with VoIP and bedside speech

The ConnectionTerminal IP is a room terminal. It is designed for connection to the IP network of a nurse call system of the Flamenco system family. The RAN room bus with the room devices is connected to the ConnectionTerminal IP .

ConnectionTerminal IP, PoE (76 0550 00) is powered via PoE+.

ConnectionTerminal IP, 24V (76 0550 10) is connected to the ward's 24 V DC power supply.

2. Safety notes

NOTICE

The complete installation of the system is described in the Technical Manual.

CAUTION

The printed circuit board includes electrostatic sensitive components. Avoid touching.

3. Product description


The ConnectionTerminal IP controls and monitors the functions of the nurse call system for one room. It has no display or control elements and is therefore installed out of arm's reach of nursing staff and patients. It must be accessible to authorised persons at all times (minimum maintenance access width of 60 cm). Fixed installation, for example in the cavity above the suspended ceiling. Sufficient heat dissipation must be ensured.

The RAN room bus is connected to the ConnectionTerminal IP. All devices in the room are connected to the RAN room bus. In addition to data transmission, the ConnectionTerminal

also supports bedside voice communication. Patient control devices are required for voice communication (Patient handset or ePat[®]lite).


4. Mounting

Figure 1. Mounting / Dismantling

- [1]^a Back box
- [2] Pattress
- [3] Four back box screws
- [4] Short patch cable
- [5] Hinge
- [6] Front panel
- [7] Catch
- [8] Security screw (M3, 10 x 3 mm, cross-head PZ)

^aNot included, order separately.

Back box installation

When installing the back box consider the space required for the ConnectionTerminal IP. Measuring from the centre of the box, the available space must be at least 120 mm above and below.

4.1. Mounting the pattress

See Figure 1 A.

Wall mounting on two-gang back box. *)

Preparation: Ideally, insert the Ethernet cable into the back box from below and use a Keystone module with RJ45 jack for the connection. The RJ45 jack must be connected to the RJ45 socket on the connection field of the Connectionterminal IP via the short patch cable. The keystone module is not included in the scope of delivery of the ConnectionTerminal IP.

1. Insert one plug of the short patch cable into the RJ45 jack of the Keystone module. Prepare the wires of the other cables gemäß Section 5: "Connecting" (page 12).
2. Feed the short patch cable and the prepared wires through outlet in the pattress.
3. Fix the pattress using the back box screws to the back box.
4. Connect the connection wires to the connection field in the pattress according to Section 5: "Connecting" (page 12).
5. Insert the second plug of the short patch cable into the RJ45 socket on the connection field.

*) **NOTE!** If mounting on a back box is not possible, the ConnectionTerminal IP may be fixed to the wall with dowels and screws.

4.2. Mounting the front panel

See Figure 1 B and C.

1. Place front panel onto the top of the pattress creating a hinge.
2. Press the front panel downwards onto the pattress with low pressing power until the catch engages.
3. If necessary, the security screw can be fitted into the hole from below (see figure).


IMPORTANT

The ConnectionTerminal IP must be configured. The configuration is carried out using the SystemOrganizer configuration software.

5. Connecting

Figure 2. Connection field

[9] Connection field

[10] IP network of the nurse call system

[11] RAN room bus

[12] 24 V DC power supply (only for 76 0550 10)

[13] LED lights up green when PoE+ voltage is applied (only for 76 0550 00)

24 V DC power supply (only for 76 0550 10)

Connector	Screw terminal, 4-pole
Screw	M2.5
Cable type	NYM 2x2.5mm ²
Wire stripping length	6 mm
Connection capacity	0.05 – 3.31 mm ²

IP network of the nurse call system

Connector	RJ45 socket
Cable type for 76 0550 00 (PoE+)	Min. CAT5e, shielded, Ø min. 0.64 mm (22 AWG)
Cable type for 76 0550 10	Min. CAT5e, shielded

The Ethernet port is internally equipped with a separation point with 2 x MOPP according to DIN EN 60601-1 (data and power supply). Therefore, no external separator is required for the terminal 76 0550 00 with PoE+ supply. The 24 V DC connection of the terminal 76 0550 10 has no such separation point. Therefore according to the German standard DIN VDE 0834, systems with 76 0550 10 must be constructed according to the principle of system isolation.

RAAN room bus	
Connector	Screw terminal, 8-pole
Screw	M3
Cable type for room bus without speech	J-Y(St)Y 2x2x0.8 or CCS32. Optional: Min. CAT5, shielded, Ø min. 0.57 mm (23 AWG).
Cable type for room bus with speech	2x J-Y(St)Y 2x2x0.8 or CCS32. Options: 2x CAT5/CAT6, shielded, Ø min. 0.57 mm (23 AWG) or 1x CAT5/CAT6 S/FTP, Ø min. 0.57 mm (23 AWG) or 1x CAT7, Ø min. 0.57 mm (23 AWG).
Wire stripping length	6 mm
Connection capacity	0.20 – 2.50 mm ² , i.e. max. 4 wires per connection point

Important information for connecting the room bus!


- The cable for speech transmission, i.e. the cable containing the wires *MIC*, *LS*, and *SCH-B*, may be stripped for a max. of 30 mm only, to ensure the protective screen function. The shield wire must be connected to the connection point *SCH-B*.
- The cable shield and the shield wire of the cable containing the wires *RAN*, *+24V-Si* and *0V* are not connected and must be removed in the back box up to the cable sheath.
- The unused wires must not be interconnected.
- Unused screws must be screwed in.

Figure 3. Connections

Power supply

(only for 76 0550 10)


Cable: NYM 2x2.5 mm²

RAAN room bus with speech

Cable: 2x J-Y(St)Y 2x2x0.8


(max. 4 wires per connection point)

RAAN room bus without speech

Cable: J-Y(St)Y 2x2x0.8

(max. 4 wires per connection point)

6. Dismantling

See Figure 1 C.

1. If a security screw is fitted into the base of the ConnectionTerminal IP, remove it.

2.  **CAUTION**
The front panel is fastened only with the catch.

Use a slotted screwdriver to push the catch upwards until the front panel loosens.

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© 2024Tunstall GmbH
Tunstall GmbH