



# Flash Receiver BE1441


# Flash Receiver BE1442/1443


# Battery backup

# BE1442

## Getting started

- Remove the battery slip on the bottom of the unit before use.
- Please charge the unit for 24 hours to ensure full backup time.


## Replacing batteries (Done by the operator)

- Remove the screw on the bottom of the unit to exchange the batteries.

### CAUTION

RISK OF EXPLOSION IF BATTERY IS REPLACED BY AN INCORRECT TYPE.

DISPOSE OF USED BATTERIES ACCORDING TO THE INSTRUCTIONS.

- Use only 4x1.2VNiMH 600AAA/900AAA Lexel rechargeable batteries to ensure full backup time.

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### Settings

No adjustments are required for normal use. The relevant descriptions are provided below, if you wish to change a setting for some reason.

### Radio key

In order to use several Bellman Visit Systems close to one another without interference, different Radio Keys can be set on the different systems. All Bellman Visit System units are supplied from the factory tuned to the same Radio Key, channel 0. This means that all Radio Key Switches on the transmitters are set to the OFF position.

A Bellman Visit transmitter is also required to alter the radio key on the Bellman Visit Flash Receiver. Proceed as follows to change the Radio Key:

- Set a Bellman Visit Transmitter to the desired Radio Key by altering its Radio Key Switch to the desired Radio Key. See the transmitter user manual for further information.
- Hold down the Test Button (9) on the Bellman Visit Flash Receiver until the green  and yellow  LEDs blink alternately. The Bellman Visit Flash Receiver will now be in programming mode for about 30 seconds.
- Press the Bellman Visit transmitter's test button.
- The LEDs     will now blink five times in quick succession to indicate a successful change of radio key.
- After changing the radio key, the Bellman Visit Flash Receiver will automatically return to normal mode.

**Please note:** all Bellman Visit products within a system must be tuned to the same Radio Key in order to operate as a group.

### Signal pattern

A Signal Pattern is the name for the way in which a receiver in the Bellman Visit System indicates activation. Changing the transmitters' Signal Switch changes the Signal Pattern which the receivers display when the transmitter is activated.

The following signal patterns are available for the Bellman Visit System: (Please note that the Flash Receiver does not emit any sound):

## 868/315MHz Europe, Asia

Type	LED pattern	Sound	Vibration	Flash
Green 1	Green is constantly lit	1 x ding dong, low-frequency tone	Separate	Yes
Green 2	Green blinks in sequences of two	2 x ding dong, low-frequency tone	Separate	Yes
Green 3	Green blinks in sequences of three	1 x ding dong, high-frequency tone	Separate	Yes
Green 4	Green blinks constantly	2 x ding dong, high-frequency tone	Separate	Yes
Yellow 1	Yellow is constantly lit	1 x ring, low-frequency tone	Short	Yes
Yellow 2	Yellow blinks in sequences of two	2 x ring ring, low-frequency tone	Short	Yes
Yellow 3	Yellow blinks in sequences of three	1 x ring, low-frequency tone	Short	Yes
Yellow 4	Yellow blinks constantly	2 x ring ring, high-frequency tone	Short	Yes
Orange 1	Orange is constantly lit	Baby	Rapid	Yes
Orange 2	Orange blinks in sequences of two	Baby	Rapid	Yes
Orange 3	Orange blinks in sequences of three	Baby	Rapid	Yes
Orange 4	Orange blinks constantly	Baby	Rapid	Yes
VMA	Red and Orange constantly blink alternately	VMA constant	Long	Yes
Fire alarm	Red blinks constantly	Fire alarm constant	Long	Yes

Hereby, Bellman & Symfon, declares that this wrist receiver and charger is in compliance with the essential requirements and other relevant provisions of directives: R&TTE 1999/5/EC, LVD 1973/23/EC, EMC 1989/336/EC and MDD 1993/42/EEC.

**For use in all EU countries**

Correct use and disposal of batteries. Replace only with the same or equivalent type recommended by the manufacturer. Please dispose of old, defective batteries in an environmentally friendly manner in accordance with the relevant legislation. Dispose used batteries at municipal collecting points or dispense them to local stores free of charge.

**Contact us**

For a complete Declaration of Conformity please contact [info@bellman.com](mailto:info@bellman.com) or call +46 31 68 28 20

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Bellman & Symfon Europe AB
Bellman & Symfon Asia Ltd.